



"VIRGO" Five Function Electric Bed

REF:0803152



Technical Specifications

- Metal Frame
- Metal Rails in Wood Colour
- Internal Wooden Covers for Panels
- Featuring Trendelenburg & Anti-Trendelenburg 12° Positions
- Product Weight: 90kg
- Max Back Tilt: 70°
- Max Leg Tilt: 30°
- Power Input 220V-50Hz, Output 24V
- Folding, Easy to Transport, Assemble & Disassemble by a Single Person
- Robust Construction
- Wooden Orthopedic Bed Frame
- Electric Elevation of Back, Knees & Feet
- Includes Lifting Pole
- High-quality Solid Casters with Brakes
- Safety Rails that can be Installed Without Tools



Int. Dimensions (LxW)	Ext. Dimensions (LxW)	Height from floor	Max Load Capacity	Back / legs / knees Elevation
203x90cm	215x102cm	Electric Elevation 39-80cm	180Kg	Electrical



Backrest section adjustment



Thigh and calf section adjustment



Trendelenburg position's angle 12



Reverse Trendelenburg position's angle 12



Height adjustment 39-80cm



MOBIAK[®]
MEDICAL • FIREFIGHTING • GASES est. 1977

www.mobiakcare.com

"LIBRA" Five Function Electric Bed

REF:0803151



Technical Specifications

- Metal Frame
- Metal Rails in Wood Colour
- External Wooden Panel Covers, covering the Elevation Mechanisms
- Featuring Trendelenburg & Anti-Trendelenburg 12° Positions
- Product Weight: 90kg
- Max Back Tilt: 70°
- Max Leg Tilt: 30°
- Power Input 220V-50Hz, Output 24V
- Folding, Easy to Transport, Assemble & Disassemble by a Single Person
- Robust Construction
- Wooden Orthopedic Bed Frame
- Electric Elevation of Back, Knees & Feet
- Includes Lifting Pole
- High-quality Solid Casters with Brakes
- Safety Rails that can be Installed Without Tools
- External Wooden Panel Covers, covering the Elevation Mechanisms



Int. Dimensions (LxW)	Ext. Dimensions (LxW)	Height from floor	Max Load Capacity	Back / legs / knees Elevation
203x90cm	203x90cm	Electric Elevation 39-80cm	180Kg	Electrical



Backrest section adjustment



Thigh and calf section adjustment



Trendelenburg position's angle 12



Reverse Trendelenburg position's angle 12



Height adjustment 39-80cm



www.mobiakcare.com

"DIAS" the Premium Home Care Bed

REF:0805600



Technical Specifications

- Metal Frame
- Metal Rails in Wood Colour
- Wooden cover head/foot, covering the Elevation Mechanisms : Optional (ref: 0805601)
- Featuring Trendelenburg & Anti-Trendelenburg 12° Positions
- Max Back Tilt: 72°
- Max Leg Tilt: 21°
- Power Input: 100–240V
- Folding, Easy to Transport, Assemble & Disassemble by a Single Person
- Safe working load for bed: 215 kg
- Robust Construction
- Wooden Orthopedic Slats
- Electric Elevation of Back, Knees & Feet
- Includes Lifting Pole
- High-quality Solid Casters with Brakes
- Safety Rails



Transport
Mechanism

Int. Dimensions (LxW)	Ext. Dimensions (LxW)	Height from floor	Max Load Capacity	Back / legs / knees Elevation
200x90cm	215x102.5cm	Electric Elevation 39-80cm	178Kg	Electrical



Backrest section
adjustment



Thigh and calf
section
adjustment



Trendelenburg
position's
angle 12



Reverse
Trendelenburg
position's
angle 12



Height
adjustment
39-80cm