

LION | Aluminium Wheelchair

Series



MOBIAK[®]
MEDICAL • FIREFIGHTING • GASES est. 1977

SPECIFICATIONS

- Aluminium - Folding Frame
- Silver
- Fabric with Lining Seat
- Detachable, Adjustable
- Folding Armrests
- 24" Solid Rear Wheels with Radial Wheels & Quick Release
- 8" x 2" Solid Front Wheels



Ref: 0810804

Ref: 0810805

Ref: 0810806

Ref: 0810830

0810804	42cm	41cm	42cm	42cm	62cm	103cm	92cm	49cm	24" 8"	135kg	16,5kg
0810805	44cm	43cm	42cm	42cm	63cm	103cm	92cm	49cm	24" 8"	135kg	16,5kg
0810806	46cm	45cm	42cm	42cm	66cm	103cm	92cm	49cm	24" 8"	135kg	16,5kg
0810830	49cm	48cm	42cm	45,5cm	67cm	99cm	90cm	46cm	24" 8"	135kg	17,5kg

2 / Years
Warranty

www.mobiakcare.com



Value in safety, invest in Quality!