

Mesh Nebulizer Ref. 0803431



The nebulizer 0803431 is a nebulization therapy device that works with new generation ultrasound technology and is designed for inhalation of almost all kinds of drugs. It is the ideal choice for adults and children, both for home and emergency treatment, travel or outdoor use.



Medical grade



Nebulizer rate adjustable



Soundless



Portable



Small nebulizer particles

Specifications:

- Input power: DC 3.0V (two "AA" 1.5V alkaline batteries)
- Operation noise: approx. 50dB(A)
- MMD: 3.49m
- Residual Volume: 0.25ml
- Weight: 108g (without batteries)
- Operation: Continuous
- Battery duration: ≥60 min
- Dimension: 67 x 48 x 125mm
- Warranty: 2 years

Package included:

- Medication Chamber
- Pediatric Mask
- Adult Mask
- Mouthpiece
- Batteries x 2
- User Manual
- Carrying Case

