

Oximeter

Ref. 0803440

○ *Finger pulse oximeter measuring O₂ (SpO₂) and heart rate.*



Automatic
shutdown



Pulse column
display

Specifications:

- Measurement range of pulse rate: 25bpm-250bpm
- Accuracy of SpO₂: 70% ~100%: ±2%
- Accuracy of pulse rate: ±1%/ ±1bpm
- Power consumption: 40mA
- Operation environment: Temperature(°C): 5~40, Relative humidity: 80%
- Boot method: Push button to boot
- Power supply: 2 × 1.5V AAA alkaline batteries
- Display: LED
- Dimension: 58 × 33 × 33.5mm
- Weight: 35g (without batteries)
- Automatic shutdown time: 8s
- Warranty: 2 Years



Package included:

- Batteries x 2
- User Manual