

Compressor nebulizer Ref. 0811850

Your Nebulizer is intended for use in the treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy.



Specifications:

- Mask with Adjustable Nebulization
- Voltage: 220V
- Power: 168Watt
- Liter flow range: 8.5 ± 0.5 lpm
- Pressure: 2.5 ± 0.5 bar
- Average nebulization rate: ≥ 0.25 ml/min
- Sound level: ≤ 65 dB
- Operation: Continuous
- Weight: 1.25kg
- Dimension: 180x150x80mm
- Warranty: 3 Years

Package included:

- Air tube
- Nebulizer Tank
- Mouthpiece
- Angular connector for mask use in supine position
- Adult Mask
- Pediatric Mask
- 5 pcs extra filters
- User Manual