

Automatic Blood-Pressure monitor

Ref. 0803441



- Digital blood-pressure monitor that measures systolic pressure (SYS), diastolic pressure (DIA) and heart rate.

Specifications:

- 80 groups of memory
- Average of three measurements
- Systolic/Diastolic and Pulse indication
- Measuring range: Diastolic: 20-210 mmHg, Cuff pressure: 0-300mmHg, Systolic: 40-260mmHg, Pulse: 40-200 bpm
- Dimension: 130 x 96 x 75mm
- Weight: 251g (without batteries)
- Power supply: 220-240V or 4 AA 1.5V batteries
- Warranty: 3 Years



Package included:

- Cuff
- USB cable
- User Manual