



The main advantage of scooters like "APOLLO" is that they allow their users to regain their independence. People who previously relied on others for transportation can now move around with comfort and autonomy. Regardless of their destination, they can go to the doctor, visit friends, go shopping and enjoy life without restrictions.



## Specifications

- Colors Available: Black-Gold
- Motor Output W: 270W,
- Motor Rated Power W: 700W
- Safe Climbing Angle: 6 degree
- Handle Bar Type Or Joystick: Delta handle bar
- Electromagnetic Brake
- Emergency Hand Brake
- Seat Height (Adjustable)
- Seat Belt
- Adjustable Seat (Forwards and Backwards)
- Light Reflectors Front And Rear
- Basket Included
- Horn
- Product Frame: Steel
- Battery Warranty: 6 Months



## Technical data



1300mm



860mm



570mm



8'' 8''




440 x 430 x 400mm



44Kg



120Kg



AGM 24V, 12 Ah



6°



8km/h



15km

Operating console



Swivel Seat

