



Mobility scooters like **"ARES"** offer an essential blend of independence, convenience, and comfort for individuals with limited mobility. Designed to empower users to move freely, they open up possibilities for navigating the world with confidence and ease.



Specifications

- Colors Available: Silver
- Motor Output W: 340W,
- Motor Rated Power W: 700W
- Turning Radius: 1400mm
- Safe Climbing Angle: 9 degree
- Handle Bar Type Or Joystick: Delta handle bar
- Electromagnetic Brake
- Emergency Hand Brake
- Headrest (Adjustable)
- Seat Height (Adjustable)
- Seat Belt
- Adjustable Seat (Forwards and Backwards)
- Turning Lights & Signal
- Reverse Signal
- Light Reflectors Front And Rear
- Alarm Signal (When Both Trurning Lights Flash Simultaneously)
- Basket Included
- Mirrors Included
- Horn
- Product Frame: Steel
- Battery Warranty: 6 Months



Technical data



1255mm



965mm



630mm



10" 10"



450 x 460 x 490mm



83Kg



150Kg



(12Vx2) 52Ah



9°



12km/h



50km

Lifting Handrests



Frontlights



Folding Backrest



Handbrake



Swivel Seat



Folding Steer



Anti-Tilt Wheels And Backlights



Front and Rear Suspensions

