

Oxygen Concentrator

Ref.2 Liters 8F-2A

Ref. 0803420

○ Small in size and weight 2 liter oxygen concentrator with modern, ergonomic design and low noise operation.



Failure detection



Medical grade



High oxygen concentration



Fashion appearance

Specifications:

- Oxygen Purity (%): $\geq 90\%$
- Oxygen Flow (Lit/min): 0-2
- Max Outlet Pressure (bar): 0,5
- Power Supply (VAC): $230V \pm 10\%$, $50Hz \pm 1Hz$
- Power (W): 150
- Noise Level dB(A): 60
- Dimensions (cm): L:30 × W:18 × H:38
- Net weight (Kg): 9
- Warranty: 3 years

Package includes:

- Power Supply Cable
- Nasal Cannula
- Humidifier
- Filter
- Firetrap
- 2x fuses
- Rubber for sealing the humidifier
- User Manual

