

# "Gemini" Wheelchairs

Ref: **0811306**

41cm



120kg

Ref: **0811301**

43cm



130kg

Ref: **0811300**

46cm



130kg

Ref: **0811610**

48cm



130kg

Ref: **0811611**

51cm



130kg



2 years  
guarantee



Double Cross Support Bar  
\*Except Ref.: 0811306



DIMENSIONS AND WEIGHT				
Measurement in cm				
0811306	43	41	60	120kg
0811301	44	43	61	130kg
0811300	47	46	64	130kg
0811610	49	48	66	130kg
0811611	52	51	68	130kg

### Technical Specifications

Frame / Color	Steel - Folding / Silver
Seat & Backrest Fabric	Fabric with Lining, with Fire Retardant Properties
Armrests	Folding
Footrests	Detachable & Adjustable



**Accessories**



Ref.: 0811357  
 Anti-Tilt Wheels



Ref.: 0811358  
 Right Footrest

Ref.: 0811359  
 Left Footrest



Ref.: 0811363  
 24'' Pneumatic  
 Spoke Wheel



Ref.: 0811369  
 Walking Stick  
 Holder



Ref.: 0810794  
 IV Stand



Ref.: 0810767  
 Fabric Bottle  
 Case



Ref.: 0811373  
 Cylinder Base



Ref.: 0811376  
 Cup Holder





Real time  
crash test  
images



**Our wheelchairs 0811300, 0811301,0811306,0811610 & 0811611 have successfully passed the crash test according to ISO 7176-19:2022.**

This means that each of this wheelchair is safe to be pushed inside a suitable vehicle and secured during travel and has anchor points, allowing the wheelchair to be secured into place with the occupant seated within the wheelchair.

The anchor points on the wheelchair are used by attaching straps to the anchor points on the wheelchair. These straps are then tightened to ensure the wheelchair will not move in the event of an accident and therefore protecting the wheelchair user from harm when traveling.

