

Auto-CPAP YH-550

Ref. 0803370

Auto-CPAP YH-550 GSM

Ref. 0803432

Remote control capability!

GSM for remote control, data export,
and setting adjustments



Product Parameters

Model: YH-550/ YH-550 GSM

Color: White

Mode: CPAP, Auto CPAP

Pressure: 4-20cm H₂O

Noise: 28dBA±2dBA

G.W./N.W.: 2.9/1.5 Kgs

Dimensions (LxHxW) cm: 41.5 × 30 × 16

○ The Auto-CPAP offers modern and flexible APAP therapy with precise event detection and a high level of comfort. Innovative design, integrated humidifier.

Technical Specifications

- Automatic turn-on/off
- Smart START/STOP
- Exhalation assistance FPS, 4 levels 0-3
- Ramp 0-45 minutes
- Leak alarm
- Silent operation 28dBA ±2 dBA
- Humidifier with 6 adjustable levels
- For detailed technical specifications, refer to the user manual

Package Contents:

- Humidifier
- Breathing tube 22mm
- Nasal mask
- Power supply
- Power cord
- Set of 2 filters
- Carrying bag
- SD card
- User manual

Warranty:

- 2 years for the device
- 1 year for the power supply
- 90 days for the humidifier