

- It is powered by a 5 V power adapter or a 3.7 V lithium battery, with the advantages of compact and portable design, easy to use, adjustable (two grades) and high in nebulization rate, which is applicable for use in hospitals, clinic and family.

Mesh Nebulizer

Ref. 0803845



Features:

- One-button operation
- Liquid nebulization function
- Power indication function
- Low-battery indication function
- Automatic shutdown without liquid
- Adjustable nebulization rate
- Charging function

Operating environment:

- Temperature: 5 °C ~ 40°C
- Relative humidity: 15%-90%
- Atmospheric pressure: 700 hPa ~ 1060 hPa

Package included:

- USB cable
- Adult mask
- Child mask
- User manual




Small in particle size



Compact and convenient to carry



One-Button operation



Low Noise



High nebulization rate



Automatic power off

Performance:

- Median particle diameter (D 0.50): 1-4µm (±25%)
- Equivalent volume particle diameter distribution: the proportion of particle diameter below 5 µm is not less than 60%
- Input power: < 5 VA

Nebulization rate:

- Grade I: 0.25 mL/min;
- Grade II: ≥ 0.15 mL/min
- Cup capacity: 10 mL
- Liquid volume: 2~ 10 mL
- Residual volume: ≤ 0.3 mL
- Noise: ≤50 dB (A)
- Working frequency: 110 KHz ± 10%
- Power supply: DC 5V or 3.7 V rechargeable lithium battery
- Safety class: class II equipment, type BF applied part
- The degree of waterproof and dustproof: IP22
- Mode of operation: continuous operation
- Dimension: 67.42 mm(L) × 48.53 mm(W) × 116.92 mm(H)
- Weight: about 0.1kg

Warranty:

- 2 years on the device
- 6 months on the drug chamber



www.mobiak.com



mobiakcaresales@mobiak.com | medicalexports@mobiak.com