

Oxygen Concentrator 5lt 8F-5A

Ref. 0803300



Technical Specifications

- Oxygen Purity (%): 87–96
- Oxygen Flow (Lit/min): 0.5–5
- Max Outlet Pressure (bar): 0.7
- Power Supply (VAC): 220
- Power (W): 400
- Noise Level dB(A): 48
- Dimensions (cm): 36×23×50 (LxWxH)
- Net Weight (Kg): 15.9

Package Contents

- Remote Control
- Power Cord (EU Type)
- Air Intake Filter
- User Manual

Alarms

- Power Interruption
- Low O₂ Purity
- Device Overheat
- High / Low Pressure