

**MOBIAK**<sup>®</sup>  
FIREFIGHTING • GASES • MEDICAL est. 1977

# Scooter PERSEUS

Ref: 0805760



The **PERSEUS** electric scooter is a reliable solution for the daily mobility of people with limited mobility. With its modern design, excellent autonomy and advanced safety features, **PERSEUS** offers independence and comfort without compromise.



2 years  
guarantee



- 1. Control console
- 2. Handle
- 3. Steering column

- 4. Batteries box
- 5. Front wheel
- 6. Foot platform

- 7. Seat and backrest
- 8. Manual lever
- 9. Rear wheel/drive wheel



Rotary switch  
for speed control

Button for releasing the  
electromagnetic brake

Power switch

2 years  
guarantee

Steering wheel  
adjustment lever

Charging socket

2 USB ports for  
charging devices  
(e.g. Smartphones)



## 2 years guarantee



Backlights

Anti-tilt wheels

### Specifications

- **Frame:** Steel
- **Colour:** White
- **Seat:** Swivel, height-adjustable and depth-adjustable (back and front)
- **Armrests:** Raised and width adjustable
- **Maximum load:** 190 kg
- **Maximum speed:** 15 km/h
- **Autonomy:** 55 km
- **Maximum obstacle height:** 60 mm
- **Maximum safe slope:** 15°
- **Steering Wheel Inclination:** up to 135°
- **Reclining Backrest:** up to 125°
- **Dimensions (overall) (L x W x H):** (1330 x 670 x 1350)mm
- **Weight of Batteries:** 36Kg
- **Net weight:** 106 kg (without batteries)
- **Effective seat depth:** 500 mm
- **Effective seat width:** 480 mm
- **Seat surface height:** 650 – 750 mm
- **Foot platform height:** 170 mm
- **Motor:** 800W
- **Batteries:** 2 pcs – 55Ah (lead acid) 6 Months Warranty
- **Charger:** 24V x 6A
- **Protection class:** IP54
- **Brakes:** Electromagnetic
- **Rear wheel diameter:** 13" with air pressure
- **Diameter of front wheels:** 11.5" with air pressure
- **Suspension:** Front & Back



Measurement  
in cm



46



48



50



65,5



58,3



133



135



65-75



11,5" 13"



6



*Value in Safety... Invest in Quality!*