

Suction Unit 7A-23D (20 L/MIN)

Ref. 0803434

2 years
guarantee



○ Effective and ergonomic suction device designed for safe removal of secretions in clinical and home environments. Ideal for healthcare professionals and home care, it offers reliable operation with easy transport and use.

Specifications:

Suction 7A-23D, 20Lt/min (Electric)

- Power supply: AC230V, 50Hz
- Max negative pressure: 90±10KPA
- Operation noise: ≤60 dB(A)
- Input power: 180VA
- Max air flow: ≥20L/min
- Liquid collection jar: 2500mL/pc*2

Warranty:

- 2 Years for the Suction

Dimensions:

- G.W.:15.2kg
- N.W.:13.5kg
- Size: 36×32×48cm

Package Contents:

- Abdominal suction tube: 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- Flexible suction tube (Φ7×Φ12, length 2 m): 1 pc
- Power cord: 1 pc
- Foot switch: 1 pc
- User manual: 1 pc