

Suction Unit 7E-A (20L/MIN)

Ref. 0803436

2 years
guarantee



Specifications:

Suction 7E-A, 20Lt/min (Electric)

- Input power: 90VA
- Power supply: AC230V, 50Hz
- Max Negative Pressure: 85±10KPA
- Max air flow: ≥20L/min (without jar)
- Operation noise: Approximately ≤65dB (A)
- Liquid collection jar: 1000mL
- Operating time: 30 minutes on, 30 minutes off

Dimensions:

- G.W : 7kg
- N.W.: 4kg
- Size: 28×19,6×28,5cm

- *Lightweight and quiet suction device, ideal for basic medical use and home care. Offers consistent performance, easy operation, and includes a fluid collection container. Suitable for daily support of patients with respiratory needs.*



Package Contents:

- Suction tube (length 2 m, Φ7×Φ12): 1 pc
- Suction catheter F8 (for children): 1 pc
- Suction catheter F12 (for adults): 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- User manual: 1 pc

Warranty:

- 2 Years for the Suction