

Suction Unit 7E-C (15 L/MIN)

Ref. 0803437

2 years
guarantee



Reliable and practical suction device, suitable for basic medical use and home care. It combines effective suction with low noise and easy transport. Ideal for daily support of patients with respiratory or other secretion removal needs.



Specifications:

- Suction 7E-C, 15Lt/min (Electric)
- Input power: 90VA
 - Max Negative Pressure: 85±10KPA
 - Max air flow: ≥20L/min (without jar)
 - Operation noise: Noise ≤65dB(A)
 - Liquid collection jar: 1000mL
 - Operating time: 30 minutes on, 30 minutes off

Dimensions:

- G.W :11kg
- N.W.:8kg
- Size: 52×46×34.5cm

Package Contents:

- Suction tube (length 2 m, Φ7×Φ12): 1 pc
- Suction catheter F8 (for children): 1 pc
- Suction catheter F12 (for adults): 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- Power cable: 1 pc
- User manual: 1 pc

Warranty:

- 2 Years for the Suction/
1 year for the Battery