

# Suction Unit 7E-H1 (28±4L/MIN)

Ref. 0803438

**2** years  
guarantee

○ Efficient and quiet suction device, suitable for medical use in clinical or home environments. It features a durable fluid collection container and offers stable operation with easy use and maintenance. Ideal for daily respiratory support needs.



## Specifications:

Suction 7E-H1 (Electric)

- Input power: 120VA
- Max Negative Pressure: 85±5kPa
- Max air flow: 28±4L/min (Device inlet)
- Operation noise: ≤65dB(A)
- Liquid collection jar: 1000mL, 1piece (PC material)
- Operating time: 30 minutes on, 30 minutes off

## Warranty:

- 2 Years for the Suction

## Package Contents:

- Suction tube (length 2 m, Φ7×Φ11): 1 pc
- Suction catheter F8 (for children): 1 pc
- Suction catheter F12 (for adults): 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- User manual: 1 pc

## Dimensions:

- G.W.: 7kg
- N.W.: 3.9kg
- Size: 48×16,5×28,5cm