

Suction Unit 7A-23A (32±4 L/MIN)

Ref. 0803464

2 years
guarantee



○ The electric suction device is ideal for medical use, offering reliable and efficient operation. Easy to use and clean, with low noise levels and consistent performance.



Specifications:

Suction 7A-23A, 32±4L (without jar)

- Input power: 180VA
- Max Negative Pressure: 90±10KPA
- Max air flow: 32±4L (without jar)
- Noise: ≤60dBA
- Operating time: 30 minutes on, 30 minutes off
- Liquid collection jar: 2000mL/pc*2

Warranty:

- 2 Years for the Suction

Package Contents:

- Suction tube (Φ7×Φ12, length 2 m): 1 pc
- Abdominal suction tube: 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- Power cord: 1 pc
- Foot switch: 1 pc
- User manual: 1 pc

Dimensions:

- G.W: 10.2kg
- N.W.: 8.5kg
- Size: 37.4×30.9×49.9cm