

○ Light weight, Flexibility and Freedom  
 ○ A private oxygen bodyguard for you and your family anytime, anywhere

# Portable Oxygen Concentrator

**Ref. 0803410**



Respond in 0.02s



Portable



## Specifications:

- Oxygen purity (%): 90-96
- Oxygen flow (Lit/min): Pulse flow / 1-4 Settings
- Maximum output pressure: 150 kPa
- Power supply AC: Input AC 100V-240V, 50Hz/60Hz, Output DC 12V, 10A
- DC Power supply: DC 12V-16V
- Noise dB(A): ≤ 55
- Battery Life: 3 hours and 20 minutes (at setting 2)
- Battery charging time: up to 5 hours from full discharge to full charge
- Weight 2.4kg (battery included)
- Dimension: 220mm x 220mm x 80mm
- Warranty: 3 years

## Package included:

- User Manual
- Battery 8 cells
- AC Power Supply
- DC Power Supply
- Carry bag
- Nasal cannula for Oxygen Therapy



Detachable molecular sieve  
 Quick and Convenient replacement



Pulse flow High efficient and Rapid  
 Oxygen delivery



Modular Lithium Batteries  
 Quick and Stable oxygen inhalation



Four-volume adjustment  
 Easy and Convenient control

