

# Commode Chair

Ref.: 0803390



## 2 years guarantee

### Features:

- Suitable for use in shower
- Tools Free, easy to assemble and disassemble
- Manufactured from plastic and Aluminum
- Its design enables the backrest and armrests to be removed according to user willing
- Extra cushions for the seat and backrest for better comfort
- Commode with capacity of 2.5 litres and lid
- Lightweight
- Easy, telescopically adjustable seat height
- Anti-slip rubber tips, on the legs



### Technical Specifications:

- Frame: Aluminum, PP (Polypropylene), Height Adjustable in 6 Levels,
- Color: Grey
- Seat/Backrest material: PE & PVC

Measurement in cm						
	46	43	54	42-54.7	5.1Kg	120Kg