

Suction Unit 7E-G1 (27±4L/MIN) (With Battery)

Ref. 0803439

2 years
guarantee



Portable suction device with rechargeable battery, ideal for use in vehicles or outdoors. Features an oil-free pump, overflow protection, and wall mounting capability (upon request). Ideal for emergencies and travel, with an optional disposable system for maximum hygiene.



Specifications:

Suction 7E-G1 (with Battery)

- DC power supply: AC 100–240V, 50/60Hz or DC 12V
- Input power: 150VA
- Max Negative Pressure: 85±5kPa
- Air flow rate: 27L±4/min (Device inlet)
- Operation noise: ≤65dB
- Operating time: 30 minutes on, 30 minutes off
- Running time of battery operation: Approx. 60 minutes

Warranty:

- 2 Years for the Suction/
1 year for the Battery

Package Contents:

- Suction tube (Φ7×Φ11, length 2 m): 1 pc
- Suction catheter F8 (children): 1 pc
- Suction catheter F12 (adults): 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- Power cable: 1 pc
- User manual: 1 pc

Dimensions:

- G.W.: 7 kg
- N.W.: 4 kg
- Size: 48×18,9×28,5cm

