

# Oxygen Concentrator

Ref. 0803302

3 years  
guarantee



Medical grade



High oxygen concentration



10L oxygen flow



Cumulative time function



Failure detection



Castors



High efficiency particulate air Filter

## Specifications:

- Flow range: 1-10L/min
- Oxygen concentration: 87% - 95.5%
- Maximum limited pressure: 100kPa
- Sound level: 55dB(A)
- Input power: 700VA
- Minimum operating time: 30minutes
- Dimension: 410mm x 375mm x 668mm
- Weight: 23.5kg
- Automatic reminder for filter cleaning every 100 hours
- Automatic reminder for filter replacing every 1800 hours

○ New 10-liter oxygen concentrator 0803302 is intended for use as an oxygen supplementation device in professional healthcare facilities and home environments. It provides a high concentration of oxygen to individuals requiring oxygen therapy.



## Package included:

- Power Supply Cable
- Humidifier
- Filter
- Fire safety valve
- Nasal cannula for Oxygen Therapy
- User Manual