

Auto-CPAP YH – 680A

Ref. 0803780

**Remote Control
Capability!**

Auto-CPAP YH – 680B

Ref. 0803777 **GSM Wifi Bluetooth**

GSM, Wifi και Bluetooth για
απομακρυσμένο έλεγχο, εξαγωγή
δεδομένων και αλλαγή ρυθμίσεων



Product Parameters

Model: YH-680A/ YH-680B

Color: Black

Mode: CPAP, Auto CPAP

Pressure: 4–20cm H₂O

Noise Level: 27dBA ± 2dBA

Dimensions (WxHxD) cm: 28.5 × 15.5 × 12.5

N.W.: 1.5Kg

○ The Auto-CPAP provides modern and flexible APAP therapy with precise event detection and a high level of comfort.

Technical Specifications:

- Auto on/off
- Auto Ramp
- Automatic humidifier temperature adjustment
- Mask fit check
- Exhalation relief FPS with 4 levels (0–3)
- Ramp time: 0–60 minutes
- Leak alarm
- Quiet operation: 27 dBA ±2 dBA
- Humidifier with 6 adjustable levels
- For detailed technical specifications, please refer to the user manual

Package Contents:

- Humidifier
- 22mm breathing tube
- Nasal mask
- Power adapter
- Power cord
- Set of 2 filters
- Carrying bag
- SD card
- User manual

Warranty:

- 2 years for the device
- 1 year for the power adapter
- 90 days for the humidifier