

# Auto-CPAP LC-APAP

Ref. 0805800

○ The Auto-CPAP LC-APAP combines advanced technology with ergonomic design, offering effective and comfortable therapy for sleep apnea.



## Technical Specifications:

- Auto on/off function
- Auto Ramp
- Ramp time: 0-60 min
- Expiratory pressure relief FPS with 4 levels (0-3)
- Leak alarm
- Quiet operation at 29 dBA
- Humidifier with 5 adjustable levels
- For detailed technical specifications, refer to the user manual

## Package Contents:

- 22mm breathing tube
- Power supply unit
- Power cable
- 2 filters
- Carrying bag
- SD card
- User manual

## Warranty:

- 2 years on the device
- 1 year on the power supply unit

## Product Parameters

Model: LC-APAP

Color: White

Mode: CPAP, Auto CPAP

Pressure: 4-20cm H<sub>2</sub>O

Noise Level: 29dBA ± 2dBA

Dimensions (W×H×L): 14 × 15,6 × 9,4cm

Net Weight: 1Kg