



The main advantage of scooters like "APOLLO" is that they allow their users to regain their independence. People who previously relied on others for transportation can now move around with comfort and autonomy. Regardless of their destination, they can go to the doctor, visit friends, go shopping and enjoy life without restrictions.



Specifications

- Colors Available: Black-Gold
- Motor Output W: 270W,
- Motor Rated Power W: 700W
- Safe Climbing Angle: 6 degree
- Handle Bar Type Or Joystick: Delta handle bar
- Electromagnetic Brake
- Emergency Hand Brake
- Seat Height (Adjustable)
- Seat Belt
- Adjustable Seat (Forwards and Backwards)
- Light Reflectors Front And Rear
- Basket Included
- Horn
- Product Frame: Steel
- Battery Warranty: 6 Months



Technical data

46cm	44cm	43cm	40cm	49cm	103cm	86cm	57cm	8" 8"
AGM 24V 12 Ah	6°	8kph	15km	120kg	44kg	270W x 24V 2 Pieces	24V x 5AMP Warranty 1 Έτος	

Operating console



Swivel Seat

