



Mobility scooters like **"ARES"** offer an essential blend of independence, convenience, and comfort for individuals with limited mobility. Designed to empower users to move freely, they open up possibilities for navigating the world with confidence and ease.



Specifications

- Colors Available: Silver
- Motor Output W: 340W,
- Motor Rated Power W: 700W
- Turning Radius: 1400mm
- Safe Climbing Angle: 9 degree
- Handle Bar Type Or Joystick: Delta handle bar
- Electromagnetic Brake
- Emergency Hand Brake
- Headrest (Adjustable)
- Seat Height (Adjustable)
- Seat Belt
- Adjustable Seat (Forwards and Backwards)
- Turning Lights & Signal
- Reverse Signal
- Light Reflectors Front And Rear
- Alarm Signal (When Both Trurning Lights Flash Simultaneously)
- Basket Included
- Mirrors Included
- Horn
- Product Frame: Steel
- Battery Warranty: 6 Months



Technical data

46cm	44cm	43cm	40cm	58,3cm	125,5cm	96,5cm	63cm	10" 10"
AGM 24V 52Ah	9°	12kph	50km	150kg	83kg	340W	24V x 6AMP	

