

BILEVEL YH-820 AUTO

Ref. 0803778

BILEVEL YH-825 AUTO

Ref. 0803779

**Remote Control
 Capability!**

BILEVEL YH-820 AUTO

Ref. 0803815 GSM

BILEVEL YH-825 AUTO

Ref. 0803816 GSM

*GSM and Wifi for remote control,
 data export, and setting adjustments*



YH - 820 AUTO



YH - 825 AUTO

The Bilevel Auto-S offers modern and flexible therapy with precise event detection and a high level of comfort.

Technical Features:

- Smart START/STOP
- Central Sleep Apnea (CSA) detection
- Leak alarm
- Quiet operation at 28dBA ±2 dBA
- Humidifier with 6-level adjustment
- For detailed technical specifications, refer to the user manual

Package Contents:

- AC Power Supply
- Humidifier
- Tube
- SD Memory Card
- Air Filters (2 pcs)
- User Manual
- Carrying Bag
- Oral-Nasal Mask
- O₂ Enrichment Connector

Warranty:

- 2 years on the device
- 1 year on the power supply
- 90 days on the humidifier

Product Parameters

Model: YH-820/ YH-825

Color: Black/White

Mode: CPAP, S, Auto S

Pressure: 4–20cm H₂O

Noise Level: 28dBA±2dBA

G.W./N.W.: 1.5Kg/1.6Kg

Dimensions (L×W×H) cm: 28,5 × 15,5 × 12,5