

Hospital Mattress

"PRO I"

Ref. 0805941

Dimensions: 191 × 82 × 14cm

"PRO II"

Ref. 0805942

Dimensions: 195 × 90x 14cm




The hospital-type mattress is designed to prevent, limit, and delay the formation of pressure ulcers in bedridden patients or individuals who remain in bed for long periods of time. It can be used in hospitals or at home.

LEVEL 1: LYING SURFACE

2 specially designed anatomical zones made of 100% PU (polyurethane) foam, density 45Kg/m³

LEVEL 2: MATTRESS BASE

100% PU (polyurethane) foam, density 45Kg/m³



Cover with zipper



Mattress of 2 layers (internally)

- Provides uniform pressure distribution on the body to avoid the formation of edema and irritation
- Suitable for patients with low to medium risk of developing bedsores (Stage I to II), provided that the patient's position is frequently changed by healthcare personnel
- Made of 100% PU foam (polyurethane), density 45Kg/m³ (internal)
- With a cover of 100% PU (polyurethane)
- With anatomical flexible zones that reduce the application of forces to the patient's abdomen
- Maximum load: 250Kg
- Softness for maximum protection of soft tissues
- High comfort and body support
- Maintenance of body temperature through special notches - channels - breathable material properties
- Ergonomic design - can also be used in beds with multi-section function following the movements of the bed
- With a waterproof & anti-slip cover, with taped seams & a cover on the zipper to prevent liquids from entering the interior
- It has a 180° zipper, for easy removal of the cover
- It has a plastic zipper and handles around the perimeter of the cover for easy transfer of the patient to the magnetic resonance imaging (MRI) scanner
- Antifungal, Antimicrobial & Antibacterial
- Cover with flame retardant properties
- The cover is removable, washable in the machine at 80°C, with the possibility of easy disinfection

*Ref number 0805942 is suitable for use with bed 0805940 HOSPITAL BED "PRO"

**Ref number 0805941 is suitable for use with bed 0805905 HOSPITAL BED

3 years
 guarantee