

COMMODE CHAIR "SEAL"

Ref. 0815017

The "SEAL" commode chair is intended for use as a toilet aid by people with limited mobility, such as the elderly, patients, or people with disabilities. It is designed for indoor use and offers safety and comfort in daily hygiene.



2 years
 guarantee

Specifications:

- Steel frame
- Suitable for all types of toilets with removable bucket
- Adjustable height in 6 levels by push button system
- Seat and backrest material: Vinyl with flame retardant properties

DIMENSIONS AND WEIGHT							
Measurement in cm							
0815017	43	46	34.5	58-62	48-63	8Kg	130Kg