



# NEW PRODUCTS

RESPIRATORY

ORTHOPEDICS

FIRST AID CARE







## CONTACT DETAILS

### Headquarters:

Kathiana Akrotiri, Chania, Crete P.C. 731 00

T.: +30 28210 63222

E: [medicalexports@mobiak.com](mailto:medicalexports@mobiak.com)

[www.mobiakcare.com](http://www.mobiakcare.com)

### Retail Shop, PYROMEDICAL Ltd:

Markou Botsari 96-98, Chania, Crete 731 36

### Design & Installation Shop

#### PYROPROSTASIA LTD:

Anapafseos 20, Chania, Crete 731 00

## DISTRIBUTION CENTERS

### Central, South Greece & Ins. Greece:

Location: Schineza, Aspropyrgos, Attica P.C. 193 00

T: +30 210 4614006, 210 5593004

### North Greece & Balkan Area:

Arkadiou 8, Kalochori, Thessaloniki P.C. 570 09

T: +30 2310 785668, 2310 778477

### Central & Eastern Crete:

P Road, Industrial Zone of Heraklion, Crete P.C. 716 01

T: +30 2810 333151

### Chania & Western Crete:

Kathiana Akrotiri, Chania, Crete P.C. 731 00

T: +30 28210 63222



## UBREATH BA200 Breath Gas Test System

Ref. **0803881**



UBREATH® Breath Gas Test System (Model BA200) quantitative measures Nitric Oxide (NO) and Carbon Monoxide (CO) in human breath, as an adjunct to establish clinical and laboratory assessments of inflammatory processes such as asthma, chronic respiratory disease, or as monitored treatment in smoking cessation clinics

### Specifications

- Screen: 10.1-inch IPS touchscreen, (1200\*1920)
- Memory: 16GB
- Operating principle: Electrochemical sensor
- Methods of use: mouthpiece or nasal
- Measurements: FeNO, CaNO, FnNO, sNO
- Measuring range: 0-4000ppb
- Exhalation time: ~10 sec. in adults, ~6 sec. in children
- Measurement appearance time: <30sec
- Accuracy: ±2.5ppb in measurement ≤50ppb, <3% in measurement >50ppb
- Repeatability: ≤1.5 ppb at measurement ≤50ppb, <3% at measurement >50ppb
- Sensor sensitivity: 0.1ppb
- Lower limit: 1ppb
- Dimensions: (259 x 115 x 200)mm
- Power supply: AC (100-240V) ~50/60Hz
- Battery: Lithium 3.7V 10000mAh, with autonomy about 6 hours
- Operating environment: (5°C - 35°C), 80% RH
- Warranty: 2 years the device, 6 months on battery & accessories



### Package contents

- Main unit
- Breath controller (change every 1500 measurements) code: 0803883
- USB cable TYPE-C
- Battery
- User manual
- Screwdriver
- OTG cable



### Accessories



Charger

**To use the device, you must have purchased the code: 0803882**

(NO sensor package + 100 patient filters /mouthpieces)



## Cough Assist

Ref. **0803284**



**PEGASO PLUS COMBI:** Cough assistant and percussor with manual, automatic and auto - adaptive modes

### FIELD OF APPLICATION

Hospitals  
Home Care  
Respiratory Therapy Clinics

### SUITABLE FOR:

Muscular dystrophy  
Spinal muscular atrophy  
Cystic fibrosis (Mucoviscidosis)  
Myasthenia gravis  
Spinal cord injury

Poliomyelitis  
Amyotrophic lateral sclerosis  
Chronic obstructive pulmonary disease (COPD)  
All conditions that includes muscular weakness

### OUTSTANDING FEATURES

Internal battery, which allows an individual therapy up to 3 days  
Portable size and lightweight (3.2 kg)  
Oscillation function during cough assistant therapy  
Inspiratory volume & Expiratory peak flow visualization  
High monitoring of respiratory parameters thanks to SpO2 oximeter  
Fast switch between Percussor & Cough assistant therapy

### TECHNICAL FEATURES

SETTINGS		
COUGH ASSISTANT		ACCURACY
Positive pressure I	From 0 to +70hPa	±10% or ±2hPa
Negative pressure E	From 0 to -70hPa	±10% or ±2hPa
Inspiration time	From 0,5 to 9,9 sec	±10% or ±2sec
Expiration time	From 0,5 to 9,9 sec	±10% or ±2sec
Pause time	From 0 to 9,9 sec	±10% or ±2sec
Trigger	Levels from 1 (most sensitive) to 9 (less sensitive)	
Inhalation Flow (Rise time)	Low-Medium-High (1,2,3) Maximum inhalation flow at level 3: 240lmp The flow depends on the generated pressure and the breathing system resistance.	
Oscillations	OFF, only on I, only on E, I and E	
Oscillations Pressure	From 1 to 15hPa	±10% or ±2hPa
Oscillations Frequency	From 1 to 20Hz	±1hz
Exhalation Flow	Maximum exhalation Flow at -70hPa: -240lpm	
Inhaled Volume reading	0cc - 2500cc (±100cc or ±15%)	
PERCUSSOR		
Inspiratory positive pressure	From 0 to +70hPa	±10% or ±2hPa
Frequency	From 50 to 600bpm	±10%
I/E Ratio	Variable from 5.0:1 to 1:5.0	±10%
OXIMETER READING		
SpO <sub>2</sub>	70%-100% accuracy= ±2%	50%-69% accuracy= ±3%
Pulse	25 to 240BPM accuracy= ±2 BPM	
OTHER SPECIFICATIONS		
Alarms	Power failure alarm, high pressures alarms, high temperature alarm, mechanical failure alarm, SpO <sub>2</sub> , Pulses	
Size	(27×26×13)cm (WxDxH)	
Weight	2.3kg	
Power supply	100/240V 50/60hz, 120VA	
Operating conditions	Temperature of functioning: from 10-40°C, Relative humidity: from 10-90%	
Safety standards	EN60601-1, EN60601-1-2	
Insulation class	Class II BF	
Risk class	2B	

### Accessories

- 0803288 Adult mask  
Replacement every 3 months
- 0803289 Paediatric mask  
Replacement every 3 months
- 0803290 Catheter mount  
Replacement every 3 months
- 0803291 Percussor expiratory adaptor  
Replacement for every new patient and every 3 months
- 0803292 Antidust filter  
(Package of 5 pieces)  
Replacement at least one a month
- 0803293 Tube 1.8m - Replacement  
for every new patient and every 3 months
- 0803294 Antibacterial filter  
Replacement for every new patient and at least once a week
- 0803295 Connector - Replacement  
for every new patient and every 3 months
- 0803296 Mouthpiece - Replacement  
for every new patient and every 3 months

## FULL FACE MASK YF-03 YUWELL

- Ref: **0803442** Small
- Ref: **0803443** Medium
- Ref: **0803444** Large



## NASAL MASK YN-01 YUWELL

- Ref: **0803445** Small
- Ref: **0803446** Medium
- Ref: **0803447** Large



### Description

#### Specifications

- Weight: 103gr
- Therapy Pressure Range: 4-30 cmH<sub>2</sub>O
- Tube Connection: φ22mm

#### Silicone Cushion

The mask is made of rebound silicone cushion which is airtight, light, and soft. It is designed to reduce facial pressure, improve fit and comfort.

#### Adaptive Silicone Cushion

The silicone cushion has an adaptive design at the nose bridge, which can be automatically adjusted according to different face shape, better relieve the pressure of the high nose and improve user experience.

#### Frame

Lightweight and Strong, Easily Adjustable & Better Hygiene

#### Swiveling Bracket

360° Swiveling Bracket. Ensures Greater Comfort of Movement during sleep.

#### Quiet Vent-hole

Honeycomb design, reducing airflow erosion intensity, gentle and quiet.

#### Soft And Light Headgear

Selected soft skin-friendly fabrics, light and stable, help patients adapt to treatment comfortably.

#### Specifications

- Weight: 63gr
- Therapy Pressure Range: 4-30 cmH<sub>2</sub>O
- Tube Connection: φ22mm

## Auto-CPAP YH-550 GSM

Ref. **0803432**

## Auto-CPAP YH-550 Wi-Fi

Ref. **0803433**

### Remote control capability !



The Auto-CPAP offers modern and flexible APAP therapy with precise event detection and a high level of comfort. Innovative design with remote control capability and an integrated humidifier.

### Technical Specifications

- GSM and Wi-Fi for remote control, data export, and setting adjustments
- Smart START/STOP
- FPS Exhalation Assist, 4 Levels 0-3
- Ramp 0-45min
- Leak alarm
- Quiet operation 28dBA  $\pm$ 2 dBA
- Humidifier with 6 level adjustment
- Detailed technical characteristics refer to the user manual

### Package Contents

- Humidifier
- Breathing tube 22mm
- Nasal Mask
- Power supply
- Power supply cable
- Set, 2 filters
- Carrying Bag
- SD card
- Instructions for use

### Warranty

- 2 years on the device
- 1 year on the power supply
- 90 days in the Humidifier

### Product Parameters

Model	YH-550
Color	White
Mode	CPAP, Auto CPAP
Pressure	4-20cm H <sub>2</sub> O
Noise	28dBA $\pm$ 2dBA
G.W./N.W.	2.9/1.5kgs
Packing	41.5 $\times$ 30 $\times$ 16cm

### Spare parts/ consumables

- Ref.: 0806206 Tube CPAP-BILEVEL
- Ref.: 0803371 Filter YH-550 YUWELL
- Ref.: 0803373 Humidifier YH-550
- Ref.: 0803372 Power supply
- Ref.: 0811590 Power Bank

## BILEVEL YH AUTO

Ref. **0803778** YH - 820

Ref. **0803779** YH - 820 (GSM)

Ref. **0803815** YH - 825

Ref. **0803816** YH - 825 (GSM)



YH - 820



YH - 825 / YH - 830

### Technical Features

- Smart START/STOP
- Central Sleep Apnea (CSA) detection
- Leak alarm
- Quiet operation at 28dBA  $\pm$ 2 dBA
- Humidifier with 6-level adjustment
- For detailed technical specifications, refer to the user manual

### Package Contents

- AC Power Supply
- Humidifier
- Tube
- SD Memory Card
- Air Filters (2 pcs)
- User Manual
- Carrying Bag
- Oral-Nasal Mask
- O<sub>2</sub> Enrichment Connector

### Warranty

- 2 years on the device
- 1 year on the power supply
- 90 days on the humidifier

### Product Parameters

Model	YH-820	YH-825
Color	Black	White
Mode	CPAP, S, Auto S	CPAP, S, Auto S
Pressure	4-20cm H <sub>2</sub> O	4-20cm H <sub>2</sub> O
Noise Level	28dBA $\pm$ 2 dBA	28dBA $\pm$ 2 dBA
Dimensions (L×W×H) cm	28,5 × 15,5 × 12,5	28,5 × 15,5 × 12,5
Net Weight	1,5 Kg	1,6 Kg

## Auto-CPAP YH - 680A

Ref. **0803780**

## Auto-CPAP YH - 680B

Ref. **0803777**

### **Remote control capability!**

GSM for remote control, data export, and setting adjustments



The Auto-CPAP provides modern and flexible APAP therapy with precise event detection and a high level of comfort.

### Technical Specifications

- Touch Screen
- With remote control capability via Wi-Fi, GSM (SIM), Bluetooth (for ref.: 0803777)
- Auto on/off
- Auto Ramp
- Automatic humidifier temperature adjustment
- Mask fit check
- Exhalation relief FPS with 4 levels (0–3)
- Ramp time: 0–60 minutes
- Leak alarm
- Quiet operation: 27 dBA ±2 dBA
- Humidifier with 6 adjustable levels
- For detailed technical specifications, please refer to the user manual

### Package Contents

- Humidifier
- 22mm breathing tube
- Nasal mask
- Power adapter
- Power cord
- Set of 2 filters
- Carrying bag
- SD card
- User manual

### Warranty

- 2 years for the device
- 1 year for the power adapter
- 90 days for the humidifier

### Product Parameters

Model	YH-680
Color	Black
Mode	CPAP, Auto CPAP
Pressure	4-20cm H <sub>2</sub> O
Noise Level	27dBA ± 2 dBA
Dimensions (W x H x D)	285 x 155 x 125mm
N.W.	1.5 Kg

## AUTO-CPAP LC-APAP

Ref. **0805800**



The Auto-CPAP LC-APAP combines advanced technology with ergonomic design, offering effective and comfortable therapy for sleep apnea.

### Technical Specifications

- Auto on/off function
- Auto Ramp
- Ramp time: 0–60 min
- Expiratory pressure relief FPS with 4 levels (0–3)
- Leak alarm
- Quiet operation at 29 dBA
- Humidifier with 5 adjustable levels
- For detailed technical specifications, refer to the user manual

### Package Contents

- 22mm breathing tube
- Power supply unit
- Power cable
- 2 filters
- Carrying bag
- SD card
- User manual

### Warranty

- 2 years on the device
- 1 year on the power supply unit



### Product Parameters

Model	LC-APAP
Color	Άσπρο
Mode	CPAP, Auto CPAP
Pressure	4-20cm H <sub>2</sub> O
Noise Level	29dBA ± 2 dBA
Dimensions (W×H×L)	14 × 15.6 × 9,4 cm
Net Weight	1 Kg

### Accessories/Spare Parts

- Code: 0806206 – LC-APAP TUBE
- Code: 0805809 – LC-APAP FILTER
- Code: 0805816 – LC-APAP Humidifier
- Code: 0805815 – LC-APAP Power Supply

## Portable Oxygen Concentrator

Ref: **0803410**



**Light weight, Flexibility and Freedom.**

**A private oxygen support for you anytime, anywhere**

**Trustworthy stable oxygen supply with 4 Core Technologies**

- Detachable molecular sieve. Quick and Convenient replacement
- Pulse flow High efficient and Rapid Oxygen delivery
- Modular Lithium Batteries
- Four-volume adjustment. Easy and Convenient control
- Battery Warranty: 1 year

### Technical Specifications

Oxygen Purity (%)	90-96
Oxygen Flow (Lit/min)	Pulse-flow / 1-4 Settings
Max Outlet Pressure	150 kPa
AC Power Supply	Input AC 100V-240V, 50Hz/60Hz, Output DC 12V,10A
DC Power Supply	DC 12V-16V
Noise Level dB(A)	≤ 55
Dimensions (cm)	L: 22 x H: 22 x W: 8
Net weight (Kg)	2.4

Battery Duration		Breath Rate, Breath/Minute (BPM)	Pulse Volume (ml) per Flow Setting			
			1	2	3	4
Setting	Hours	15	14,0	28,0	42,0	50,0
		20	10,5	21,0	31,5	37,5
1	5	25	8,4	16,8	25,2	30,0
		30	7,0	14,0	21,0	25,0
2	3:20	35	6,0	12,0	18,0	21,4
		40	5,3	10,5	15,8	18,8
3	2:20					
4	2					
		Total Volume per Minute (ml/min)	210	420	630	750

### Package Contents

- User Manual
- Battery 8 cells
- AC Power Supply
- DC Power Supply
- Carry bag
- Nasal cannula for Oxygen Therapy

## Oxygen Concentrator 10Lt

Ref: **0803302**



- High oxygen concentration
- Cumulative time function
- Failure detection
- High efficiency particulate air Filter

### Technical Specifications

Oxygen Purity (%)	87-96
Oxygen Flow (Lit/min)	0-10
Max Outlet Pressure	0,5
Power Supply (VAC)	230V±10%, 50Hz±1Hz
Power (W)	700
Noise Level dB(A)	≤55
Dimensions (cm)	L: 41 x W: 37.5 x H: 66.8
Net weight (Kg)	23,5

- Automatic reminder for filter cleaning every 100 hours
- Automatic reminder for filter replacing every 1800 hours

### Package Contents

- Humidifier
- Fire safety valve
- Nasal cannula for Oxygen Therapy
- Filter
- Power Supply Cable
- User Manual

## 9F-5B Oxygen Concentration 5Lt

Ref: **0803430**

The 9F-5B oxygen concentrator provides quiet operation, ensuring that you can relax without annoying noises.

It features advanced fault detection to alert users to any operational problems.

Designed with a sleek, modern look that fits into any home décor.



### Technical Specifications

Oxygen Purity (%)	87-96
Oxygen Flow (Lit/min)	0-5
Max Outlet Pressure (bar)	0,5
Power Supply (VAC)	230V±10%, 50Hz±1Hz
Power (W)	350
Noise Level dB(A)	43
Dimensions (cm)	40×30×55
Net weight (Kg)	17,5

### Package Contents

- Cannula
- Humidifier
- Humidifier tube
- 2x inlet filters
- Dast filter
- 2x fuses
- Firetrap
- User Manual

## Compressor Nebulizer

Ref: **0811850**



Your Nebulizer is intended for use in the treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy.

- Suitable for hospital and home use
- Ideal for all ages
- Modern and Ergonomic Design with carrying handle
- Compact & lightweight design

### Technical Specifications

Voltage	220 V AC / 50 Hz
Power	168 Watt
Medication Capacity	8 ml
Particle Size	0.5-10 $\mu$ m
Sound Level	<65 dB(A) (1 meter away)
Average Nebulization Rate	>0.25 mL/min (0.9% saline water)
Compressor Pressure Range	2.5 $\pm$ 0.5 bar
Operating Pressure Range	0.65-1.31 bar
Liter Flow Range	8.5 $\pm$ 0.5 lpm
Operating Temperature	5 °C – 40 °C
Operating Humidity Range	15-80% (RH)
Storage Temperature Range	-25 °C – 60 °C
Dimensions	18x15x8 cm
Weight	1.25 Kg
Warranty	3 years

### Package Contents

- Air tube, Nebulizer Tank, Mouthpiece, Angular connector for mask use in supine position, Adult Mask, Pediatric Mask, 5 pcs extra filters, User Manual

## Mesh Nebulizer with Batteries

Ref: **0803431**



The nebulizer 0803431 is a nebulization therapy device that works with new generation ultrasound technology and is designed for inhalation of almost all kinds of drugs.

It is the ideal choice for adults and children, both for home and emergency treatment, travel or outdoor use.

### Technical Specifications

- Input power: DC 3.0V
- (two "AA" 1.5V alkaline batteries)
- Operation noise: approx. 50dB(A)
- MMD: 3.49m
- Residual Volume: 0.25ml
- Weight: 108g (without batteries)
- Operation: Continuous
- Battery duration: 60 min
- Dimension: 67 x 48 x 125mm
- Warranty: 2 years

### Package Contents

- Medication Chamber
- Pediatric Mask
- Adult Mask
- Mouthpiece
- Batteries x 2
- User Manual
- Carrying Case



## Mesh Nebulizer CONTEC - Rechargeable

Ref: 0803845



It is powered by a 5 V power adapter or a 3.7 V lithium battery, with the advantages of compact and portable design, easy to use, adjustable (two grades) and high in nebulization rate, which is applicable for use in hospitals, clinic and family.

### Specifications

#### Performance

- Median particle diameter (D 0.50): 1–4 $\mu$ m ( $\pm$ 25%)
- Equivalent volume particle diameter distribution: the proportion of particle diameter below 5  $\mu$ m is not less than 60%
- Input power: < 5 VA

#### Nebulization rate:

- Grade I: 0.25 mL/min.
- Grade II:  $\geq$  0.15 mL/min
- Cup capacity: 10 mL
- Liquid volume: 2 - 10 mL
- Residual volume:  $\leq$  0.3 mL
- Noise:  $\leq$ 50 dB (A)
- Working frequency: 110 KHz  $\pm$  10%
- Power supply: DC 5V or 3.7 V rechargeable lithium battery
- Safety class: class II equipment, type BF applied part
- The degree of waterproof and dustproof: IP22
- Mode of operation: continuous operation

#### Operating environment:

- Temperature: 5 °C ~ 40°C
- Relative humidity: 15% - 90%
- Atmospheric pressure: 700 hPa ~ 1060 hPa

#### Package included:

- Usb cable
- Adult mask
- Child mask
- User manual

**Dimension:** 67.42 mm(L)  $\times$  48.53 mm(W)  $\times$  116.92 mm(H)

**Weight:** 100g

# Oxygen Therapy Consumables



## Multi-vent Mask Adjustable Venturi, elongated, with seven connectors

Ref: **0803202**

**Adult**

Ref: **0803203**

**Pediatric**



- Tube length 2.1m
- Single venturi valve is available with fixed concentration
- Choice of 7 different valves, with different concentrations: 24% (blue), 28% (yellow), 31% (white), 35% (green), 40% (pink), 50% (orange), 60% (red)

## Non-Rebreath Oxygen Mask with tubing 2,10 m

Ref: **0803204**

**Adult**

Ref: **0803205**

**Pediatric**



- Tube length 2.1m
- Supply of high concentration oxygen
- Offered with head strap and adjustable nose clip



## Oximeter

Ref: **0803440**



**Finger pulse oximeter measuring Oxygen (SpO<sub>2</sub>) and heart rate.**

- Automatic shutdown
- Pulse column display

### Technical Specifications

- Measurement range of pulse rate: 25bpm-250bpm
- Accuracy of SpO<sub>2</sub>: 70% ~100%: ±2%
- Accuracy of pulse rate: ±1%/ ±1bpm
- Power consumption: 40mA
- Operation environment: Temperature(°C): 5~40,
- Relative humidity: 80%
- Boot method: Push button to boot
- Power supply: 2 × 1.5V AAA alkaline batteries
- Display: LED
- Dimension: 58 × 33 × 33.5mm
- Weight: 35g (without batteries)
- Automatic shutdown time: 8s
- Warranty: 2 Years

### Package Contents

- Batteries x 2
- User Manual



## Infrared Forehead Thermometer YT-1

Ref: **0803309**

The device is used to measure the forehead temperature of the human body.

### Package Included:

- Thermometer case
- Batteries



## Portable Vein Finder

Ref: **0815740**

**2** Years  
Warranty



The portable vein finder provides a timely and accurate image of blood vessels on the skin. It has an ergonomic design and advanced infrared imaging technology to provide an accurate digital image of the veins on the patient's skin in real time. It also helps doctors and nurses easily locate veins.

### Indications:

Must be used exclusively by specialized and trained personnel. It is adapted to different colors for different ages, body types, skin colors, weights and various operating environments. It is suitable for adults, children and infants. In addition, it allows the observation of the thickening, length, curvature and branch information of the veins.

### Technical Specifications:

- Projection Colors: 12
- Light: Infrared Light
- Projection Size: 3
- Brightness Levels: 6
- Enhancement Mode: Yes
- Inversion Mode: Yes
- Sleep Mode: Yes
- Battery Capacity Monitoring: Yes
- Infrared Wave: 760-940nm
- Detection Depth: 0-10 mm
- Detection Distance: 18-30 cm
- Accuracy Error:  $\leq 0.3$  mm
- Noise Level:  $\leq 35$  dB
- Rechargeable Lithium Battery: 3400 mAh
- Continuous Operation Time: 3 Ωπεç
- Power Supply : 100V-240V - 50/60Hz AC
- Weight: 280 gr
- Product Dimensions (L x W x H): 20 x 6.2 x 5.5 cm

### Package Included:

- Vein Finder Device
- Tourniquet
- Charger Cable
- Calibration Card
- Instructions for Use



## Automatic Blood-Pressure monitor

Ref: **0803441**



**Digital blood-pressure monitor that measures systolic pressure (SYS), diastolic pressure (DIA) and heart rate.**



### Technical Specifications

- 80 groups of memory
- Average of three measurements
- Systolic/Diastolic and Pulse indication
- Measuring range:  
Diastolic: 20-210 mmHg,  
Cuff pressure: 0-300mmHg,  
Systolic: 40-260mmHg,  
Pulse: 40-200 bpm
- Dimension: 130 x 96 x 75mm
- Warranty: 2 Years
- Weight: 251g (without batteries)
- Power supply: 220-240V or 4 AA 1.5V batteries

### Package Contents

- Cuff
- USB cable
- User Manual

## PULSE OXIMETER

ref. **0803847**



The Pulse Oximeter 0803847 adopts Photoelectric Oxyhemoglobin inspection Technology integrated with Capacity Pulse Scanning & Recording Technology, which can be used to measure human oxygen saturation and pulse rate through finger. The device is suitable for use in the user's home environment, hospitals, and healthcare facilities.

### Features

- Display of SpO2 value
- Display of heart rate value and bar graph
- Display of pulse waveform
- Display of PI value
- Language selection: English
- Low battery indicator
- Alarms: finger detection failure, limit exceedance, low battery level
- Sound volume adjustment
- Screen brightness adjustment
- Memory for storing SpO2 & PR
- Data export capability to a computer
- For detailed technical specifications, refer to the user manual.

### Package Contents

- USB Cable
- SpO2 Oximeter Sensor
- Charger
- Screwdriver
- User Manual

**Dimension: 159 mm(L) x 72 mm(W) x 31mm(H)**

**Weight: about 250 g (including battery)**

### Accessories/Spare Parts

- Code:0803903 Oximetry Sensor

## AMBULATORY 24 HOURS BLOOD PRESSURE MONITOR

ref: **0803846**



### Operational environment:

- Temperature: 5°C~40°C
- Relative humidity: 15%~85%
- Atmospheric pressure: 700hPa ~ 1060hPa

Specific EMC, climate, mechanical environment: Do not use a mobile phone in the vicinity of the device, for the strong radiation field produced by mobile phone would interfere with the normal function of use.

### Storage environment:

- Temperature: -20°C~+55°C
- Relative humidity: ≤95%

Specific EMC, climate, mechanical environment: The device after packaged should be stored in a environment of -20°C~+55°C, relative humidity no more than 95%, and well-ventilated room with no corrosive gas. Strong shock, vibration and snow and rain should be avoided during transportation.

### Package included

- Cuff for adults
- BP extending tube
- Usb cable
- Case
- 2x screws
- User manual

**Dimension: 128mm(L) × 69mm (W) x 36mm(H) (without packaging)**

**Weight: <300 g (with battery)**

### Accessories / Spare Parts

- Code: 0803908 Cuff for adults
- Code: 0803909 BP extending tube

The ABPM 0803846 is a handheld ambulatory blood pressure monitor, which is designed according to oscillography theory. The device could monitor human body blood pressure for up to 24 hours continuously and dynamically, providing an accurate basis for the diagnosis. It is applicable for using in hospital, clinics and other medical institutions.

### Features

#### Performance:

##### NIBP:

- Measure Method: Oscillometry
- Measure Mode: The upper arm measure
- Automatic Measure Interval: 5, 10, 15, 20, 30, 40, 60, 90, 120, 180, 240 minutes
- Measure range: Pressure: 0kPa (0mmHg) -39.6kPa (297mmHg)
- Resolution: 1mmHg
- Accuracy: ±3mmHg
- Prompt parameter: SYS, DIA
- Inflation: automatic inflation by force pump
- Deflation: automatic deflation

##### PR:

- Measure range: 40bpm~240bpm
- Resolution: 1bpm

##### Safety:

Power supply: DC 3V (2 × 1.5V AA alkaline dry battery)

Safety type: internally powered device, type BF applied part with defibrillation protection

##### Product features:

- Compact and portable, user-friendly interface, easy to use.
- Patient range: adult, pediatric, neonate.
- 24-hour ambulatory NIBP monitoring function, up to 350 groups of ambulatory NIBP data can be recorded for once.
- Perfect combination of automatic and manual measurement method, up to 300 groups of data can be recorded for once by manual measure.
- High-definition color TFT display, strong visibility.
- By data review interface such as "data list", "trend graph, big font", NIBP data is clear at a glance.
- Display of low power prompt, error message and time.
- Supply two kinds of unit: mmHg/kPa.
- English menu.
- Over-limit prompt function for parameters.
- Communicate with PC, PC software can achieve data review, measured results analysis, view of trend graph, reports printing and other functions.

##### Software features:

- Connect to the device by USB interface.
- Download NIBP measure results from the device.
- Display of scoop-shape trend graph, filling-type trend graph, histo-gram, pie chart, correlation line graph.
- Edit every piece of NIBP data and add annotation to it.
- Edit basic information, doctor's advice, NIBP status instruction, current medicine-taken information, etc.
- Support report printing and print preview.

## DIGITAL BLOOD PRESSURE MONITOR WITH OXIMETRY

ref: **0803871**

**New Product**



The blood pressure monitor is used for non-invasive measurement of arterial blood pressure and SpO<sub>2</sub> oximetry in humans (adults, children, neonates). Each measurement includes detailed timing, systolic pressure, diastolic pressure, mean pressure, and pulse rate.

**2** / **4** Years  
Warranty

### Features

- 2.8" Color LCD screen
- Cuff pressure range: 0-297 mmHg
- Measurement range:
  - Adults: SYS: 30-270 mmHg, DIA: 10-220 mmHg
  - Pediatric: SYS: 30-235 mmHg, DIA: 10-195 mmHg
  - Neonatal: SYS: 30-135 mmHg, DIA: 10-100 mmHg
- Pulse measurement range: 40-240 bpm
- Overpressure protection:
  - Adults: 295±5 mmHg
  - Pediatric: 240±5 mmHg
  - Neonatal: 145±5 mmHg
- Measurement analysis: Pressure: 1 mmHg, Pulse: 1 bpm
- Measurement accuracy: Pressure: ±3 mmHg, Pulse: ±5 bpm
- Ability to adjust blood pressure and oximetry limits
- Ability to activate an audible alarm when limits are exceeded
- Operating temperature: 5-40°C
- Operating atmospheric pressure: 700 hPa – 1060 hPa
- Storage environment: Temperature: -20 to +55°C, Humidity ≤ 95%
- Power supply: 4 alkaline batteries AA (110-240V transformer optional)
- Power consumption: ≤ 600 mA
- IP22
- For detailed technical specifications, refer to the user manual.

### Package Contents:

- Adult cuff
- 4 x Alkaline batteries AA
- USB cable for data transfer
- Transportation case
- User manual
- 0803872 Adult SpO<sub>2</sub> sensor

**Dimensions: 130 mm (W) X 110 mm (L) X 80 mm (H)**

**Weight: approximately 300 g (without batteries)**

### Accessories

- 0803872 Adult SpO<sub>2</sub> sensor
- 0803873 Children's SpO<sub>2</sub> sensor
- 0803874 Neonatal SpO<sub>2</sub> sensor
- 0803875 Adult cuff (22-43 cm)
- 0803876 Neonatal cuff (6-11 cm)
- 0803877 Infant cuff (10-19 cm)
- 0803878 Child's cuff (18-26 cm)
- 0803879 Power supply



## Electrocardiograph 12 Channel with Printer

ref. **0803840**



The device is an electrocardiograph which can collect 12-lead ECG signal simultaneously and print ECG waveform with thermal printing system. It features in, recording and displaying ECG waveform in auto/manual mode, measuring, analyzing and diagnosing ECG waveform parameters automatically, prompting for "Lead off" and "Lack of paper", English menu, built-in lithium battery, AC/DC, selecting rhythm lead freely to observe abnormal heart rhythm, case database management.

### Functions:

- Adopt 10.1" high-resolution color LCD, operate by full-touch screen and function buttons, convenient to use
- Sync collection for 12-lead ECG, support 12-lead and Cabrera-lead waveform display. Get high-quality ECG waveform via power frequency filter (50/60 Hz), baseline filter and EMG filter (25/35 Hz) of ECG signal.
- Display 3/6/12-lead ECG waveform, HR value, print mode, sensitivity, paper speed, filter state, clock, battery power, background grid, measured parameter and interpretation, etc. on one screen, and with prompt functions for lead-off, overload and system working state
- Power supply: AC/DC (which can adapt AC 50/60Hz), with built-in rechargeable lithium battery, charging circuit, perfect over-current and over-voltage protection circuit
- In optimal DC state, up to 10-hour standby time, continuous print for more than 3-hour, which meets the requirements of visiting a patient at home and body examination
- English operating interface, the report can be printed, built-in virtual keyboard
- Built-in thermal printer, support multiple printing modes and formats, such as, auto M x N, M x N+1, M x N + 2, rhythm M line and manual, etc. The contents including time, paper speed, sensitivity, calibration signal, lead name, filter state and patient information can be printed. The waveform length printed, measurement parameter output, diagnosis conclusion, superimposed QRS waveform, histogram, trend chart and interval list, etc. can be set. With the print functions (timing print and auto-print for arrhythmia, etc.), which meet different requirements.
- With the functions of auto-measurement and auto interpretation for routine ECG parameters, provide measured results and auto-diagnosis conclusions for HR, PR interval, P Duration, QRS Duration, QT/QTc interval, P/QRS/T Axis, R (V5), S (V1), R (V5) + S (V1) amplitude, etc. which reduces the doctor' s burden.
- Built-in memory stores 4000 cases (data saving time: 10 s) at least, convenient for case review and statistics.
- Historical cases can be reviewed, inquired, modified, transmitted, printed, lead corrected, the cases can also be exported to other electronic file formats (dat, pdf, xml, bmp, jpeg, etc)
- Support external USB standard keyboard, mouse, scanner and Printer

### Package included:

- Cable ECG
- Ground wire
- 6x Suction cups for ECG
- 4x limb clamp ECG
- Printer paper
- Power supply cable
- User manual

• **Dimensions:** 340mm(L) x 320mm(W) x 86mm(H)

• **Net weight:** 5Kg

### Accessories / Spare Parts

- Code: 0803910 Cable ECG
- Code: 0803911 Ground wire
- Code: 0803912 6x Suction cups for ECG
- Code: 0803913 4x limb clamp ECG
- Code: 0803842 Printer Paper
- Code: 0803841 Trolley base

## ECG HOLTER 12 CHANNEL

ref. **0803844**



Dynamic ECG Systems adopt 12-lead system of international standard, which can continuously record ECG waveform for 24 hours and analyze ECG waveform by the PC software. It is applicable for use in medical institutions and the community.

### Function

- The recorder can store real-time ECG data for 24 hours.
- Arrhythmia analysis based on MCSSTM and TMCATM can greatly reduce doctor's work.
- Sync analysis for 12-lead, which ensures that QRS waveform can be found accurately, and no distortion.
- With more than 10 templates (such as atrial fibrillation, ect.) and many user-defined modules, which can almost distinguish any kind of pathologic waveforms.
- Flexible analysis channel selection function, which ensures that any channel can be selected as the main analysis channel.
- Flexible atrial fibrillation analysis, which ensures that physicians can use the whole or segmented Auto/manual analysis, more accurately in atrial fibrillation analysis.
- All pacemakers (such as AAI, VVI, DDD etc.) can be analyzed by the powerful pacemaker analysis function.
- With many pacing analysis templates, such as "Dual chamber pacing", "Ventricular pacing", "Ventricular Pseudo fusion" and "Ventricular asynchronous pacing", etc.
- Single of full-lead ECG in any time bucket can be reviewed by the fast review analysis function.
- With analysis functions for 5-minute, 1-hour and 24-hours heart rate variability.
- One-Key print, convenient and fast to print the reports.
- Perfect case management function.
- Abnormal ST can be analyzed according to the whole and event by the advanced ST segment analysis and total load of myocardial ischemia, which can help the doctors comprehensively judge myocardial ischemia.
- The risk of sleep breath pause can be predicted by the unique "Sleep breath pause syndrome" analysis function.
- The death risk in patients with myocardial infarction can be predicted by the "HRT" analysis function.
- "T-Wave alternation" is an important index to predict the malignant arrhythmia and sudden cardiac death.
- With QTD (QT dispersion), VCG (vector cardiogram), VLP (ventricular late potential) and TVCG (time vector cardiogram) analysis modules, more valuable in analysis report.

### Accessories / Spare Parts

- Code: 0803907 Electrode cable

### Package Contents:

- ECG cable
- 2x Packs of 10 ECG electrode pads
- Device case
- Usb cable
- Usb stick
- 2x Screws
- User manual

**Dimensions:** 69mm (L) x 50mm (W) x 17.3mm (H)

**Weight:** 100g



## Sleep Apnea Screen Meter

Ref: **0803843**



### Package included:

- Usb cable
- Nasal cannula
- 2x Cannula adaptor
- 2x SpO<sub>2</sub> sensor (Adult & Pediatric)
- User manual

**Dimension:** 69(L) x 50(W) x 17.3(H)mm

**Weight:** 100g

The product is suitable for sleep breathing, SpO<sub>2</sub> monitoring. Suitable for family, community medical institutions and hospitals.

### Specifications

#### Performance:

- Nose air flow measurement range: 0~40 rpm, accuracy:  $\pm 2$  rpm.
- SpO<sub>2</sub> measurement range: 0~100% (resolution 1%), accuracy: 70~100%,  $\pm 2$ ; Less than 70%, unspecified.
- PR measurement range: 30~250 bpm, accuracy:  $\pm 2$  or  $\pm 2\%$  (whichever is greater).
- Resistance to surrounding light: The deviation between the value measured in the condition of man-made light or indoor natural light.
- Accuracy under low perfusion: low perfusion 0.4%, SpO<sub>2</sub> accuracy:  $\pm 4\%$  PR accuracy  $\pm 2$ bpm or  $\pm 2\%$ , whichever is greater.
- Power consumption:  $\leq 100$ mA
- Power supply: 3.7V rechargeable lithium battery.
- Safety type: internally powered equipment, type BF applied part.

#### Features:

- SpO<sub>2</sub> and PR value display, pulse sound and bar graph indication.
- Pulse waveform and nose flow waveform display.
- Alarms for low power
- Figer out and exceeding limits
- Battery status indication
- Real-time clock
- Timing power on/off
- Storage of multiple cases
- Continuous 12 hours data record
- Uploading data via USV port
- PC analysis software

## PATIENT MONITOR WITH PRINTER

ref: **0803834**



The monitor has abundant functions that can be used for clinical monitoring with adults, pediatrics and neonate. Users may select different parameter configurations according to different requirements. The monitor, power supplied by 100-240V~, 50/60Hz, adopts 12.1" color TFT LCD displaying real-time, date and waveform. It can synchronously display eight-channel waveform and full monitoring parameters equipped with an optional 48mm thermal recorder. The monitor can be connected to the central monitoring system via wire or wireless network to form a network monitoring system.

This device can monitor such parameters as ECG, RESP, NIBP/IBP, SpO<sub>2</sub>, and dual-channel TEMP, etc. It integrates parameter measurement module, display and recorder in one device to form compact and portable equipment. Its replaceable internal battery brings lot of convenient for patient moving.

### Functions

- 12.1" TFT color LCD
- Fanless design, quiet, energy-saving and clean, which reduce the possibility of cross-infection
- All-round monitor for adults, pediatric and neonate
- With standard interface, oxygen graph, trend graph, big character interface and view bed, convenient to observe
- Finish all operations with keys and knobs
- Maximum 8-channel waveform synchronous display
- Display 7-lead ECG waveform on one screen, cascade ECG waveform display
- Adopt digital SpO<sub>2</sub> technology, anti-motion and anti-ambient light interference, and measurement can be performed under the circumstance of weak filling
- Heart rate variability (HRV) analysis function
- NIBP measurement mode: Manual/AUTO/STAT, storage for 4800-group NIBP data
- Review for 71 alarm events of all parameters and 60 arrhythmia alarm events
- Drug concentration calculation and titration table functions
- One-touch printing of trend graph
- Connect to Central Monitoring System with Wi-Fi or wired mode
- AC/DC, built-in rechargeable lithium battery achieves uninterrupted monitoring, autonomy of the battery: up to 3 hours
- Anti-high frequency surgical unit, defibrillation-proof (special leads are necessary)

### Functions:

#### Standard parameters:

ECG, RESP, SpO<sub>2</sub>, PR, NIBP, dual-channel TEMP

#### -ECG:

Heart rate (HR)  
ECG waveform  
Arrhythmia and ST-segment analysis

#### -RESP:

Respiration rate (RR)  
Respiration waveform

#### -SpO<sub>2</sub>:

Oxygen saturation (SpO<sub>2</sub>)  
Plethysmogram (PLETH) waveform  
Pulse rate (PR)  
Bar graph

#### -NIBP:

Systolic pressure (SYS),  
Diastolic pressure (DIA)  
Mean pressure (MEAN)

#### -TEMP: T1, T2, TD

#### -IBP (optional):

CH1: SYS, DIA, MAP

CH2: SYS, DIA, MAP  
IBP waveform

#### -CO<sub>2</sub> (optional):

ELCO<sub>2</sub>, InsCO<sub>2</sub> Inspired  
Minimum CO<sub>2</sub>,

#### AwRR: Air Way Respiration

It has abundant functions, such as audible and visual alarm, trend data storage and output, NIBP measurement, alarm event marking and drug concentration calculation, etc.

### Package included

- 2x Packs 10pcs Pads ECG electrodes
- ECG cable
- SpO<sub>2</sub> Sensor
- BP extending tube
- Cuff for adults
- Temperature sensor
- Printer paper 5pcs
- Power cable
- User manual

**Dimension: 319mm(L) x 167mm(W) x 268mm(H)**

**Weight: < 3kg (excluding accessories)**

### Accessories / Spare Parts

- 0803908 Cuff for adults
- 0803909 BP extending tube
- 0803914 ECG cable
- 0803915 SpO<sub>2</sub> Sensor
- 0803916 Temperature sensor
- 0803917 Printer paper
- 0803837 Trolley base
- 0803838 Bag
- 0803839 Packs 10pcs Pads ECG electrodes

## Manual Suction Unit 7B-1 (20 L/MIN)

Ref. 0803435



The Manual Suction Device is a portable and easy-to-use solution for immediate suction of fluids. It features an operating lever, ideal for emergency situations or use outside the doctor's office.

### Suction 7B-1 (20 L/MIN) (Manual)

- Max Negative Pressure:  $\geq 40$ kpa
- Max flow:  $\geq 20$ L/min
- Liquid Collection Jar: 200ml

### Dimensions

- G.W.: 1.5kg
- N.W.: 0.202kg
- Size: 17×18×5cm

### Warranty

- 2 Years for the Suction

### Package Contents

- Suction Tube 8mm: 1pc
- Suction catheter F8 (children): 1 pc
- Suction catheter F12 (adults): 1 pc
- User Manual

## Suction Unit 7E-A (20 L/MIN)

Ref. 0803436



Lightweight and quiet suction device, ideal for basic medical use and home care. Offers consistent performance, easy operation, and includes a fluid collection container. Suitable for daily support of patients with respiratory needs.

### Suction 7E-A, 20L/min (Electric)

- Input power: 90VA
- Power supply: AC230V, 50Hz
- Max Negative Pressure:  $85 \pm 10$ KPA
- Max air flow:  $\geq 20$ L/min (without jar)
- Operation noise: Approximately  $\leq 65$  dB (A)
- Liquid collection jar: 1000mL
- Operating time: 30 minutes on, 30 minutes off

### Dimensions

- G.W: 7kg
- N.W.: 4kg
- Size: 28×19,6×28,5cm

### Package Contents

- Suction tube (length 2 m,  $\Phi 7 \times \Phi 12$ ): 1 pc
- Suction catheter F8 (for children): 1 pc
- Suction catheter F12 (for adults): 1 pc
- Fuses ( $\Phi 5 \times 20$ ): 2 pcs
- Air filter: 2 pcs
- User manual: 1 pc

### Warranty

- 2 Years for the Suction

## Suction Unit 7E-C (20 L/MIN)

Ref. 0803437



Reliable and practical suction device, suitable for basic medical use and home care. It combines effective suction with low noise and easy transport. Ideal for daily support of patients with respiratory or other secretion removal needs.



### Suction 7E-C, 20L/min (Electric)

- Input power: 90VA
- Power supply: AC230V, 50Hz
- Max Negative Pressure: 85±10KPA
- Max air flow: ≥20L/min (without jar)
- Operation noise: Noise ≤65dB(A)
- Liquid collection jar: 1000mL
- Operating time: 30 minutes on, 30 minutes off

### Dimensions

- G.W.: 7kg
- N.W.: 3.7kg
- Size: 37×18×28cm

### Package Contents

- Suction tube (length 2 m, Φ7×Φ12): 1 pc
- Suction catheter F8 (for children): 1 pc
- Suction catheter F12 (for adults): 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- Power cable: 1 pc
- User manual: 1 pc

### Warranty

- 2 Years for the Suction

## Suction Unit 7E-H1 (28±4L/MIN)

Ref. 0803438



Efficient and quiet suction device, suitable for medical use in clinical or home environments. It features a durable fluid collection container and offers stable operation with easy use and maintenance. Ideal for daily respiratory support needs.



### Suction 7E-H1 (Electric)

- Input power: 120VA
- Max Negative Pressure: 85±5kPa
- Max air flow: 28±4L/min(Device inlet)
- Operation noise: ≤65dB(A)
- Liquid collection jar: 1000mL, 1piece(PC material)
- Operating time: 30 minutes on, 30 minutes off

### Dimensions

- G.W :7kg
- N.W.:3.9kg
- Size: 48×16,5×28,5cm

### Package Contents

- Suction tube (length 2 m, Φ7×Φ11): 1 pc
- Suction catheter F8 (for children): 1 pc
- Suction catheter F12 (for adults): 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- User manual: 1 pc

### Warranty

- 2 Years for the Suction

## Suction Unit 7E-G1 (27±4L/MIN) (with battery)

Ref. 0803439



Portable suction device with rechargeable battery, ideal for use in vehicles or outdoors. Features an oil-free pump, overflow protection, and wall mounting capability. Ideal for emergencies and travel, with an optional disposable system for maximum hygiene.

### Suction 7E-G1 (with Battery)

- DC power supply: AC 100-240V,50/60Hz or DC 12V
- Input power: 150VA
- Max Negative Pressure: 85±5kPa
- Air flow rate: 27L±4/min(Device inlet)
- Operation noise: ≤65dB
- Operating time: 30 minutes on, 30 minutes off
- Running time of battery operation: Approx. 60 minutes

### Dimensions

- G.W.: 7kg
- N.W.: 4kg
- Size: 48×18,9×28,5cm

### Package Contents

- Suction tube (Φ7×Φ11, length 2 m): 1 pc
- Suction catheter F8 (children): 1 pc
- Suction catheter F12 (adults): 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- Power cable: 1 pc
- User manual: 1 pc

### Warranty

- 2 Years for the Suction/ 1 year for the Battery

## Suction Unit 7A-23D (20 L/MIN)

Ref. 0803434



Effective and ergonomic suction device designed for safe removal of secretions in clinical and home environments. Ideal for healthcare professionals and home care, it offers reliable operation with easy transport and use.

### Suction 7A-23D, 20L/min (Electric)

- Power supply: AC230V, 50Hz
- Max negative pressure: 90±10KPA
- Operation noise: ≤60 dB(A)
- Input power: 180VA
- Max air flow: ≥20L/min
- Liquid collection jar: 2500mL/pc\*2

### Warranty

- 2 Years for the Suction

### Dimensions

- G.W :15.2kg
- N.W.:13.5kg
- Size: 36×32×48cm

### Package Contents

- Abdominal suction tube: 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- Flexible suction tube (Φ7×Φ12, length 2 m): 1 pc
- Power cord: 1 pc
- Foot switch: 1 pc
- User manual: 1 pc



## Suction Unit 7A-23A (32±4 L/MIN)

Ref. 0803464



The electric suction device is ideal for medical use, offering reliable and efficient operation. Easy to use and clean, with low noise levels and consistent performance.

### Suction 7A-23A 32±4L (without jar)

- Input power: 180VA
- Max Negative Pressure: 90±10KPA
- Max air flow: 32±4L (without jar)
- Noise: ≤ 60dBA
- Operating time: 30 minutes on, 30 minutes off
- Liquid collection jar: 2000mL/pc\*2

### Warranty

- 2 Years for the Suction

### Dimensions

- G.W :10.2kg
- N.W.:8.5kg
- Size: 37.4×30.9×49.9cm

### Package Contents

- Suction tube (Φ7×Φ12, length 2 m): 1 pc
- Abdominal suction tube: 1 pc
- Fuses (Φ5×20): 2 pcs
- Air filter: 2 pcs
- Power cord: 1 pc
- Foot switch: 1 pc
- User manual: 1 pc



# Wheelchairs with Medium Wheels

## Wheelchair "AEGEUS" 24"

ref. **0815003**

43cm

ref. **0815002**

46cm

ref. **0815001**

50cm



The "AEGEUS" wheelchair is designed to transport people with limited mobility and can be used by the patient themselves or with the help of a companion. They offer autonomy and ease of movement, depending on the model and the user's needs.



### Technical Specifications

FRAME	Steel - Foldable
COLOR	Blue
SEAT	Fabric
ARMRESTS	Foldable
FOOTRESTS	Detachable, Rotating, Swing Away with adjustable height
REAR WHEELS	24" MAG wheel with solid tire
FRONT WHEELS	8" x 2" Solid (200x25mm)
BRAKES	Attendant brakes

- Double cross bar
- Quick Release
- Option to add bottle holder 0815035
- Option to add I.V stand 0815033

## Wheelchair "AEGEUS" 12"

ref. **0815006**

43cm

ref. **0815005**

46cm

ref. **0815004**

50cm



Single cross support bar

<b>0815001</b>	52cm	50cm	40cm	44cm	70cm	110cm	94cm	51cm	24" 8"	125kg	19.5*kg
<b>0815002</b>	47cm	46cm	40cm	44cm	65cm	115cm	94cm	51cm	24" 8"	125kg	19.5*kg
<b>0815003</b>	44cm	43cm	40cm	44cm	63cm	112cm	94cm	51cm	24" 8"	125kg	19.5*kg

\* With footrests & rear wheels. The wheelchair without footrests and rear wheels weighs 13kg.

### Technical Specifications

FRAME	Steel - Foldable
COLOR	Blue
SEAT	Fabric
ARMRESTS	Foldable
FOOTRESTS	Detachable, Rotating, Swing Away with adjustable height
REAR WHEELS	12" Solid
FRONT WHEELS	8" x 2" Solid (200x25mm)
BRAKES	Attendant brakes

- Option to add bottle holder 0815035
- Option to add I.V stand 0815033

<b>0815004</b>	51cm	50cm	41cm	43cm	63cm	102cm	93cm	51cm	12" 8"	125kg	16*kg
<b>0815005</b>	47cm	46cm	41cm	43cm	61cm	100cm	93cm	51cm	12" 8"	125kg	15.5**kg
<b>0815006</b>	45cm	43cm	41cm	43cm	59cm	98cm	93cm	51cm	12" 8"	125kg	15.5**kg

\* With footrests. The wheelchair without footrests weighs 14.5kg.

\*\* With footrests. The wheelchair without footrests weighs 14kg.

## Commode chair "CYAN" (U-Shape)

Ref: **0805410**



### Technical Specifications

FRAME	Steel
COLOUR	Grey-Blue
SEAT	PVC
ARMRESTS	Adjustable-Foldable
FOOTRESTS	Detachable Adjustable Length Rotating
REAR CASTORS	4" Compact
FRONT CASTORS	4" Compact
BRAKES	On the Rear Castors

44,5cm	43cm	44cm	46cm	55cm	77cm	97cm	53cm	4" 4"	136kg	10,3Kg	

### Adjustable bathtub board with rotating base (110 kg)

Ref: **0812202**



75cm	34,5 cm	45-67cm	-	110Kg

- Weight: 1,40 kg
- Swivel base with diameter: 9,5 cm

## Commode Wheelchair

Ref: **0806777**



### Technical Specifications

FRAME	Steel
SEAT	PVC Leather, "U" Shape, with Commode
ARMRESTS	Foldable
FOOTRESTS	Detachable Length Adjustable
REAR WHEELS	5"x1" (127x25mm) Solid with Brake
FRONT WHEELS	5"x1" (127x25mm) Solid
ATTENDANT BRAKES	No
MAX LOAD	120 Kg

0806777	50cm	44cm	44cm	38cm	60cm	91cm	90cm	49cm	5" 5"	120kg	12kg

## Wheelchair “JADE”

Ref: **0815007**



“JADE” wheelchairs are designed for transporting people. They facilitate the movement of people with limited mobility and can be used with the assistance of a companion.

### Technical Specifications

FRAME	Steel - Foldable with folding backrest
COLOR	Green
SEAT	Fabric
ARMRESTS	Fixed
FOOTRESTS	Detachable, Rotating, Swing Away with adjustable height
REAR WHEELS	12"x 1.8" (305x45mm) Solid
FRONT WHEELS	6" x 1" (152x25mm) Solid

### Accessories (not included)

- Ref.: 0815008 Transfer Bag

0815007	46.5cm	47cm	40cm	40cm	65.5cm	90cm	92cm κλειστό 72cm	47cm	12" 6"	125kg	13.5*kg	

\*With footrests. The wheelchair without footrests weighs 12kg

## Foldable Transport Wheelchair “I-TRAVEL”

Ref: **0805662**



Adjustable suitcase-style handle

The “I-TRAVEL” stroller offers flexibility & modern design. It has a detachable backrest, a lightweight and foldable frame that makes it easy to transport and store.

0805662	43cm	39.5cm	36cm	36cm	51cm	83cm	86.5cm	48cm	8" 6"	136kg	7.5kg	

## Commode Wheelchair “RHINO” 24”

Ref: **0815018**



136kg



The “RHINO” wheelchairs of this type are designed to transport people, facilitating daily care and hygiene. They are ideal for use indoors, allowing the patient to be moved to the toilet or bathroom, either independently or with the help of an attendant. They also function as a toilet seat and are suitable for use in the shower, offering safety, comfort, and convenience for the user’s personal hygiene.



### Technical Specifications

FRAME	Aluminum - Non Foldable
COLOR	Blue
SEAT	PU seat, waterproof with commode
ARMRESTS	Flip up with waterproof PU pad
FOOTRESTS	Detachable, rotating, Swing away with adjustable height
REAR WHEELS	24" MAG solid mag wheel with Quick release
FRONT WHEELS	5" solid

- Suitable for use in the shower

<b>0815018</b>	49cm	47cm	43cm	42.5cm	66cm	98.5cm	95.6cm	53cm	24" 5"	136kg	17*kg

\* With footrests & rear wheels. The wheelchair without footrests and rear wheels weighs 11kg.

## Commode Wheelchair U-Shape “AEGEUS”

Ref: **0815020**



125kg



The wheelchair “AEGEUS” is designed to transport people, offering support and ease of movement in everyday life. It is suitable for people with mobility problems or limited mobility and can be used with the help of an attendant. They can also be used as a toilet seat, facilitating personal hygiene, especially indoors.



### Technical Specifications

FRAME	Steel – foldable
COLOR	Blue
SEAT	PVC with commode – “U” Shape
ARMRESTS	Flip back
FOOTRESTS	Detachable, rotating, swing away with adjustable height
REAR WHEELS	12" Solid
FRONT WHEELS	8" x 2" Solid (200x50mm)

- Single Crossbar

<b>0815020</b>	46cm	44cm	41cm	40.5cm	62cm	100cm	91cm	49cm	12" 8"	125kg	17*kg

\* With footrests. The wheelchair without footrests weighs 16kg.

## Commode Wheelchair “LAPIS”

Ref: **0815000**



The “LAPIS” wheelchair is designed to transport people with limited mobility, offering safety and comfort during their daily commute. Suitable for use by both the user themselves and with the help of a companion. In addition, it has a built-in container, allowing it to be used as a toilet seat.

### Technical Specifications

FRAME	Steel – foldable
COLOR	Blue
SEAT	PVC with commode – “U” Shape
ARMRESTS	Folding-Flip back
FOOTRESTS	Detachable, Rotating, Swing away with adjustable height
REAR WHEELS	24" MAG Solid & Quick release
FRONT WHEELS	8" x 1" Solid (200x25mm)
BACKREST	With folding backrest mechanism



<b>0815000</b>	48cm	46cm	40cm	46cm	64cm	109cm	102cm	51cm	24" 8"	125kg	19*kg	

\* With footrests & rear wheels. The wheelchair without footrests and rear wheels weights 11,5kg.

## Wheelchair “ONYX”

Ref: **0815010**



The “Onyx” wheelchair features a reinforced frame and high-strength materials, making it ideal for users with high demands!

### Technical Specifications

FRAME	Steel - Foldable
COLOUR	Black
SEAT	Fabric
SEAT & BACKREST FABRIC	With fire retardant properties
ARMRESTS	Detachable
FOOTRESTS	Detachable, rotating, swing away with adjustable height
REAR WHEELS	24" MAG Wheel solid
FRONT WHEELS	8" x 2" Solid (200 x 50mm)

- Double cross bar
- Option to add Reclining footrests 0815011,0815012

<b>0815010</b>	69cm	67cm	49cm	43cm	89cm	127cm	91cm	50cm	24" 8"	272kg	28*kg	

\* With footrests. The wheelchair without footrests weighs 26kg.



## Wheelchair AQUA-RECLINING


Ref. **0805498**

The **AQUA-RECLINING** wheelchair is ideal for transporting people with mobility problems or limited mobility. Suitable for indoor use, it facilitates daily care, access to the toilet and personal hygiene, with the help of a caregiver.



### Technical Specifications

FRAME	Aluminum - Non foldable
SEAT	PU Seat- Waterproof with Commode
HEIGHT ADJUSTMENT	Yes
SUITABLE FOR USE IN THE SHOWER	Yes
SEAT AND BACK TILT	Adjustable seat and backrest tilt from 0° to 28°
HEADREST	Yes, height adjustable
ARMRESTS	Flip up with PU Waterproof Pad
FOOTRESTS	Aluminium, detachable and rotating with adjustable length
REAR WHEELS	5"x 1" (127x25mm) solid with brake
FRONT WHEELS	5"x 1" (127x25mm) solid with brake
ATTENDANT BRAKES	No
MAX LOAD	136Kg



											GROSS WEIGHT
<b>0805498</b>	43.5cm	43cm	45cm	55cm	104.5 cm	118 - 123cm	56.5 61.5cm	5" 5"	136kg	14.5kg	16.5kg

## Lift chair electrical, “ESTIA IV”

Ref: **0811956**

**2** / **4** Years  
Warranty



FRAME	STEEL
SEAT & BACKREST FABRIC MATERIAL	NYLON WITH COTTON LINING
REAR CASTORS	3"
FRONT CASTORS	5"
BRAKES	IN FRONT AND REAR CASTORS
FOOTREST	FOOTREST PLATFORM
WITH HYDRAULIC LIFTING SYSTEM	THE HYDRAULIC PUMP CONTROLS THE LIFT OF THE SEAT BASE AND THE LEG LIFT SYSTEM.

<b>0811956</b>	46cm	43cm	43cm	37cm	59cm	88cm	93-109cm*	46-62cm*	3" 5"	120kg	28.5kg		
							102-119cm**	56-72cm**					

\* Measurement at the lowest position of the orange locking pin.  
\*\* Measurement at the highest position of the orange locking pin.

## Lift chair electrical, “ESTIA III”

Ref: **0811955**

**2** / **4** Years  
Warranty



### Technical Specifications

FRAME	Aluminum
FRAME	Foldable & the Seat Section is Detachable
BACKREST AND ARMREST MATERIAL	ABS
SEAT MATERIAL	Steel, welded around with high- strength tubes
SEAT CUSHION	Removable with washable cover
REAR WHEELS	3". With brake & 360o rotation capability
FRONT WHEELS	5", With brake
BRAKE	On the front & rear wheels
FOOTREST	Fixed Footplate (Platform)
MINIMUM RADIUS	77 Cm
WITH ELECTRIC LIFTING SYSTEM	The electric lifting system controls the lifting of the seat base from the vertical part of the frame
NOISE	≤ 65 dB
BATTERY	Lithium 1 pc X 24V X 4 AH (battery pack is hidden in bottom tube)
DURATION OF BATTERY	Battery life is estimated at 800 charge cycles.
POWER SUPPLY	IN: 220V AC, 50Hz, 30W / OUT: 24V,1A
IP GRADE	IP25
MAX LOAD	120 Kg
NET WEIGHT - WITHOUT THE SEAT CUSHION	23 Kg

46cm	40cm	57cm	77cm	87-115cm	46-76cm	3" 5"	120kg	23kg*				

\*Without the seat cushion

## Reclining Wheelchair “DIONE I”

Ref: **0803400**



Solid with spoked wheels & Quick Release

Seat Tilt



The DIONE wheelchair is a special type of wheelchair, with a comfortable ergonomic seat, and provides all the following adjustments & advantages for greater flexibility and safety.

- Roof spacer
- With anti-tilt system
- Escort bar
- Escort brakes
- Safety belt

### Technical Specifications 0803400 / 0803401

FRAME	Steel
RECLINING BACKREST	Adjustable from 95° to 145°
BACKREST HEIGHT	Adjustable
SEAT	Comfortable with adjustable seat depth
SEAT TILT	With adjustable seat tilt from 0° to 32°
ARMRESTS & SIDE PANELS	Height & width adjustable
CHEST SUPPORT	Detachable, adjustable
FOOTRESTS	Detachable, elevating & length adjustable with supporting pad
KNEE SUPPORT	Rotating knee support, adjustable with cushion
REAR WHEELS	<b>0803400</b> 24" Solid with spoked wheels & Quick Release <b>0803401</b> 12" Mag solid & Quick Release
FRONT WHEELS	8"x2" (200x50mm) solid, adjustable
MAX USER WEIGHT	150kg
HEAD SUPPORT	Height, width and angle adjustable
MATERIAL OF UPHOLSTERY, SEAT CUSHION & HEADREST CUSHION	Polyurethane (PU) foam (inside) with nylon lining

## Reclining Wheelchair “DIONE II”

Ref: **0803401**



12" Mag solid & Quick Release

Seat Tilt



MEASUREMENT WITH HEADREST IN MAXIMUM HEIGHT POSITION		
<b>0803400 / 0803401</b>	84cm	134cm

<b>0803400</b>	40cm	*49.5 - 51.5cm	**52.5 - 59cm	56-72cm	121cm	114cm	56cm	24" 8"	150kg	40kg	
<b>0803401</b>	40cm	*49.5 - 51.5cm	**52.5 - 59cm	56-72cm	121cm	114cm	56cm	12" 8"	150kg	38kg	

\*The depth of the seat is adjustable. The seat cushion has a depth of 49.5cm  
\*\* The backrest height is adjustable. The backrest cushion has a height of 56cm

## Reclining Wheelchair "PROTEUS I"

Ref. **0811980**

**2** / 4 Years  
Warranty



*The Reclining wheelchair "Proteus I" is designed for the transportation of people with mobility problems or limited mobility. It can be used both by the user and with the help of an attendant. The ergonomic design and reclining features offer increased comfort and support, making it ideal for everyday use.*

### Technical Specifications

FRAME	Steel
RECLINING BACKREST	Reclining backrest with tilt from 90° έως 160°
SEAT	PVC Leather
ARMRESTS	Foldable
SIDE PANELS	With side panel
FOOTRESTS	Detachable & elevating & length adjustable
REAR WHEELS	24" Spoke solid with Quick Release
FRONT WHEELS	8"x2" (200x50mm) PU solid
MAX USER WEIGHT	120kg



- Head support
- Anti-tilt wheels
- Attendant brakes
- Safety belt



											GROSS WEIGHT
0811980	46.5cm	44cm	74*cm	66cm	117*cm	125**cm	55cm	24" 8"	120kg	23.5kg	26.5kg

\* Measurement including footrests. The measurement without the footrests is 92cm.

\*\*Measurement including headrest. The measurement without the headrest is 98cm.





## Power wheelchair foldable “ODYSSEUS”

Ref: **0805700**



### Technical Specifications

FRAME	Aluminum- Folding (once battery is removed)
COLOUR	Black
SEAT	Detachable seat cushion
ARMRESTS	Flip back armrest with PU cushion
FOOTREST	Foldable
REAR WHEELS	12" with air pressure
FRONT WHEELS	8" Solid
MAX SPEED	6 Km/h
AUTONOMY*	20 Km
MAX TILT	10°
MOTOR	250W x 24 V x 2 Pieces / 1 year Warranty
BATTERY	20Ah x 24V x 1 Piece / 6 Month Warranty
CHARGER	24V x 3AMP / 1 year Warranty
BRAKES	Electromagnetic

\* Approximative autonomy as also determined by the weight of occupant, the battery condition, the ambient temperature, the way of driving, the terrain and the wind speed

43cm	41cm	41cm	46cm	61cm	95cm	103cm	46cm	12" 8"	135kg	28kg	

## Power Wheelchair Foldable "FELIX"

Ref: **0805762**



*The FELIX electric wheelchair is the best choice for people with limited mobility who are looking for a combination of durability, flexibility and ergonomics in their daily mobility. Made of high-quality lightweight aluminum, the FELIX offers durability and at the same time easy folding, ideal for storage and transport in a car. It is suitable for indoor or outdoor use*



### Technical Specifications

Frame	Aluminum, folding
Controller	Right or Left Adjustable Controller
Motor	250W x 24V x 2 with electromagnetic brake
Battery	1 piece Lithium 24V 12Ah
Charger	24V 2A
Max speed	6 km/h
Travel mileage	20-25 km
Gradeability	13°
Front wheel	8" PU solid
Rear wheel	12" PU solid
Armrest	Flip back armrest with PU cushion
Footrest	Foldable – Height adjustable
Seat	Detachable seat cushion with fire retardant properties



- Control Lever That Can Be Adjusted To The Right Or Left
- Anti-Tip Casters
- Safety Belt
- Mud Fenders
- Easy To Transport

### Warranty

- 2 years for the frame
- 1 year for the electric parts of the wheelchair
- 6 months for the battery

<b>0805762</b>	47cm	46cm	43cm	50-60cm	64cm	93cm	95-105cm	51cm	12" 8"	130kg	25kg*

\*Net weight without batteries. The weight with batteries is 29kg.



## Scooter-Type "ARES"

Ref: **0811820**

*Mobility scooters like "Ares" offer an essential blend of independence, convenience, and comfort for individuals with limited mobility. Designed to empower users to move freely, they open up possibilities for navigating the world with confidence and ease.*



FRAME	Steel
COLOUR	Silver
SEAT	With fire retardant properties- Rotating seat
ARMRESTS	Folding
FOOTREST	Foldable
REAR WHEELS	10" with air pressure
FRONT WHEELS	10" with air pressure
MAX SPEED	12 Km/h
AUTONOMY*	50 Km
MAX TILT	9°
MOTOR	340W
BATTERY	(12V x 2 pieces) 52 Ah
CHARGER	24V x 6AMP
BRAKES	Electromagnetic

\* Approximative autonomy as also determined by the weight of occupant, the battery condition, the ambient temperature, the way of driving, the terrain and the wind speed



0811820	1255 x 583 x 965mm	10" x 10"	450 x 460 x 490mm	83Kg	150Kg	(12Vx 2 ) 52Ah	9°	12km/h	50km

## Foldable Scooter Wheelchair "AEGLE"

Ref: **0805761**

AEGLE is the ideal folding electric scooter for people with limited mobility who seek independence, comfort and style in their everyday life. With a modern design and lightweight construction, it folds easily in a few seconds thanks to its ergonomic design, offering ease of transport and storage. It is suitable for indoor or outdoor use.



Rotary switch for speed control

Foldable

### Warranty

- 2 years for the frame
- 1 year for the electric parts of the wheelchair
- 6 months for the battery



Seat tilt 4°



### SPECIFICATIONS

Frame	Aluminum - Foldable
Maximum speed (driving forward)	6 km /h
Autonomy *	15km
Maximum obstacle height	75 mm
Maximum safe slope	12 °
Dimensions (overall) (L x W x H)	(97 x 52 x 95) cm
Net weight including batteries	25 kg
Seat tilt	4°
Motor	180W
Batteries	1 pc - 24V 10Ah (lithium)
Battery charger	DC output 24 V, 2A
Battery charger protection class	IP20
Reversing alarm	Yes
Brakes	Electromagnetically
Rear wheel diameter	8 " solid
Diameter of front wheels	7 " solid
Control console	Rotary switch for speed control
Recommended Storage & Use Temperature	-30°C to +45°C



- Fold and unfold effortlessly with just the push of a button on the remote control
- Designed for easy transport and convenient storage
- Option to change the plastic trim to blue or red



## Scooter PERSEUS

Ref: **0805760**

The **PERSEUS** electric scooter is a reliable solution for the daily mobility of people with limited mobility. With its modern design, excellent autonomy and advanced safety features, PERSEUS offers independence and comfort without compromise.

**2** / **4** Years  
Warranty



Rotary switch for speed control



Steering Wheel Inclination up to 135°



Basket



Reclining  
Backrest  
up to 125°

Suspension  
Front & Back



Anti-tilt wheels

Backlights

FRAME	Steel
COLOUR	White
SEAT	Swivel, height-adjustable and depth-adjustable (back and front)
ARMRESTS	Raised and width adjustable
REAR WHEELS	13" with air pressure
FRONT WHEELS	11.5" with air pressure
MAXIMUM SPEED	15 Km/h
AUTONOMY*	55 km
MAXIMUM OBSTACLE HEIGHT	60 mm
LEG PLATFORM HEIGHT	170 mm
MAXIMUM SAFE SLOPE	15 °
MOTOR	800W
BATTERIES	2 тμχ 24V -55Ah (lead acid) 6 Months Warranty
CHARGER	24V x 6AMP
BRAKES	Electromagnetic

\* Approximative autonomy as also determined by the weight of occupant, the battery condition, the ambient temperature, the way of driving, the terrain and the wind speed



# PERSEUS



Power switch

Steering wheel adjustment lever

Charging socket

2 USB Ports for charging devices (e.g. Smartphones)



<b>0805760</b>	1330 x 670 x 1350mm	11.5" 13"	650-750mm	480 x 500 x 655 mm	106Kg*	190Kg	(24V x 2) 55Ah	15°	15km/h	55km

\*Net weight without batteries. The weight including batteries is 142 Kg.

## Electric Wheelchair "ROBODRIVE"

Ref. **0805720**



## Ergonomically Crafted and Visually Appealing Design



### SPECIFICATIONS

<b>FRAME</b>	Folding
<b>COLOUR</b>	White
<b>SEAT</b>	Fabric Cushion
<b>ARMRESTS</b>	Foldable
<b>FOOTREST</b>	Fixed on the frame & foldable
<b>REAR WHEELS</b>	10" Pneumatic Gear tires
<b>FRONT WHEELS</b>	10" Off-road omnidirectional wheels
<b>MAX SPEED</b>	6Km/h
<b>ENDURANCE MILEAGE*</b>	25Km
<b>MAX TILT</b>	10°
<b>MOTOR</b>	200W x 24 V x 2 Pieces
<b>BATTERY</b>	20Ah x 24V x 1 Piece / 6-Month Warranty
<b>CHARGER</b>	29.4V x 4AMP
<b>BRAKES</b>	Electromagnetic

\*Approximative autonomy as also determined by the weight of occupant, the battery condition, the ambient temperature, the way of driving, the terrain and the wind speed.

- Anti-tilt Casters.
- Safety Belt.
- Lithium Battery.
- Mud fender.
- Autonomous shock absorption system
- Seat upholstery made of memory foam.
- Automatic folding mechanism.
- Intelligent management featuring application.
- Smart Monitoring with error notifications.
- Folding Dimensions (L x W x H): (83 x 63 x 37) cm.
- Easy application (app) for managing functions via mobile.

<b>0805720</b>	47cm	43cm	43cm	45cm	65cm	108cm	95cm	55cm	10" 10"	120kg	44kg	



## Electric Wheelchair "TECHRIDE"

Ref. **0805721**



*Sophisticated and functional,  
tailored specifically for your needs*

### SPECIFICATIONS

FRAME	Folding
COLOUR	Green
SEAT	Fabric Cushion
ARMRESTS	Foldable
FOOTREST	Fixed on the frame & foldable
REAR WHEELS	10" Pneumatic Gear tires
FRONT WHEELS	8" Solid
MAX SPEED	8Km/h
ENDURANCE MILEAGE*	23Km
MAX TILT	15°
MOTOR	200W x 24 V x 2 Pieces
BATTERY	20Ah x 24V x 1 Piece / 6-Month Warranty
CHARGER	29.4V x 4AMP
BRAKES	Electromagnetic

\*Approximative autonomy as also determined by the weight of occupant, the battery condition, the ambient temperature, the way of driving, the terrain and the wind speed.

- Anti-tilt Casters.
- Safety Belt.
- Lithium Battery.
- Mud fender.
- Autonomous shock absorption system.
- Seating made from memory foam.
- Manual folding mechanism.
- Intelligent management featuring application.
- Smart Monitoring with error notifications.
- Folding Dimensions (L x W x H): (70 x 63 x 45) cm.

0805721	49cm	42cm	43cm	47cm	62cm	100cm	94cm	50cm	10" 8"	150kg	27kg



## Electric Wheelchair "TURBOTECH"

Ref. **0805722**



### SPECIFICATIONS

FRAME	Folding
COLOUR	Grey
SEAT	Fabric Cushion
ARMRESTS	Foldable-Detachable
FOOTREST	Fixed on the frame-Foldable & Height Adjustable
REAR WHEELS	12" Pneumatic Gear tires
FRONT WHEELS	12" Off-road omnidirectional wheels
MAX SPEED	10Km/h
ENDURANCE MILEAGE*	31Km
MAX TILT	10°
MOTOR	200W x 24 V x 2 Pieces
BATTERY	20Ah x 24V x 1 Piece / 6-Month Warranty
CHARGER	29.4V x 4AMP
BRAKES	Electromagnetic

\*Approximative autonomy as also determined by the weight of occupant, the battery condition, the ambient temperature, the way of driving, the terrain and the wind speed.

- Anti-tilt Casters.
- Safety Belt.
- Lithium Battery.
- Mud fender.
- Autonomous shock absorption system.
- Seating made from memory foam.
- Manual folding mechanism.
- Intelligent management featuring application.
- Smart Monitoring with error notifications.
- Folding Dimensions (L x W x H): (86 x 64 x 53) cm.

<b>0805722</b>	47cm	42cm	43cm	48cm	62cm	105cm	120-130cm	50cm	12" 12"	150kg	36kg

*Designed to deliver the pinnacle of  
off-road driving experience!*



## Stand-Up Wheelchair "MOVEUP PRO MANUAL"

Ref. **0805723**





*Simple to maintain an upright & comfortable position.  
Self-sufficiency in everyday living.*



### SPECIFICATIONS

Color	White
Maximum Speed	4.5 Km/h
Driving Distance Range	Under 25°C, 15 Km with an even speed, on even ground
Motor	250W x 24 V x 2 Pieces
Maximum safe slope	5°
Height of obstacle	≤ 4 cm
Width of trench	≤ 5 cm
User's height	1.5m - 1.9m
Minimum turning Diameter	0.63 m
Bluetooth	2.4GHz, 1mW
Battery Capacity/Voltage	15AH/25.6V / 6 months warranty
Charger	DC 24V 4A
Input of the charger	AC 100V-240V/50-60Hz



- Lithium Battery
- Remote control through application.
- Manual lift mechanism

<b>0805723</b>	36cm	22cm	29cm	53cm	78cm	120cm	70cm	5" 10"	120kg	86kg

## Stand-Up Wheelchair "MOVEUP PRO ELECTRICAL"

Ref. **0805724**





*Simple to maintain an upright & comfortable position.  
Self-sufficiency in everyday living.*



### SPECIFICATIONS

Color	White
Maximum Speed	4.5 Km/h
Driving Distance Range	Under 25°C, 15 Km with an even speed, on even ground
Motor	250W x 24 V x 2 Pieces
Maximum safe slope	5°
Height of obstacle	≤ 4 cm
Width of trench	≤ 5 cm
User's height	1.5m - 1.9m
Minimum turning Diameter	0.63 m
Bluetooth	2.4GHz, 1mW
Battery Capacity/Voltage	15AH/25.6V / 6 months warranty
Charger	DC 24V 4A
Input of the charger	AC 100V-240V/50-60Hz



- Lithium Battery
- Remote control through application.
- Electrical lift mechanism

<b>0805724</b>	36cm	22cm	29cm	53cm	78cm	120cm	70cm	10" 5"	120kg	86kg

## Stand-Up Wheelchair "JOY"

Ref. **0805725**




*Simple to maintain an upright & comfortable position.  
Self-sufficiency in everyday living.*



### SPECIFICATIONS

Color	White
Maximum Speed	4.5 Km/h
Driving Distance Range	Under 25°C, 15 Km with an even speed, on even ground
Motor	250W x 24 V x 2 Pieces
Maximum safe slope	5°
Height of obstacle	≤ 4 cm
Width of trench	≤ 5 cm
User's height	1.5m - 1.9m
Minimum turning Diameter	0.63 m
Bluetooth	2.4GHz, 1mW
Battery Capacity/Voltage	15Ah/25.6V / 6 months warranty
Charger	DC29.4V, 4A
Input of the charger	AC220V



- Lithium Battery
- Remote control through application.
- Electrical lift mechanism

0805725	27cm	20cm	20cm	20cm	53cm	78cm	107cm	44cm	5" 10"	60kg	91kg



# Reclining Armchairs

## Reclining Armchair "CINNAMON I"

Ref: **0815700**



COVER MATERIAL	Soft chenille fabric
INTERNAL LINING MATERIAL	High density foam -the seat has a density of 35D, hardness 70
SIDE POCKET	Yes
CONTROLLER	Yes
USB PORT ON THE CONTROLLER	Yes
MAX LOAD	135kg
MOTORS	1
CASTORS	2 without brakes
MOVEMENTS	Raising armchair and simultaneously reclining backrest and raising footrest

• Option to add covers Code: 0815709

## Reclining Armchair "CLOUD I"

Ref: **0815701**



<b>0815700, 0815701</b>	52cm	50cm	67,5cm	70cm	73cm	110cm	126cm	166cm	113cm	135kg 45kg

• Option to add covers Code: 0815707



## Reclining Armchair "CINNAMON II"

Ref: **0815702**



• Option to add covers Code: 0815711

COVER MATERIAL	Soft chenille fabric
INTERNAL LINING MATERIAL	High density foam -the seat has a density of 35D, hardness 70
SIDE POCKET	Yes
CONTROLLER	Yes
USB PORT ON THE CONTROLLER	Yes
MAX LOAD	150kg
MOTORS	2
CASTORS	2 with brakes
MOVEMENTS	Raising armchair and simultaneously reclining backrest and raising footrest ,raising backrest, reclining backrest

## Reclining Armchair "CLOUD II"

Ref: **0815703**



• Option to add covers Code: 0815710

<b>0815702, 0815703</b>	50cm	54cm	67,5cm	80cm	74cm	104cm	126cm	166cm	114cm	150kg	51kg			

## Reclining Armchair "CLOUD XXL"

Ref. **0815704**

**3** Years Warranty

"CLOUD XXL" is designed to offer comfort and safety, and is ideal for use at home, in clinics, or in elderly care facilities. The reclining function" allows for rest and relaxation, while its sturdy construction supports different weights, even users with increased physical support needs. With an ergonomic design and discreet appearance that suits any space



COVER MATERIAL	Soft chenille fabric
INTERNAL LINING MATERIAL	High density foam -the seat has a density of 35D, hardness 70
SIDE POCKET	Yes
CONTROLLER	Yes
USB PORT ON THE CONTROLLER	Yes
MAX LOAD	250kg
MOTORS	2
CASTORS	2 with brakes
MOVEMENTS	Simultaneously raising armchair and footrest lowering ,simultaneously lowering chair and raising footrest, raising backrest, reclining backrest



• Option to add covers Code: 0815712



<b>0815704</b>	75cm	55cm	67,5cm	112cm	92cm	110cm	127cm	167cm	132cm	250kg	64kg

# Reclining Armchairs

## Reclining Armchair with Vibration Massage


Ref. **0815708**

**COMING SOON**

The RECLINING ARMCHAIR WITH VIBRATION MASSAGE is designed to provide comfort and safety, it is ideal for use at home, in clinics or in elderly care facilities. The reclining function allows for rest and relaxation, while the integrated vibration massage function offers additional relief. Its sturdy construction supports different weights, even users with increased physical support needs. With an ergonomic design and a discreet appearance that fits into any space



COVER MATERIAL	Leather, waterproof & anti-stain
INTERNAL LINING MATERIAL	High density foam -the seat has a density of 35D, hardness 70
SIDE POCKETS	Yes
CONTROLLER	Yes
USB PORT ON THE CONTROLLER	Yes
MAX LOAD	135kg
MOTORS	1
FUNCTIONS	Massage and Heating
CASTORS	2 without brakes
MOVEMENTS	Raising armchair and simultaneously reclining backrest and raising footrest



0815708	52cm	50cm	67,5cm	70cm	73cm	110cm	126cm	166cm	113cm	135kg	45kg

## Electric Patient Lift “ELYROS”


Ref: **0811880**



The “ELYROS” heavy-duty lift is specially designed for the safe and comfortable lifting and moving of people with limited mobility and increased body weight. It is suitable for use in hospitals, care units, and even at home, offering convenience and comfort to both patients and caregivers.



- Rechargeable battery
- Electric lift
- Steel frame
- Brakes on rear wheels
- Net weight: 45 kg
- Input Voltage: 220V – 50Hz
- Low battery indicator
- Safe lift stop if the patient or caregiver accidentally blocks the arm
- Battery warranty: 6 months



- A. Maximum lifting arm height (upper end): 200cm
- B. Minimum and maximum lifting height (distance of lifting arm hooks from the ground): 90 - 172cm
- C. Height of the bottom surface of the base from the ground: 4.5cm
- D. Height of the top surface of the base from the ground: 11cm
- E. Minimum depth: 125cm
- F. Internal length of base: 80cm
- G. Minimum inside width of closed base: 49cm
- H. Minimum outer width of closed base: 60cm
- I. Bend radius range: 131cm
- L. Internal width of open base: 80cm
- M. External width of open base: 91cm
- Diameter of front wheels: 7.5cm
- Diameter of rear wheels with brakes: 10cm



## Bath Chair “WAVE” Rose

Ref: **0803851** **Small**

Ref: **0803852** **Medium**

Ref: **0803853** **Large**



The “WAVE” bath chair is designed to provide stable and comfortable support for children and teenagers during bath or shower time. It is ideal for indoor use, inside the bathroom, and is intended for users who have difficulty sitting or maintaining their balance.



### Technical Specifications

FRAME	Aluminum, Foldable
COLOR	Rose
FABRIC	Polyester mesh fabric with fire retardant properties
BASE	Aluminum, with anti-slip support section

	0803851	0803852	0803853
MAX LOAD	100Kg	100Kg	100Kg
SEAT DEPTH	28cm	35cm	40cm
USEFUL SEAT WIDTH	43cm	43cm	43cm
BACKREST HEIGHT	49cm	60cm	71.5cm
CALF SECTION LENGTH	25cm	29.5cm	34cm
OVERALL LENGTH	100*cm	123*cm	142.5*cm
OVERALL WIDTH	43.8cm	43.8cm	43.8cm
SEAT HEIGHT (MEASUREMENT FROM THE FLOOR)	15-40cm	15-46cm	15-51cm
AGE	1-6 years	4-12 years	9-16 years

\*In full horizontal extent





## Bath Chair "WAVE" Blue

Ref: **0803848** **Small**

Ref: **0803849** **Medium**

Ref: **0803850** **Large**



The "WAVE" bath chair is designed to provide stable and comfortable support for children and teenagers during bath or shower time. It is ideal for indoor use, inside the bathroom, and is intended for users who have difficulty sitting or maintaining their balance.

### Technical Specifications

FRAME	Aluminum, Foldable
COLOR	Blue
FABRIC	Polyester mesh fabric with fire retardant properties
BASE	Aluminum, with anti-slip support section

	0803851	0803852	0803853
MAX LOAD	100Kg	100Kg	100Kg
SEAT DEPTH	28cm	35cm	40cm
USEFUL SEAT WIDTH	43cm	43cm	43cm
BACKREST HEIGHT	49cm	60cm	71.5cm
CALF SECTION LENGTH	25cm	29.5cm	34cm
OVERALL LENGTH	100*cm	123*cm	142.5*cm
OVERALL WIDTH	43.8cm	43.8cm	43.8cm
SEAT HEIGHT (MEASUREMENT FROM THE FLOOR)	15-40cm	15-46cm	15-51cm
AGE	1-6 years	4-12 years	9-16 years

\*In full horizontal extent



## Base with Wheels for Bath Chair “WAVE”

Ref: **0803854**

Small

Ref: **0803855**

Medium

Ref: **0803856**

Large



*Combines with WAVE bath chairs for use in the shower*



### Technical Specifications

<b>FRAME</b>	Aluminum, with anti-slip support section
<b>REAR CASTORS</b>	3" swivel & anti-slip
<b>FRONT CASTORS</b>	3" swivel & anti-slip
<b>BRAKES</b>	In all castors

	<b>0803854</b>	<b>0803855</b>	<b>0803856</b>
<b>MAX LOAD</b>	100Kg	100Kg	100Kg
<b>LENGTH</b>	112cm	115cm	117,5cm
<b>WIDTH</b>	49,8cm	49,8cm	49,8cm
<b>HEIGHT</b>	43cm	43cm	43cm
<b>AGE</b>	1-6 years	4-12 years	9-16 years

## Commode Chair “SEAL”

Ref: **0815017**



The “SEAL” commode chair is intended for use as a toilet aid by people with limited mobility, such as the elderly, patients, or people with disabilities. It is designed for indoor use and offers safety and comfort in daily hygiene.



- Steel frame
- Suitable for all types of toilets with removable bucket
- Adjustable height in 6 levels by push button system
- Seat and backrest material: Vinyl with flame retardant properties

43cm	46cm	58-62cm	48-63cm	34.5cm	8kg	130kg

## 360° ROTATING SHOWER CHAIR

Ref: **0805660**



### Technical Specifications

• Frame: Aluminum

48-51cm	48-51cm	40cm	34cm	43-55cm	34cm	136kg	4,9kg

## 360° ROTATING BATHTUB TRANSFER BENCH

Ref: **0805661**



*Horizontal and 360° movement for easy transport in and out of the bath*

### Technical Specifications

• Frame: Aluminum

87cm	87cm	40cm	34cm	48-58cm	34cm	136kg	6,2kg

### Adjustable bathtub board with rotating base (110 kg)

Ref: **0812202**



75cm	34,5 cm	45-67cm	-	110Kg

- Product weight: 1,40 kg
- Swivel base with diameter: 9,5 cm

## Electric Bath Lift Chair “SPLASH”

Ref: **0815720**



Adjustable height

### Package Contents:

- Controller
- Charger
- Backrest
- Base (cushion and 6 suction cups)
- User Manual

The electric bath lift chair “SPLASH” offers comfort, safety, and independence in the bathroom, making it easier for people with limited mobility to use the bathtub. With a compact, foldable design and battery operation, it can be easily installed in any bathtub, without cables or external power supply. It features smart safety functions, non-slip suction cups, and an articulated seat system for maximum stability and ease of use.

### Specifications

OPERATING TEMPERATURE	5-40°C
OPERATING HUMIDITY	Below 40°C 20%-90%
STORAGE TEMPERATURE	-10-50°C
STORAGE HUMIDITY	20%-90%
CHARGING VOLTAGE	100-240 V, 50/60 Hz - 0,8 A
BATTERY	14.8 V, 2.2 Ah
IP RATING	IP68



- Easy seat height adjustment using a control lever (6-48 cm)
- Easy to install in almost any bathtub — no plumbing or electrical work required
- Removable backrest with 50° tilt (90°-140°)
- Easy to transport and store due to its small size when folded
- Equipped with 6 suction cups on the base for non-slip action
- Waterproof frame that can be raised and lowered during bathing

107cm	70cm	6-48cm	61cm	36.6cm	45cm	64cm	29cm	0°-140°	136kg	13.5kg

## Inflatable Bath with Pump

Ref: **0804616**



- Bathtub dimensions (LxWxH): 190x90x22cm
- Material: PVC
- Bathtub weight: 2Kg
- Pump weight: 200gr
- Capacity: 35lt
- Pump power consumption: 150W
- Voltage: 220V
- Double air chambers provide protection against water and foam leakage
- The water drain allows for quick and safe emptying of the bathtub after use

### Warranty

- 1 year for the inflatable bathtub
- 1.5 years for the pump

### Package Contents:

- Inflatable bathtub
- Electric pump
- Repair kit

## Bed Urine Bag with T-Valve

Ref: **0809544**



- With non-return valve
- With T-tap valve
- Capacity 2 lt
- Internal pipe diameter: 6 mm
- Pipe length 90 cm
- Step connector with protective cap
- Distributed in packs of 25 pieces

## Disposable urine bag without outlet

Ref. **0809012**



- With non-return valve
- Without outlet valve
- Capacity 2 lt
- Internal pipe diameter: 6 mm
- Pipe length 90 cm
- Step connector with protective cap
- Distributed in packs of 25 pieces

## Bed Pan


Ref: **0810224**



- Length: 37cm
- Width: 23.8cm
- Height: 8.5cm
- Material: PP
- Net Weight: 110gr
- Capacity: 1L

## Reaching Aid with Grabber

Ref: **0812216** 80cm



- 1. Head
- 2. Main Tube
- 3. Trigger
- 4. Handle

### Technical specifications

- Tube: anodized aluminum
- Header: ergonomic PP (polypropylene) and TPR
- Handle: ergonomic PP (polypropylene) and TPR
- Trigger: PP (polypropylene)



## Kid's wheelchair "HELIOS"

Ref. **0815009**



The "HELIOS" children's wheelchair is designed for transporting children. It facilitates the movement of people with limited mobility and can be used either with the help of an attendant or by the child themselves.

<b>FRAME</b>	Steel – foldable
<b>SEAT</b>	Fabric - Nylon
<b>ARM REST</b>	Flip back
<b>FOOTREST</b>	Detachable, Rotating, Swing away with adjustable height
<b>REAR WHEELS</b>	22" MAG wheel with solid tire & Quick release
<b>FRONT WHEELS</b>	6" x 1.25" Solid (150x32mm)
<b>SAFETY BELT</b>	Yes
<b>HANDLES</b>	With Adjustable Height


• With optional reclining footrests, 0815031, 0815032



<b>0815009</b>	36cm	35cm	36cm	34cm	57cm	94cm	87-96*cm	44cm	22" 6"	68Kg	16**Kg

\*With the handles in the lowest position, the total height is 87cm. With the handles in the highest position, the total height is 96cm.

\*\*With footrests & rear wheels. The wheelchair without footrests and rear wheels weights 12.5kg.



## Knee Walker “STROLLY”

Ref: **0815019**



The “STROLLY” walker is designed to support the knee and improve reduced walking ability. It is intended for users with sufficient mobility to be able to move outside the home for daily activities and social participation without risk. It provides stability when walking and has a specially designed basket for carrying small loads, offering practicality and independence in everyday use.



### Specifications:

- Frame: Steel - Foldable
- Handles: Height-Adjustable
- Brakes: Safety brakes on handles that can lock
- Front and Rear Wheels: 7.5", Solid
- Seat for user: Height adjustable, with PU pad
- With steel storage basket - removable
- Total Dimensions (folded) (L x W): 13.5 x 34 x 5 cm
- Max Load for Basket: 5 Kg



40cm	20cm	43 – 58cm	82-106cm	92cm	42cm	10*Kg	136Kg

\*The measurement is with the basket. Without the basket the Weight of the product is 9.5 Kg.


## Rollator Walker "MIGHTY"

Ref: **0812215**



The Rollator walker MIGHTY features a reinforced frame and high-strength materials, making it ideal for users with high demands!

FRAME	Aluminum
FRAME	Foldable
COLOUR	Black
HANDLES	Height adjustable, in 6 positions
ATTENDANT BRAKES	Safety brakes on handles that can lock
REAR WHEELS	7.3"x1.4" (188x36mm) Solid
FRONT WHEELS	7.3"x1.4" (188x36mm) Solid
SEAT FOR USER	Yes
CANE HOLDER	Yes
STORAGE BAG	Yes
SUSPENTION	No
MAX LOAD	180 Kg
MAX LOAD FOR BAG	2.5 Kg



56cm	62cm	62cm	55cm	86-99cm	8.3Kg	180Kg

## Rollator Walker "INDIGO"

Ref: **0805411**



- Height of Handles Adjustable, in 6 Positions
- Safety Brakes on Handles that can Lock
- Includes Cane Holder
- Includes Storage Bag



### Technical Specifications

FRAME	Aluminum
FRAME	Foldable for easy storage
COLOUR	Blue
HANDLES	Height adjustable, in 6 positions
ATTENDANT BRAKES	Safety brakes on handles that can lock
REAR WHEELS	9" Solid
FRONT WHEELS	9" Solid
SEAT FOR USER	Yes
CANE HOLDER	Yes
STORAGE BAG	Yes
SUSPENTION	No
MAX LOAD	136 Kg
MAX LOAD FOR BAG	5 Kg

45.5cm	57cm	64,5cm	58cm	85-98cm	8Kg	136Kg

\*The measurement is without the bag.

## Rollator Walker "RUBY"

Ref: **0805499**



- Height of Handles Adjustable, in 6 Positions
- Safety Brakes on Handles that can Lock
- Includes Cane Holder
- Includes Storage Bag

### Technical Specifications

FRAME	Aluminum
FRAME	Foldable for easy storage
COLOUR	Red
HANDLES	Height adjustable, in 6 positions
ATTENDANT BRAKES	Safety brakes on handles that can lock
REAR WHEELS	8" Solid
FRONT WHEELS	8" Solid
SEAT FOR USER	Yes
CANE HOLDER	Yes
STORAGE BAG	Yes
SUSPENTION	Yes
MAX LOAD	136 Kg
MAX LOAD FOR BAG	5 Kg

45.5cm	64cm	67cm	60cm	80-93cm	8.5Kg	136Kg

\*The measurement is without the bag.

## Heavy duty foldable walker

Ref: **0812200**



60cm	48cm	46cm	82-107cm	3Kg	225Kg

- Frame: Aluminum
- Folding for easy use
- Height-Adjustable
- PVC handgrips
- Compatible with Legs with castor (Ref: 0812207)



## Bariatric tetrapod cane Large base (225 kg)

Ref: **0812205**

30x21cm	74-97cm	225Kg

- Stainless Steel



## Bariatric tetrapod cane Small base (225 kg)

Ref: **0812206**

21x16cm	74-97cm	225Kg

- Stainless Steel



## Legs with castor for bariatric walker

Ref: **0812207**

Compatible with: 0812200




## Walker with 360° Swivel Castors

Ref: **0810643**



- Aluminum Frame
- Foldable
- Height Adjustable
- PVC handgrips
- Swivel castors with 360-degree rotation.
- Compatible with wheel conversion accessory 0815815

55cm	44,5cm	51cm	84-102cm	3.9Kg	136Kg



## Walker 360o accessory with double castors brake (pair)

Ref: **0815815**

Accessory for converting the wheels of the walker 0810643 in dual wheels with rubber tips that can act as brakes (self-locking tips)



**New Product**

## Bariatric underarm crutches Large (250 kg)

Ref: **0812203**



140-160cm	250Kg

- Frame: Steel


## Heavy Duty Forearm Crutch

Ref: **0812217**



- Frame: Steel
- Adjustable Height
- Net Weight: 1.1Kg

107.5 -137.5 cm	23-30.5cm	220Kg



## Premium Four-Leg Red Cane with Closed Handle


Ref: **0805663**



- Frame: Aluminum with plastic handle (TPR)
- Ergonomic handle
- Adjustable height
- Non-slip base
- Net Weight: 0.82Kg
- Can be used on either the left or the right hand (ambidextrous)



26x20cm	71-97cm	136Kg



## Bariatric underarm crutches Medium (250 kg)

Ref: **0812204**



120-140cm	250Kg

- Frame: Steel



## Underarm Crutches Universal Height-Adjustable

Ref: **0804594**



- Aluminum frame
- Universal size ( Fits sizes Small, Medium, Large)
- Adjustable Height 98-158cm
- Tube thickness 1,15mm

98-158cm	136Kg

## Premium Four-Leg Blue Cane with Closed Handle


Ref: **0805664**



- Frame: Aluminum with plastic handle (TPR)
- Ergonomic handle
- Adjustable height
- Non-slip base
- Net Weight: 0.82Kg
- Can be used on either the left or the right hand (ambidextrous)



26x20cm	71-97cm	136Kg



## Chest Support Belt

Ref. **0803650**

UNIVERSAL  
- ONE SIZE



A: 43cm  
B: 16cm

Total max length  
(with straps): 203cm

- breathable material
- padded fabric for maximum comfort
- adjustable straps to fit all body sizes

A

Ideal for comfortable sitting posture on chairs and wheelchairs.

## Body Support Belt

Ref. **0803651**

UNIVERSAL  
- ONE SIZE



A: 40cm  
B: 28cm  
C: 15cm

- breathable material
- padded fabric for maximum comfort
- adjustable straps to fit all body sizes



## Heel Protector with Cushion

Ref. **0803652**

UNIVERSAL  
- ONE SIZE



- Polyester Fibers
- Prevents Bed Sores, Suitable for Pressure & Friction Relief
- Dimensions (LxWxH): (26 x 21 x 11) cm

## Electric Multi Crank Bed “LIBRA II”

Ref: **0811241**

**2** / **4** Years  
Warranty



### Technical Specifications

- Metal frame
- Aluminum rails in wood color
- External wooden panel covers covering the lifting mechanisms
- Features Trendelenburg & Anti-Trendelenburg 12° function
- Product weight 96kg
- Maximum backrest angle 70°
- Maximum leg angle 30°
- Input power supply 220V-50Hz, output 24V
- Includes 3 motors and 1 remote control
- Foldable, transportable, easily and quickly assembled & disassembled by 1 person
- Metal orthopedic slats (bed with wooden slats 0803152)
- Electric back, knee, and leg elevation
- Hanger included
- Compact wheels with brakes
- Safety rails that can be easily attached or removed (no tools required)
- Warranty for electrical parts: 1 year

	Int. Dimensions (LxW)	Ext. Dimensions (LxW)	Height from floor	Safe Working Load (SWL)*	Back / legs / knees elevation
<b>0811241</b>	200x90cm	215x103cm	Electric Elevation 39-80cm	180Kg	Electrical

*This includes the combined weight of the patient, mattress, and any\* accessories or equipment placed on the bed*

## Electric Multi Crank Bed “VIRGO II”

Ref: **0811240**



### Technical Specifications

- Metal frame
- Aluminum rails in wood color
- Internal wooden covers for the headrests
- Features Trendelenburg & Anti-Trendelenburg 12° function
- Product weight 90kg
- Maximum backrest angle 70°
- Maximum leg angle 30°
- Input power supply 220V-50Hz, output 24V
- Includes 3 motors and 1 remote control
- Foldable, transportable, easily and quickly assembled & disassembled by 1 person
- Metal orthopedic slats (bed with wooden slats 0803152)
- Electric back, knee, and leg elevation
- Hanger included
- Compact wheels with brakes
- Safety rails that can be easily attached or removed (no tools required)
- Warranty for electrical parts: 1 year

	Int. Dimensions (LxW)	Ext. Dimensions (LxW)	Height from floor	Safe Working Load (SWL)*	Back / legs / knees elevation
<b>0811240</b>	200x90cm	212x103cm	Electric Elevation 39-80cm	180Kg	Electrical

*This includes the combined weight of the patient, mattress, and any\* accessories or equipment placed on the bed*

## Wooden folding massage bed - 2 Sections - With accessories - BLUE (Width 70)

Ref: **0812300**



### SPECIFICATIONS

- Dimensions - Open (LxWxH): **185 x 70 x 61-86 cm**
- Dimensions - Closed (LxWxH): **93 x 71 x 16.5 cm**
- Net Weight: **15.5 kg**
- Material: Beech wood, PVC leather, Polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Bed
- Side armrest
- Front armrest
- Adjustable headrest
- Headrest cushion with hole
- Carrying bag



Ref: <b>0812300</b>	<b>Blue</b>	15,5Kg	200Kg
Ref: <b>0812301</b>	<b>Black</b>	15,5Kg	200Kg
Ref: <b>0812302</b>	<b>Beige</b>	15,5Kg	200Kg

## Wooden folding massage bed - 2 Sections With accessories - BLACK (Width 60)

Ref: **0812303**



### SPECIFICATIONS

- Dimensions - Open (LxWxH): **180 x 60 x 61-86 cm**
- Dimensions - Closed (LxWxH): **90 x 61 x 16.5 cm**
- Net Weight: **14.3 kg**
- Material: Beech wood, PVC leather, Polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Bed
- Side armrest
- Front armrest
- Adjustable headrest
- Headrest cushion with hole
- Carrying bag

## Accessories

### Semi-cylindrical Cushion

Ref: **0812314**



- Dimensions (LxWxH): **66 X 22 X 11.5 cm**

### Semi-cylindrical Cushion

Ref: **0812315**



- Dimensions (LxWxH): **66 X 15 X 15 cm**

### Paper Holder (massage bed)

Ref: **0812316**



- Dimensions (LxWxH): **72 X 15 X 3 cm**

## Wooden folding massage bed - 2 Sections Without accessories - BLACK (Width 60)

Ref: **0812304**



### SPECIFICATIONS

- Dimensions - Open (LxWxH): **180 x 60 x 61-86 cm**
- Dimensions - Closed (LxWxH): **90 x 61 x 16.5 cm**
- Net Weight: **12.3 kg**
- Material: Beech wood, PVC leather, Polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Bed
- Carrying bag

# Massage Beds

## Aluminum folding massage bed 3 Sections - With accessories - BLACK (Width 60)

Ref: **0812305**



### SPECIFICATIONS

- Dimensions - Open (LxWxH): **180 x 60 x 56.5-78 cm**
- Dimensions - Closed (LxWxH): **90 x 61 x 16.5 cm**
- Net Weight: **15.9 kg**
- Material: Aluminum, PVC leather, Polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Bed
- Side armrest
- Front armrest
- Adjustable headrest
- Headrest cushion with hole
- Carrying bag

## Wooden Bed - BEIGE (Width 70)

Ref: **0812307**



### SPECIFICATIONS

- Dimensions - Open (LxWxH): **185 x 70 x 65-85 cm**
- Net Weight: **30.6 kg**
- Material: Beech wood, PVC leather, Polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Bed
- Adjustable headrest-backrest
- Headrest cushion
- Bed frame legs
- Shelf for placing items
- Bag with screws

## Aluminum folding massage bed 2 Sections- With accessories - BLACK (Width 70)

Ref: **0812306**



### SPECIFICATIONS

- Dimensions - Open (LxWxH): **185 x 70 x 56.5-78 cm**
- Dimensions - Closed (LxWxH): **93 x 71 x 16.5 cm**
- Net Weight: **15.2 kg**
- Material: Aluminum, PVC leather, Polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Bed
- Side armrest
- Front armrest
- Adjustable headrest
- Headrest cushion with hole
- Carrying bag

## Foot stool (158 kg)

Ref: **0812201**



- Frame: stainless steel
- Surface: anti-slip & mat, made of rubber
- Maximum User Weight: 158Kg
- Total dimensions: 41 x 30 x 23 cm
- Base surface dimensions: 37 x 26 x 2 cm

## Folding massage chair – BLACK

Ref: **0812308**



### SPECIFICATIONS

- Dimensions - Open (LxWxH): **75 x 49 x 104-124 cm**
- Dimensions - Closed (LxWxH): **98 x 20.5 x 50**
- Net Weight: **15,9 kg**
- Material: Iron structure, PVC leather, polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Massage chair
- Headrest
- Headrest Cushion
- Chest support cushion

## Height adjustable “PONY” type stool with backrest & wheels - BLACK

Ref: **0812313**



### SPECIFICATIONS

- Seat Diameter: **34.5 cm**
- Adjustable base height: **45-60 cm**
- Net Weight: **5,6 kg**
- Material: Iron structure, PVC leather, nylon, polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Seat
- Casters
- Air pole lifting base
- Adjustable headrest
- Backrest

## Height adjustable stool with backrest & wheels - BLACK / BEIGE

Ref: **0812311**

Ref: **0812312**



### SPECIFICATIONS

Ref: **0812311, 0812312**

- Seat Diameter: **34.5 cm**
- Adjustable base height: **45-60 cm**
- Net Weight: **7,2 kg**
- Material: Iron structure, PVC leather, nylon, polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Seat
- Casters
- Air pole lifting base
- Leg ring - Footrest
- Backrest

## Height adjustable stool with backrest & wheels - BLACK / BEIGE

Ref: **0812309**

Ref: **0812310**



### SPECIFICATIONS

Ref: **0812309, 0812310**

- Seat Diameter: **34.5 cm**
- Adjustable base height: **45-60 cm**
- Net Weight: **3,7 kg**
- Material: Iron structure, PVC leather, nylon, polymeric foam
- For internal Use

### PACKAGE CONTENTS

- Seat
- Casters
- Air pole lifting base

## Air Cell Cushion "ANESIS I" TPU

Ref: **0804514**

Ref: **0804515**



• Dimensions (WxLxH): 42 x 42 x 5cm

• Dimensions (WxLxH): 40 x 46 x 5cm

Ref: **0804516**

• Dimensions (WxLxH): 46 x 46 x 5cm



- Silky feel and softness for maximum soft tissue protection
- Cover with air permeability, easy to wash & dry, with fire retardant properties (6 months guarantee)
- Lightweight and easy to carry
- Relieves and Relieves User Pressure from the Sitting Position
- Internal airflow automatically adjusts during sitting position, depending on body structure and weight, for a better fit
- Easy adjustment of the cushion to the desired position (with special straps)
- Protective TPU (thermoplastic polyurethane) cushion coating for easy cleaning
- Manual Air Piston (6 months warranty)
- Ergonomic design



CUSHION PRODUCT CODE	CUSHION DIMENSIONS (LxWxH)	SUITABLE FOR WHEELCHAIRS WITH A SEAT WIDTH:
0804514	42 x 42 x 5cm	41cm - 43cm
0804515	40 x 46 x 5cm	43cm - 46cm
0804516	46 x 46 x 5cm	46cm



## Air Cell Cushion "ANESIS II" TPU

Ref: **0804517**

Ref: **0804518**



• Dimensions (WxLxH): 42 x 42 x 5cm

• Dimensions (WxLxH): 40 x 46 x 5cm

Ref: **0804519**

• Dimensions (WxLxH): 46 x 46 x 5cm



- Silky feel and softness for maximum soft tissue protection
- Cover with air permeability, easy to wash & dry, with fire retardant properties (6 months guarantee)
- Lightweight and easy to carry
- Relieves and Relieves User Pressure from the Sitting Position
- Internal airflow automatically adjusts during sitting position, depending on body structure and weight, for a better fit
- Easy adjustment of the cushion to the desired position (with special straps)
- Protective TPU (thermoplastic polyurethane) cushion coating for easy cleaning
- Ergonomic design
- Cushions offer the opportunity of automatically adjustment of air flow, according to user needs based on innovative design.



CUSHION PRODUCT CODE	CUSHION DIMENSIONS (LxWxH)	SUITABLE FOR WHEELCHAIRS WITH A SEAT WIDTH:
0804517	42 x 42 x 5cm	41cm - 43cm
0804518	40 x 46 x 5cm	43cm - 46cm
0804519	46 x 46 x 5cm	46cm - 48cm

## Wall-mounted metal First Aid box with one shelf and lock

Ref: **0805407**

- With 1 Shelf
  - Lock with Key
  - Material: Metal
  - For internal use, avoid contact with water
  - Dimensions (LxWxH): 30 x 12 x 40 cm
- Available without Contents**



## Wall-mounted metal First Aid box with two shelves and lock

Ref: **0805406**

- With 2 Shelves
  - Lock with Key
  - Material: Metal
  - For internal use, avoid contact with water
  - Dimensions (LxWxH): 30 x 12 x 47 cm
- Available without Contents**



## First Aid Kit MOBIAK - Large

Ref: **0804602**

- **Contains:** Scissors, Plastic Tweezers, Sterile Gauze Sponges, PBT Conforming Bandages, Non-woven Tape, Surgical Mask, Vinyl Gloves, Cleansing Wipe, NaCl 0.9%, Adhesive Plaster, Transparent PU Wound Dressing, Triangular Bandage, Safety Pins, Burn Dressing, Emergency Blanket, Rain Poncho, Pen Flashlight, Notebook, EVA Case L, List of contents, Elastic Strap Tourniquet, CPR Face Shield, First Aid Booklet

- **Dimensions (LxWxH):** 22.5x16.5x5.5cm



## First Aid Kit MOBIAK - Small

Ref: **0804600**

- **Contains:** Scissors, Plastic Tweezers, Sterile Gauze Sponge, PBT Conforming Bandages, Non-woven Tape, Vinyl Gloves, Cleansing Wipe, NaCl 0.9%, Adhesive Plaster, Transparent PU Wound Dressing, Triangular Bandage, Safety Pins, EVA Case S, List of contents, First Aid Booklet

- **Dimensions (LxWxH):** 17x11x4cm



## AED Cabinet with Alarm

Ref: **0805495**

- **Dimensions (LxWxH):** 42,5 x 18 x 42,5cm



## First Aid Kit MOBIAK - Medium

Ref: **0804601**

- **Contains:** Scissors, Plastic Tweezers, Sterile Gauze Sponge, Pbt Conforming Bandages, Non-Woven Tape, Vinyl Gloves, Cleansing Wipe Nacl 0.9%, Adhesive Plaster, Transparent Pu Wound Dressing, Triangular Bandage, Safety Pins, Ischemic Tourniquet, CPR Mask, EVA Case-M, List Of Contents, First Aid Booklet

- **Dimensions (LxWxH):** 19 x 14 x 5cm



## Transport stretcher chair (2 in 1)

Ref: **0805919**



### Technical Specifications

- Frame: Aluminium
- Fabric: PVC (orange color), waterproof
- Dimensions open (LxWxH) : 192 x 50 x 20cm
- Dimensions in sitting position of the stretcher:
  - Backrest height: 89 cm
  - Seat height (from floor): 52cm
  - Seat surface width: 39cm
  - Backrest surface width: 48cm
- Dimensions folded (LxWxH) : 106 x 51 x 13cm
- Product Weight: 9.6Kg
- Weight Tolerance: 159Kg
- Lightweight, easy to carry.
- Safe to use and easy to sterilize.
- Includes two orange safety belts. (code 0805927)
- With two wheels of diameter 150mm
- Foldable
- Can be converted from a transport chair to a wheeled stretcher for patient transport and in horizontal position

## Folding Stretcher with wheels

Ref: **0805918**



### Technical Specifications

- Frame: Aluminum
- Fabric: CANVAS (blue color), waterproof
- Dimensions open (LxWxH): 185 x 50 x 20cm
- Dimensions closed (LxWxH): 92.5 x 50 x 9cm
- Product Weight: 6.5Kg
- Maximum Load: 159Kg
- Lightweight, foldable, small size, easy to carry
- Safe to use and easy to sterilize
- Includes two black safety belts (code 0804595)
- Features two 100mm diameter wheels and stabilizers

## Folding Aluminium Stretchers

Ref: **0804591**

Ref: **0805996**

Available Upon Request



### Technical Specifications

- Frame: Aluminium
- Material: PVC
- Includes: 2 Safety Belts
- Packing Dimensions (LxWxH) : 104 x 18 x 11 cm
- Dimensions Open (LxWxH) : 208 x 55 x 14cm
- Product Weight: 4,78Kg
- Weight Tolerance: 159Kg

### Safety belt, black, with plastic buckle (1 pair)

Ref: **0804595**

- Compatible with: 0804593, 0805910, 0805918, 0805907



Length: 170 cm  
Width: 5 cm



## Spine Board - Floating stretcher

Ref: **0805907**



### Spider straps system (Child)

Ref: **0804597**

- Compatible with: **0805907**



### Technical Specifications

- The stretcher is made of high-density polyethylene – HDPE. This makes it rigid, lightweight and reliable for use in emergency situations.
- It has 18 handles around the perimeter for operator convenience.
- Its material is durable with high strength and non-absorbent. It is also easy to clean and disinfect.
- The maximum Load is 170Kg.
- It Includes four (4) safety belts (**Ref: 0804595**), quick-release, with plastic buckle and 5 cm width for the restraint of the injured person.
- It can be used on boats, ambulances, hospitals, in sports and mountaineering clubs and generally outdoors for the transport of patients and injured people.
- It can float on water
- It can be used in combination with the Head Immobilization Splint 0805908.
- It can be used in combination with the inflatable Set (leg and hand) splints with set **Ref: 0805909**
- It is compatible for use with X-ray and MRI devices. (X-Ray & MRI compatible)
- It can be used in combination with Spider straps system **0804597 & 0804599**
- Product Net Weight: 7 Kg

## Basket Stretcher

Ref: **0805906**



### Technical Specifications

- The stretcher is made of high-density polyethylene – HDPE and has aluminum tubes around the perimeter. This makes it robust, flexible and reliable for use in emergency situations.
- Its material is durable with high strength and non-absorbent. It is also easy to clean and disinfect.
- The maximum Load is 250Kg.
- It includes four (4) safety belts, with a car-type lock, 5cm width for tying up the injured and a lifting system with four lashing straps with safety rings, which fit into the special round holes on the stretcher.
- It has an adjustable footrest mechanism, which helps with both horizontal and vertical transport.
- Contains internally, an integrated soft, non-absorbent coating 1cm thick for the comfort of the wearer, which is easily cleaned and disinfected.
- It is foldable, you can divide it into two parts so that it can be stored and transported easily.
- It is suitable and functional for rescue needs either in the mountains or at sea.
- It is not compatible for use with X-ray and MRI devices. (X-Ray & MRI non-compatible)
- Dimensions: 218x62x16 cm
- Dimension folded: 124x62x20 cm
- Product Net Weight: 20Kg with all accessories on and 15Kg without

## Scoop Stretcher

Ref: **0805910**



### Technical Specifications

- The stretcher is made of anodized aluminum with aluminum locking mechanisms, which makes it lightweight, durable, flexible, and reliable for use in emergency situations.
- Its material is with high strength and non-absorbent. It is also easy to clean and disinfect.
- It is foldable, splits into two parts so that it can be stored and transported easily.
- It has 6 handles around the perimeter for operator convenience, but both the front and rear can be used as handles.
- The stretcher frame is made in two halves that can be detached and attached under the patient without requiring the patient to be moved.
- The length is adjustable in four positions #0:168cm / #1: 178cm / #2: 189cm / #3: 200cm and can be adapted to most body types.
- The maximum Load is 170Kg.
- It Includes four (3) safety belts (**Ref: 0804595**), quick-release, with plastic buckle and 5 cm width for the restraint of the injured person
- It can be used on boats, ambulances, hospitals, in sports and mountaineering clubs and generally outdoors for the transport of patients and injured people.
- It is not compatible for use with X-ray and MRI devices. (X-Ray & MRI non-compatible)
- Dimensions: 168-200x44x6 cm
- Dimension folded: 120x44x8 cm
- Product Net Weight: 7.3Kg

## Adjustable cervical collar first aid

Ref: **0805920**



### Technical Specifications

- X-ray, CT and MRI suitable
- Suitable for easy access for pulse, airway and tracheal measurement
- 4 positions of height adjustment
- Dimensions open in full extent (LxH) : 56.5 x 10-15 cm
- Internal Lining: Foam material with a 100% cotton stockinet
- Product Weight: 140gr

## Adjustable cervical collar first aid (child)

Ref: **0805925**



### Technical Specifications

- X-ray, CT and MRI suitable
- Suitable for easy access for pulse, airway and tracheal measurement
- 3 positions of height adjustment
- Dimensions open in full extent (LxH): 48 x 8-11 cm
- Internal Lining: Foam material with a 100% cotton stockinet
- Product Weight: 70gr

## Head immobilizer adult

Ref: **0805908**



### Technical Specifications

- It consists of two sections made of foam material with vinyl lining, which are fixed to a base with straps so that it can be adapted to different types of stretchers.
- The Head Immobilizer 0805908 can be used with the Spine Board 0805907 and the Scoop Stretcher 0805910.
- The head strap is designed to conform patient's forehead.
- There are holes for monitoring the patient's ear.
- It is designed to support the injured person's head naturally, avoiding as much as possible extra pressure on the head.
- Minimal interference with X-ray, MRI, or CT procedures.
- Waterproof plastic coating makes it easy to prevent bacterial growth within components. Will not absorb blood or bodily fluids. Durable and easy to clean.
- Direct contact with skin, can be used in very cold weather.
- It can be used on boats, ambulances, hospitals, in sports and mountaineering clubs and generally outdoors for the transport of patients and injured people
- Product Net Weight: 1 Kg
- Dimensions: 40,5x27x14,5 cm
- Product parts: Base Board with eight "D" shape loop x1, Block with round hole x 2, Immobilization Belt x2, Head Strap x2
- Temperature range: -20°C ~ 50°C

## Soft Stretcher

Ref: **0804592**



### Technical Specifications

- Frame: Water-resistant And Washable Nylon
- Includes: 2 Safety Belts
- Dimensions when deployed (LxW) : 205 x 70 cm (measurement without the handles)
- Dimensions when closed (LxW) : 74 x 35 cm
- Product Weight: 1.72Kg
- Weight Tolerance: 159Kg

## Air splint set

Ref: **0805909**



### Technical Specifications

- Air splints are designed for fractures, broken, crushed, sprains or burned extremities.
- Air splints are provided to people who need protection and support for painful, swollen, or weak joints and their surrounding structures.
- Immobilizes, cushions, and protects injured limbs so the user can be transported safely & easily.
- The Air Splint are made of a soft and transparent plastic material, that is non-toxic, very resistant and easily adapts to limbs.
- The splints are made of two linings and between them there is a space which can contain the air that enters through a valve.
- Ready to use-inflates in seconds. Application is faster, simpler, and more comfortable than conventional bandaging/splinting techniques.
- The air fills the splint till creating a strong tension that avoids any movement of the limb. The inflatable splints have an outside zip to make the application to the limb easier.
- It is compatible for use with X-ray and MRI devices. (X-Ray & MRI compatible)
- Can be used in combination with any stretcher (0805906, 0805907, 0805910, 0804591, 0804592, 0804593) and or the Head Immobilizer 0805908.
- Temperature range: 0°C ~ 50°C
- Product Net Weight: 1,2 Kg
- It is compatible for use with X-ray and MRI devices. (X-Ray & MRI compatible)

## First Aid Medi Bag

Ref: **0805970**



- Material: 600D oxford fabric, plastic (bottom of bag)
- Dimensions: (53 x 32 x 24) cm
- Capacity: 40L

## First Aid Medi Mini Bag

Ref: **0805971**



- Material: 420D PVC, plastic (Handle)
- Dimensions: (25x16x10) cm
- Capacity: 4L

## First Aid Medi Backpack

Ref: **0805972**



- Material: 1200D oxford fabric, plastic (buckle of shoulder straps)
- Dimensions: 48x31x11cm
- Capacity: 16L

## First Aid Kit DIN13164 MOBIAK

Ref: **0804605**



### DESCRIPTION

- Adhesive Tape
- Adhesive Bandage 10x6cm
- Fingertip Bandage 4x5cm
- Finger Bandage 12x2cm
- Plaster Strip 7,2x1,9cm
- Plaster Strip 7,2x2,5cm
- First Aid Dressing DIN13151-G
- First Aid Dressing DIN13151-M
- First Aid Dressing DIN13151-K
- First Aid Dressing DIN13152-A
- Elastic Bandage DIN 61634-FB 6
- Elastic Bandage DIN 61634-FB 8
- Wound Compresses
- Triangular Bandage
- Emergency Blanket
- Scissor
- Disposable Gloves
- Face Mask
- Cleansing Wipe
- Contents List
- General Guide for First Aid Procedure

## Aluminum Stair Chair Climber

Ref. **0805976**




The aluminium Stair Chair Climber is specially designed for transporting patients up and down stairs, allowing quick and safe removal in case of emergency. An essential tool for ambulances, hospitals and sports facilities, it is only used by specially trained personnel. Its robust construction and ergonomic design make it a reliable solution in critical moments.



### Technical Specifications

- Frame Material: Aluminum
- Designed with special tracks, for convenient and safe transport in stairs
- Dimensions open (WxLxH): 163 × 51 × 118cm
- Dimensions folded (LxWxH): 106 × 51 × 30cm
- Carton Dimensions (LxWxH): 109 × 55 × 31 cm
- Lightweight and foldable
- Equipped with four wear-resistant rubber wheels
- Can also be used as a wheeled stretcher on level ground
- Product Gross Weight: 24Kg
- Product Net Weight: 21Kg



## Electric Aluminum Stair Chair Climber

Ref. **0807110**



The aluminium electric Stair Chair Climber is specially designed for transporting patients up and down stairs, allowing quick and safe removal in case of emergency. An essential tool for ambulances, hospitals and sports facilities, it is only used by specially trained personnel. Its robust construction and ergonomic design make it a reliable solution in critical moments.

### Technical Specifications

- Designed with special tracks, for convenient and safe transport in stairs
- Easy cleaning and disinfection
- Lightweight and foldable
- Equipped with four wear-resistant rubber wheels
- Length of rare foldable handle: 26cm
- Front adjustable retention handle with button: 15cm-35cm
- Height of adjustable grab bar: 108cm, 127cm, 146cm, 164cm
- Elevating footrest
- Safety belt
- Charger: 29.4V, 2A
- Battery: 24V / 13.2Ah
- Working power: AC 100-240V 50/60Hz
- Motor: DC 24V / 120W
- IP protection class: IP33
- Frame Material: Aluminum
- Max user weight: 180Kg
- Product gross weight: 39Kg
- Product net weight: 34Kg



## Electric Stair Climbing System for Wheelchair

Ref: **0807109**


The electric Stair Climbing System with a wheelchair crawler is specially designed for transporting patients upstairs, allowing quick and safe evacuation in case of emergency. It is intended for wheelchair users. An essential tool for ambulances, hospitals and sports facilities, used only by specially trained personnel. Its robust construction and ergonomic design make it a reliable solution in critical moments.



### Technical Specifications

- Frame material: Aluminium
- Tracks incline : 20°
- Designed with special tracks, for convenient and safe transport in stairs
- Easy cleaning and disinfection
- Lightweight and foldable
- Equipped with four wear-resistant wheels
- Height adjustable grab bar: 110-165CM with 3 height positions
- Adjustable width of the wheelchair placement section: 40-80CM
- Charger: 29.4V 2A
- Battery: 24V 13AH Lithium-ion battery
- Working power: AC 100-240V 50/60HZ
- Motor: 24V 200W
- LCD Display on the control panel
- Tracks speed: 3-8m/min
- LED Light for easy transport in the dark
- Dimensions closed (H\*L\*W): 109\*45\*30cm
- Product gross weight: 48Kg
- Product net weight: 36Kg

### Requirements for stairs



## Resuscitation Bag PVC AMBU Type

Ref: **0805921** **Infant**

Ref: **0805922** **Child**

Ref: **0805923** **Adult**



AMBU resuscitation bags are used to provide artificial respiration to patients who are not breathing or are breathing inadequately to assist breathing in emergency situations.

### PVC Infant Resuscitation Bag

- Suitable for infants weighing 5-10 kg
- Infant silicone mask
- PVC manual ventilation balloon with capacity of 280 ml
- Special safety valve to control the pressure up to 40 cm H<sub>2</sub>O
- 1600 ml oxygen bag
- Instructions for use (GR-EN)

### PVC Resuscitation Bag for Children

- Suitable for children weighing 10-40 kg
- Silicone mask for children
- PVC manual ventilation balloon with capacity of 550 ml
- Special safety valve to control the pressure up to 40 cm H<sub>2</sub>O
- 1600 ml oxygen bag
- Instructions for use (GR-EN)

### PVC resuscitation bag for adults

- Suitable for adults weighing 40 kg and above
- Adult silicone mask
- PVC manual ventilation balloon with capacity of 1500 ml
- Special safety valve to control the pressure up to 60 cm H<sub>2</sub>O
- 2000ml oxygen bag
- Instructions for use (GR-EN)

## Guedel Type Oropharyngeal Airway Set (Set of 8)

Ref: **0807104**



### SPECIFICATIONS

- The airways are color coded to help recognize the sizes; the length of the tube is marked on the airway with a number.
- The flange limits the depth of insertion -so can't accidentally lose the tube in patient's airway.
- The reinforced bite block prevents occlusion of the airway by clenching of the teeth.
- The rigid construction is kink resistant.

## Box for AMBU Resuscitation Bags

Ref: **0805926**



- Dimensions (LxWxH): 27x16x10.5cm

## Wall Bracket First Aid Box MOBIAK - Medium

Ref: **0804603**



Dimensions (WxLxH): 29 x 20 x 9cm

Available without contents

## Wall Bracket First Aid Box MOBIAK - Large

Ref: **0804604**



Dimensions (WxLxH): 42 x 32 x 15cm

Available without contents

## CPR Pocket Resuscitator Mask

Ref: **0805924**



CPR mask is designed for CPR, provides a physical barrier between the rescuer and victim.

## Compression Band for Bleeding Control Tourniquet

Ref: **0804611**



The Compression Band is intended for the stabilization of bleeding.

## CPR Keychain

Ref: **0805973**



- Ideal for a gift to friends and partners.
- With certified CPR mask.
- Dimensions: WxLxH (7 x 2 x 7)cm

## Resque Tube

Ref: **0807100**



- Material: Polyurethane Foam & Vinyl
- Dimensions (LxWxH): (124 x 16.3 x 7)cm
- Net weight: 1 Kg

## Hemostatic Tourniquet

Ref: **0807102**



- Material: Nylon & polyester
- Dimensions (LxW): 3.8 x 96cm
- Net Weight: 84 gr
- It is easy to use with one hand and it has a recording section for time monitoring, for preventing tissue damage due to prolonged hemostasis
- It can be applied by a rescuer or by the patient himself, as it allows easy adjustment of the required pressure
- It is also effectively applied over clothing

## Adjustable Blood Collection Band

Ref: **0807103**



- Material: Nylon & polyester
- Dimensions (LxW): 2.5 x 45cm
- Net Weight: 23gr
- The blood collection belt is quickly placed with one hand and securely locked
- With adjustable tightness & fixation clip
- Made of durable & comfortable strap
- Suitable for use in various medical settings, such as hospitals, clinics and blood donation centers

## Flexible Roll Splint

Ref: **0807106**



- Material: EVA Foam & Aluminum
- Dimensions (LxW): 100 x 11cm
- Net weight: 162 gr
- Small in size and weight, it can be placed in any first aid kit
- Reusable and waterproof
- Can be easily washed and disinfected with common cleaning and disinfection products
- It is placed quickly and easily, using a bandage or tape, without requiring any other auxiliary equipment
- Requires minimal training to use, is simple to apply and effective
- Does not contain latex

## First Aid Scissors

Ref: **0807107**



- Material: Ανοξειδωτο Ατσάλι και PVC
- Length: 19cm
- Net weight: 50 gr
- Easily and accurately cuts bandages
- With finger-shaped plastic handle for maximum comfort
- With serrated blade
- Blunt edge to minimize the risk of patient injury
- Round blade for easy cutting

## Body Bag

Ref: **0807105**



- Ideal for safe transport of corpses.
- Lightweight and durable construction with high resistance to tearing, puncturing, and leakage.
- Zipper along the entire length of the bag.
- 6 Handles for easy transport.
- Transparent pocket for placing the deceased's details, fixed on the zip.
- Material: Waterproof double-layer canvas.
- Color: Black
- Dimensions: 200 x 50 x 25 cm
- Maximum User Weight 136 kg
- Net Weight: 1.4 Kg

**CAUTION: For safe transport, the bag should be lifted simultaneously by all handles.**

## Bed Pan

Ref: **0810224**



- Length: 37cm
- Width: 23.8cm
- Height: 8.5cm
- Material: PP
- Net Weight: 110gr
- Capacity: 1L

## ABENA Bedpan Liner 20pcs (PACK)

Ref: **0809860**



- Disposable
- Absorbing Capacity  
660 ml of liquid

## ABENA Commode Liner 20pcs (PACK)

Ref: **0809861**



- Disposable
- Absorbing Capacity  
660 ml of liquid







*Thank you for your trust!*

