

ORTHOPEDICS CATALOGUE



CONTACT DETAILS

Headquarters:

Kathiana Akrotiri, Chania, Crete P.C. 731 00

T.: +30 28210 63222

E: medicalexports@mobiak.com

www.mobiakcare.com

Retail Shop, PYROMEDICAL Ltd:

Markou Botsari 96-98, Chania, Crete 731 36

DISTRIBUTION CENTERS

Central, South Greece & Ins. Greece:

Location: Schineza, Aspropyrgos, Attica P.C. 193 00

T: +30 210 4614006, 210 5593004

North Greece & Balkan Area:

Arkadiou 8, Kalochori, Thessaloniki P.C. 570 09

T: +30 2310 785668, 2310 778477

Central & Eastern Crete:

P Road, Industrial Zone of Heraklion, Crete P.C. 716 01

T: +30 2810 333151

Chania & Western Crete:

Kathiana Akrotiri, Chania, Crete P.C. 731 00

T: +30 28210 63222





contents

neck

arm

thoracolumbar

leg

footcare

kinesiology

cold-hot therapy



neck



SOFT CERVICAL COLLAR

- Made of medium density foam, radiolucent, hypoallergenic, antiperspirant with cotton lining

Indications:

- Neck injuries
- Cervical syndrome
- Neck muscle injuries

REF	H2421	H2422	H2423
NECK HEIGHT	6cm	8cm	10cm
SIZE/ PERIPHERY NECK		S 32 - 38cm	M 38 - 44cm L 44 - 50cm





PHILADELPHIA COLLAR



- Made of plastazot material ventilated and hypoallergenic
- It is divided into two sections that have Velcro parts for easy application
- The front part has a tracheal hole and support for the lower jaw to offer better stabilization

Indications:

- Severe neck injuries
- Traumatic treatment in neck injuries
- Occipital condyle fracture
- Herniated intervertebral cervical disc

REF	SIZE/ CIRCUMFERENCE NECK
	S 25 - 33cm
H2419	M 33 - 41cm
	L 41 - 49cm





ADJUSTABLE NECK BRACE

- It consists of two plastic sections, radiolucent, covered with foam, self-adhesive and washable linings
- In the front part it has a tracheal hole, lower jaw immobilization support and neck height adjustment system

Indications:

- Intervertebral disc herniation
- Rheumatoid arthritis neck
- Traumatic neck injuries



REF

SIZE/
CIRCUMFERENCE NECK

H2420

ONE SIZE (adjustable)





arm



SHOULDER BRACE

Indications:

- Intended for Shoulder injuries, Shoulder dislocation

- Made of neoprene fabric for compression and maintaining warmth in the shoulder area
- With transverse support strap with Velcro for easy application
- Ideal for re-joining sports after injury
- Can be used on either left or right hand (ambidextrous)

REF	SIZING / HUMERUS CIRCUMFERENCE
	S 34 - 36cm
	M 36 - 38cm
A2413	L 38 - 40cm
	XL 40 - 42cm
	XXL 42 - 44cm





DOUBLE ARM STRAP

- With adjustable cotton strap
- Easy to use
- Ambidextrous

Indications:

- Conservative rehabilitation of shoulder, arm, elbow injuries
- Aid for patients with a plaster cast

COMING SOON

REF	SIZE	STRAP LENGTH	STRAP WIDTH
A2462	ONE SIZE	119cm	3cm





BREATHABLE ARM SLING

- Made of cotton
- With adjustable hanging strap and protective shoulder pad

Indications:

- Postoperative use
- Conservative rehabilitation in shoulder, arm, elbow injuries
- Aid to a patient with a plaster cast

REF	FOREARM LENGTH
	S 32 - 36cm
A2434	M 37 - 41cm
	L 42 - 46cm





ARM SLING

- Made of cotton
- With adjustable hanging strap and protective shoulder pad

Indications:

- Postoperative use
- Conservative rehabilitation in shoulder, arm, elbow injuries
- Aid to a patient with a plaster cast

REF	FOREARM LENGTH
A2433	XS 28 - 33cm
	S 33 - 38cm
	M 38 - 43cm
	L 43 - 48cm





ARM SLING

- Arm support with elastic shoulder support for more comfort
- Height Adjustment with Buckle, offering different positions
- Lightweight, easy to install and very comfortable
- Can be used on either left or right hand (ambidextrous)

Indications:

- Providing arm support in a fixed position to reduce pain and discomfort
- Post-operative use
- Conservative rehabilitation in shoulder, arm, elbow injuries
- Aid to a patient with a plaster cast

REF	SIZE/ CIRCUMFERENCE ARM
0806935	ONE SIZE (adjustable)



SHOULDER IMMOBILIZER

- Made of hypoallergenic, ventilated foam
- It has a perimeter chest strap and a hand suspension strap with Velcro for easy length adjustment and fit
- Sets the elbow at 90 degrees firmly immobilized

Indications:

- Shoulder dislocation
- Sub cephalic humerus fracture
- Postoperative use

REF

SIZING / HUMERUS
CIRCUMFERENCE

A2432

ONE SIZE (adjustable)





SHOULDER IMMOBILIZER VELPEAU TYPE



- Made of breathable fabric
- Special design for shoulder-arm immobilization
- Easy to use
- Ambidextrous

Indications:

- Shoulder dislocation
- Subcranial humeral dislocation
- Fixed non-displaced scapula dislocation

COMING SOON

REF	SIZE/ CIRCUMFERENCE*
A2467	S/M 70 - 105cm
	L/XL 105 - 135cm

*Under-bust measurement





PADDED ARMSLING 15/30 DEGREES



- Construction of high density foam covered with hypoallergenic fabric
- It has a divided pillow with an anatomical curve for proper fit
- It has fabric arm sling and fabric immobilization straps
- Can be used on either left or right hand (ambidextrous)

Indications:

Preservative in partial rupture of the rotator cuff, post-operative for healing of irreparable fractures, small-medium tension rupture of the rotator cuff.



REF

SIZING / HUMERUS
CIRCUMFERENCE

A2445

ONE SIZE (adjustable)



PADDED ARMSLING 30 DEGREES

- Provides abductor immobilization of the shoulder for pain reduction and recovery
- Made of high-density foam covered with hypoallergenic fabric
- Includes a sleeve and stabilization strap
- Ambidextrous

Indications:

- Conservative in partial rupture of rotator cuff tendons
- Postoperative treatment of non-displaced fractures
- Small or medium-sized suturing of rotator cuff tendons

COMING SOON

REF	SIZE	SPLINT LENGTH	SPLINT WIDTH
A2459	ONE SIZE	38cm	17cm





SHOULDER – ARM ABDUCTION IMMOBILIZER 45/60 DEGREES WITH SLING

- Construction of high density foam covered with hypoallergenic fabric
- The pillow has an anatomical curve for proper fit
- It has fabric arm sling and fabric immobilization straps
- Can be used on either left or right hand (ambidextrous)

Indications:

- After surgeries on the tendinous shoulder rotator cuff
- Postoperative immobilization after implantation of prosthetic shoulders
- In case of impingement syndrome (shoulder tendonitis)
- In case of acute and chronic folliculitis
- In habitual dislocations and subluxations

REF

SIZING / HUMERUS
CIRCUMFERENCE

A2444 **ONE SIZE** (adjustable)





HUMEROUS FRACTURE BRACE

- Made of polyethylene
- It has an inner lining made of plastazot, which makes it extremely light and flexible
- It has stabilization straps for equal distribution of pressure on the humerus

Indications:

- Humerus diaphysis fracture

REF	SIZE/ CIRCUMFERENCE ARM
	S 23 – 30cm
A2435	M 30 – 37cm
	L 37 – 44cm





ELASTIC ELBOW SUPPORT WITH STRAP



- Made of high-strength elastic fabric with 3D knitting
- Features a stabilization strap
- Provides gentle, gradual compression to tissues and soft tissues
- Suitable for exercise and daily activities
- Ambidextrous

Indications:

- Elbow epicondylitis recovery
- Elbow ligament-tendon injuries

COMING SOON

REF	ELBOW CIRCUMFERENCE	HEIGHT
A2477	S 18 - 25cm	24cm
	M 25 - 30cm	24cm
	L 30 - 35cm	24cm
	XL 35 - 40cm	24cm
	XXL 40 - 45cm	24cm



AEROPRENE ELBOW SUPPORT



- Made of Aeroprene fabric
- Easy to use
- Suitable for sports

Indications:

- Tendon-ligament injuries
- Ideal for sports reintegration

COMING SOON

REF	SIZE	HEIGHT
A2505	ONE SIZE	18.5cm





EPICONDYLITIS ELBOW STRAP



Indications:

- Elbow epicondylitis recovery

- Features a high quality silicone pad
- With adjustable velcro strap
- Ambidextrous

COMING SOON

REF	SIZE/ ELBOW CIRCUMFERENCE
A2472	S/M 30 - 40cm
	L/XL 40 - 50cm





CUBITAL ELBOW BRACE

- Made of hypoallergenic fabric
- Places the elbow joint in a resting position, preventing its involuntary bending
- Lightweight
- Ambidextrous

Indications:

- Conservative in ulnar neuritis
- Elbow fracture

COMING SOON

REF SIZE

A2468 ONE SIZE





ELBOW BRACE ROM

- Adjustable Range 34 to 43.5cm
- Equipped with mechanism for Immobilization and Adjustment of the Range of Motion, Graduated in Degrees (Range Of Motion)

REF

0808306 Left

0808307 Right





TELESCOPIC ROM ELBOW BRACE



- Made of lightweight metal alloy. It has telescoping system to cover all lengths
- Easy adjustment system from 0–90 degrees extension and 0–120 degrees flexion
- It has quick lock system and anatomically shaped pads internally, hypoallergenic, detachable & washable

Indications:

- Conservatives in radius head fracture
- Medial or lateral condyle fractures
- Reversed elbow dislocation
- Tendon ligament rupture
- Postoperative transcondylar and supracondylar fractures
- Biceps fixation
- Medial collateral ligament rupture
- Compression fractures with unstable pulley
- Radius head fracture
- Olecranal fracture

REF	SIZE/ CIRCUMFERENCE ARM
A2453	ONE SIZE
Left	
Right	



ELASTIC WRIST BAND

- Provides support and stability to the wrist through gentle, effective compression
- Features a velcro adhesive for better stabilization & retention
- Elastic & comfortable to use during sports and daily activities

Indications:

- Mild wrist injuries

COMING SOON

REF	SIZE	HEIGHT
A2478	ONE SIZE	7.7cm





ELASTIC WRIST BAND REINFORCED

- Provides support and stability of the wrist through gentle effective compression
- Features a velcro adhesive for better stabilization
- Elastic & comfortable to use as it rests and soothes pain, maintaining the range of motion of the wrist
- Ambidextrous

Indications:

- Mild wrist injuries

COMING SOON

REF	SIZE	HEIGHT
A2480	ONE SIZE	7.7cm



WRIST SUPPORT NEOPRENE



- Made of neoprene fabric to prevent moisture
- With Velcro for better fixation
- Ideal for reintegration and use during sports
- Ambidextrous

Indications:

- Minor wrist injuries
- Return to sports after injury

COMING SOON

REF	SIZE	HEIGHT
A2463	ONE SIZE	14cm





WRIST SUPPORT WITH THUMB NEOPRENE

- Made of neoprene fabric to prevent moisture
- Equipped with Velcro for better fixation
- Ideal for reintegration and use during sports
- Ambidextrous

Indications:

- Minor wrist-thumb injuries
- Return to sports after injury

COMING SOON

REF	SIZE	HEIGHT
A2511	ONE SIZE	14cm





THUMB SPLINT UNIVERSAL

- Made of hypoallergenic, breathable fabric
- With metal immobilization bar support along the thumb area
- Easy application with two Velcro straps



Indications:

- Thumb injuries
- De Quervain syndrome
- Thumb arthritis

REF	SIZE	LENGTH	WIDTH
A2454	ONE SIZE	40.5cm	5cm



WRIST BRACE

- Made of hypoallergenic, anti-sweat, non-stretch fabric
- With aluminum alloy palm panel
- Easy installation and application with 3 straps and Velcro

Indications:

- Carpal tunnel syndrome
- Wrist injuries and contusions
- Postoperative and conservative rehabilitation in colles fracture



REF		WRIST CIRCUMFERENCE	
A2440	Right	S	14 - 17cm
		M	18 - 21cm
A2450	Left	L	22 - 25cm



WRIST BRACE UNIVERSAL

- Made of hypoallergenic, anti-sweat, non-stretch fabric
- With aluminum alloy palm panel
- Easy installation and application with 3 straps and Velcro

Indications:

- Carpal tunnel syndrome
- Wrist injuries and contusions
- Postoperative and conservative rehabilitation in colles fracture



REF		SIZE	LENGTH	WIDTH
A2451	Right Left	ONE SIZE	20cm	10.5cm



FOREARM SPLINT AMBIDEXTROUS 25CM

- Features 2 metal bars that prevent dorsiflexion and palmar flexion
- Provides wrist support to reduce pain and aid recovery
- Easy to apply, made of hypoallergenic & antiperspirant fabric
- Ambidextrous

Indications:

- Carpal Tunnel Syndrome
- Wrist Injuries and Sprains
- Postoperative and Conservative Rehabilitation of Colles Fracture

COMING SOON

REF	SIZE	LENGTH
A2455	ONE SIZE	25cm





THUMB-WRIST BRACE SPICA

- Made of hypoallergenic, anti-sweat, non-stretch fabric, 20cm length
- With aluminum alloy palm panel and thumb immobilization panel
- Easy installation and application with 3 straps, Velcro and thumb immobilization strap

Indications:

- Thumb injuries
- Scaphoid fracture
- De Quervain syndrome
- Postoperative thumb tendon surgery

REF		WRIST CIRCUMFERENCE
A2437	Right	S 18 - 25cm
		M 25 - 30cm
A2438	Left	L 30 - 35cm





THUMB-WRIST BRACE SPICA UNIVERSAL



- Made of hypoallergenic, anti-sweat, non-stretch fabric, 20cm length
- With aluminum alloy palm panel and thumb immobilization panel
- Easy installation and application with 3 straps, Velcro and thumb immobilization strap

Indications:

- Thumb injuries
- Scaphoid fracture
- De Quervain syndrome
- Postoperative thumb tendon surgery

REF	SIZE	LENGTH	WIDTH
A2439	Right Left	ONE SIZE	20cm 11.5cm





RESTRAINER (PAIR) MOBIAK



- Made of lightweight and breathable material with an internal anatomical, foam lining
- Ability to adjust the length via an adjustable strap & velcro closure
- Available as a pair
- Ambidextrous

Indications:

- Immobilization of patients' limbs

COMING SOON

REF	SIZE	LENGTH	STRAP WIDTH
A2479	ONE SIZE	28.5cm	7.5cm





FINGER SPLINT OVAL 8



- Protects and supports the smaller joints of the fingers
- Comfortable & lightweight for everyday use
- Can be used without restriction in the bathroom and the sea

Indications:

- Used only by patients with orthopedic problems
- **Conservative treatment in:** Fractures, Rheumatoid arthritis, Dismorphism or hyperextension, boutonnière deformity

COMING SOON

REF	CIRCUMFERENCE (MEASURED IN THE MIDDLE SEMI-CIRCLE)			
	No1: 4.3cm	No4: 5.1cm	No7: 6cm	No10: 6.9cm
A2461	No2: 4.6cm	No5: 5.4cm	No8: 6.3cm	No11: 7.2cm
	No3: 4.9cm	No6: 5.7cm	No9: 6.6cm	No12: 7.5cm



FINGER SPLINT FOLD OVER

- Made of aluminum
- Lined with foam which makes it very comfortable for the patient
- It has Velcro straps for its immediate application

Indications:

- Finger ligament and tendon injuries
- Injuries and fracture of phalangophalangeal and interphalangeal joint

REF	SIZE/ FINGER LENGTH
	S 4 – 5cm
A2452	M 5 – 6.5cm
	L 6.5 – 8cm



FINGER SPLINT TYPE BOUTONNIERE

- Made of high quality plastic, stainless steel & soft EVA inner lining
- Immobilizes the finger while allowing active flexion
- Easy to use

Indications:

- Conservative treatment of finger fractures
- Rheumatoid arthritis
- Ligamentous injuries of fingers
- Boutonniere deformity

COMING SOON

REF	SIZE/ FINGER LENGTH
-----	------------------------

S: 5cm

A2471

M: 7.5cm

L: 10cm





WATERPROOF ARM COVER (ADULT)

REF	LENGTH
-----	--------

0808511	56cm
----------------	------

- Easy to apply
- Reusable waterproof PVC covering with plastic and elastic ring
- Convenient for use in shower or bath
- Intended for single patient use



WATERPROOF ARM COVER (CHILD)

REF	LENGTH
-----	--------

0808512	45.72cm
----------------	---------





**thoraco-
lumber**

thoracolumbar



POST- OPERATIVE BELT 16CM

- Made of high-strength elastic fabric, breathable and hypoallergenic
- Perfect fit to the body
- Offers light compression, useful for mild pain or long-term treatment

Indications:

- Mild to moderate back pain
- Supports the abdominal area after surgeries
- Prevents sagging in the abdominal area after pregnancy

COMING SOON

REF	SIZE/CIRCUMFERENCE	HEIGHT
W2456	S 70 - 90cm	16cm
	M 90 - 110cm	16cm
	L 110 - 130cm	16cm





POST- OPERATIVE BELT 20CM



- Made of high-strength elastic fabric, breathable and hypoallergenic
- Perfect fit to the body
- Offers light compression, useful for mild pain or long-term treatment

Indications:

- Mild to moderate back pain
- Supports the abdominal area after surgeries
- Prevents sagging in the abdominal area after pregnancy

COMING SOON

REF	SIZE/CIRCUMFERENCE	HEIGHT
	S 70 - 90cm	20cm
W2457	M 90 - 110cm	20cm
	L 110 - 130cm	20cm





POST - OPERATIVE ELASTIC SUPPORT

- Made of high resistance, elastic ventilated and hypoallergenic fabric
- With large divided Velcro for better application
- Soft lining in the abdominal area for greater comfort

Indications:

- Mild or moderate lower back pain
- Supports abdominal area following abdominal surgeries
- Prevents sagging in the abdominal area following pregnancy

REF	SIZE/ CIRCUMFERENCE	HEIGHT
	S 70 - 90cm	26cm
W2408	M 90 - 110cm	26cm
	L 110 - 130cm	26cm





POST-OPERATIVE ABDOMINAL (UMBILICAL) SUPPORT



- Made of high resistance, elastic ventilated and hypoallergenic fabric
- With Velcro closure
- With anatomical cushion that fits the area of the umbilical hernia

Indications:

Supports the abdominal area around the hernia. Abdominal hernias, including postoperative, umbilical, periumbilical hernias, lumbar hernias. Prevention of compression of severe hernias during intense physical exertion.

REF	SIZE/CIRCUMFERENCE	HEIGHT	
W2443	S	80 - 90cm	16cm
	M	90 - 100cm	16cm
	L	100 - 110cm	16cm
	XL	110 - 120cm	16cm
	XXL	120 - 135cm	16cm





UMBILICAL HERNIA BELT

- Made of high resistance, elastic ventilated and hypoallergenic fabric.
- It's extremely small size makes it very discreet and comfortable for the patient.

Indications:

Conservative treatment
of umbilical hernia

REF	SIZE/CIRCUMFERENCE
W2448	ONE SIZE (adjustable)





HERNIA BELT SINGLE MOBIAK

- Easy and controlled application with velcro
- With anatomical pressure pad for the inguinal hernia area

Indications:

- Conservative Treatment of Inguinal Hernia

COMING SOON

REF

SIZE/PELVIC
CIRCUMFERENCE

W2474 Right

S: 80 – 85cm

M: 85 – 95cm

W2473 Left

L: 95 – 105cm



thoracolumbar



AEROPRENE BACK SUPPORT



- Made of Aeroprene fabric
- Easy to use
- Suitable for sports

Indications:

- Return to sports after injuries
- Sciatica, back pain

COMING SOON

REF	SIZE	BACK HEIGHT	FRONT HEIGHT
W2504	ONE SIZE	19cm	13.5cm





LUMBAR SUPPORT (LOW)

- Made of high-strength elastic fabric, breathable and hypoallergenic
- Perfect fit to the body
- Offers restriction of lumbar movements to reduce pain and discomfort
- Suitable for use during work

Indications:

- Back pain
- Sciatica
- Spondyloarthropathy
- Muscle-ligament injuries
- Intervertebral disc herniation

COMING SOON

REF	SIZE/CIRCUMFERENCE	BACK HEIGHT	FRONT HEIGHT
W2458	S/M: 80 - 100cm	16.2cm	13.5cm
	L/XL: 100 - 130cm	16.2cm	13.5cm





LUMBOSACRAL CORSET

- Made of high resistance, elastic ventilated and hypoallergenic fabric
- With four plastic panels in the back part
- With Velcro closure in the front part for easy application

Indications:

Back pain, Sciatica,
Spondyloarthropathy,
Muscle-ligament injuries,
Intervertebral disc herniation.



REF	SIZE/CIRCUMFERENCE
W2447	S 80 – 90cm
	M 90 – 100cm
	L 100 – 110cm
	XL 110 – 120cm
	XXL 120 – 135cm
FRONT HEIGHT	13cm
BACK HEIGHT	26cm



LUMBAR SUPPORT WITH STAYS AND TENSIONER

- Made of high resistance, elastic ventilated and hypoallergenic fabric
- Easy closure with Velcro and extra straps offering greater stability
- With four plastic panels and removable support pad

Indications:

Back pain, Sciatica,
Spondyloarthropathy,
Muscle-ligament injuries,
Intervertebral disc herniation.



REF	SIZE/CIRCUMFERENCE
W2449	S 80 – 90cm
	M 90 – 100cm
	L 100 – 110cm
	XL 110 – 120cm
	XXL 120 – 135cm
FRONT HEIGHT	13cm
BACK HEIGHT	26cm

ELASTIC 3D KNITTED LUMBAR SUPPORT WITH STAYS



- Made of high resistance, elastic ventilated and hypoallergenic fabric
- With four plastic panels in the back part
- With Velcro closure in the front part for easy application

Indications:

Mild or moderate back pain, mild post-traumatic injuries, prophylaxis in sports activities and muscle weakness, sciatica, humeral spondyloarthropathy, muscle-ligament injuries, lumbar intervertebral disc herniation.

REF	SIZE/CIRCUMFERENCE
W2442	S 80 - 90cm
	M 90 - 100cm
	L 100 - 110cm
	XL 110 - 120cm
	XXL 120 - 135cm
FRONT HEIGHT	17cm
BACK HEIGHT	27cm



LUMBAR SUPPORT WITH STAYS AND TENSIONER (32cm)

- Made of high resistance, elastic ventilated and hypoallergenic fabric
- Easy closure with Velcro and extra straps offering greater stability to the lumbar area
- With 2 metal and 2 plastic panels at the back

Indications:

Supports the back muscles in low back problems (Lumbago) in combination with zygapophyseal joint disorders. Severe back pain, sciatica, degenerative processes, weakness and muscle weakness, pre- and post-operative treatments, herniated disc, arthritis of the spine, arthritis of the spine and spondylolisthesis of the lumbar region, osteomalacia.

REF	SIZE/ CIRCUMFERENCE
	S 70 - 90cm
	M 90 - 110cm
W2446	L 110 - 130cm
	XL 110 - 130cm
	XXL 110 - 130cm
FRONT HEIGHT	16cm
BACK HEIGHT	32cm



LUMBAR SUPPORT BRACE WITH LACING COMPRESSION SYSTEM

- Made of non-stretch ventilated hypoallergenic fabric
- System with compression cords for maximum pressure distribution

Indications:

- Back pain and sciatica of the lumbar spine
- Lumbar spine intervertebral disc herniation
- Spondyloarthropathy of the lumbar spine
- Traumatic injuries of muscles and ligaments of the lower thoracic spine or lumbar spine

REF	SIZE	CIRCUMFERENCE
W2401	ONE SIZE (adjustable)	110cm - 165cm (including the extension straps 25cm)
FRONT HEIGHT		17cm
BACK HEIGHT		22.5cm



LUMBAR SUPPORT BRACE WITH LACING COMPRESSION SYSTEM



- Made of non-elastic breathable fabric
- Features a compression cord system
- Perfect fit and support

Indications:

- Back pain and sciatica of the lumbar spine
- Herniated intervertebral disc of the lumbar spine
- Spondyloarthropathy of the lumbar spine
- Traumatic injuries of muscles and ligaments of the lower thoracic spine or lumbar spine
- Spondylolisthesis

COMING SOON

REF	SIZE/CIRCUMFERENCE
W2460	S/M: 70 – 90cm
	L/XL: 90 – 110cm
	XXL/XXXL: 110 – 140cm
FRONT HEIGHT	14.7cm
BACK HEIGHT	26cm

thoracolumbar



POSTURE HARNESS

Indications:

Initial stages of spinal kyphosis, Improving spinal posture.

- Made of high resistance, elastic ventilated and hypoallergenic fabric
- With two plastic panels in the back part
- Belt and shoulder straps close at the front part for greater patient comfort

REF	SIZE/CIRCUMFERENCE
-----	--------------------

S	80 - 90cm
----------	-----------

M	90 - 100cm
----------	------------

T2407	L	100 - 110cm
--------------	----------	-------------

XL	110 - 120cm
-----------	-------------

XXL	120 - 135cm
------------	-------------





POSTURE CORRECTOR



- Made of breathable, hypoallergenic fabric
- Easy to use

Indications:

- Spinal kyphosis

COMING SOON

REF	SIZE	FRONT HEIGHT	BACK HEIGHT
W2509	ONE SIZE	11cm	35cm

thoracolumbar



UNIVERSAL THORACOLUMBAR SPINE SUPPORT BRACE WITH LACING COMPRESSION SYSTEM

- Construction of non-elastic ventilated hypoallergenic fabric
- System with compression cords that help in the correct application of pressure exerted by the splint
- On the back part, it has a plastic panel, detachable, lined with foam

Indications:

- Back pain and sciatica of the lumbar spine
- Lumbar spine intervertebral disc herniation
- Spondyloarthropathy of the lumbar spine
- Traumatic injuries of muscles and ligaments of the lower thoracic spine or lumbar spine
- Spondylolisthesis

REF	SIZE	CIRCUMFERENCE
T2402	ONE SIZE (adjustable)	80cm - 132cm
FRONT HEIGHT		16.5cm
BACK HEIGHT		22.5cm







SPINAL ORTHOSIS TAYLOR TYPE

- Made of high resistance, elastic ventilated & hypoallergenic fabric
- Easy closure with Velcro and extra straps that offer greater stability to the lumbar area
- It has 2 flexible metal panels at the back that cover the entire spine area
- It has two straps for lifting support of the thoracic part of the spine which close with Velcro in the front part

Indications:

It helps to correct the dorsal spine and helps to treat kyphosis. Provides stabilization for hyperextension in vertebral fractures of the lumbar and thoracic spine. It is also recommended for any trauma to the thoracic and lumbar vertebrae and in the early stages of osteoporosis.

REF	SIZE/CIRCUMFERENCE
T2405	S 80 - 90cm
	M 90 - 100cm
	L 100 - 110cm
	XL 110 - 120cm
	XXL 120 - 135cm



SPINAL ORTHOSIS TAYLOR TYPE (Y-SHAPE)

- Made of aluminum alloy lined with breathable, hypoallergenic, cotton fabric.
- It has T-shape metal panels on the back.
- The front part closes like a belt with Velcro and extra Velcro straps for maximum stability.
- The upper part has underarm braces which place the patient's torso in the correct position.

Indications:

It helps to correct the dorsal spine and helps to treat kyphosis. Provides stabilization for hyperextension in vertebral fractures of the lumbar and thoracic spine. It is also recommended for any trauma to the thoracic and lumbar vertebrae and in the early stages of osteoporosis.

REF	SIZE/CIRCUMFERENCE	FRONT HEIGHT	BACK HEIGHT	
T2406	S	80 - 90cm	51cm	
	M	90 - 100cm	52cm	
	L	100 - 110cm	21.5cm	53cm
	XL	110 - 120cm		54cm
	XXL	120 - 135cm		55cm





TALL SPINAL SUPPORT "SPINE ORTHOSIS"

- Construction of aluminum alloy frame
- Frame easily configurable, integrated in hypoallergenic antiperspirant fabric that offers easy closure at the front like a lumbar belt
- It also has torso stabilization straps in the upper part

Indications:

- Osteoporosis
- Osteoporotic fractures
- Kyphosis spine
- Prevention of osteoporotic fractures
- SCHEURMANN disease

REF	SIZE/CIRCUMFERENCE
-----	--------------------

T2424	S 85 - 100cm
	M 100 - 115cm
	L 115 - 130cm
	XL 130 - 145cm





THORACOLUMBAR BRACE WITH LACING COMPRESSION SYSTEM AND SHOULDER STRAPS



- Construction of non-elastic ventilated hypoallergenic fabric
- System with compression cords that help in the correct application of pressure exerted by the splint
- In the rear part the height of the splint is easily adjusted with Velcro
- The splint leaves the upper abdomen part free

Indications:

- Type I osteoporosis of the spine
- Prevention of Osteoporotic Vertebral Fractures type II
- Spondyloarthropathies of thoracic lumbar spine
- Scheuermann's disease
- Kyphosis



REF	SIZE	CIRCUMFERENCE
T2403	ONE SIZE (adjustable)	80cm - 132cm
FRONT HEIGHT		16.5cm
BACK HEIGHT		45cm

thoracolumbar



3 POINT HYPER-EXTENSION BRACE

- Made of lightweight aluminum alloy
- Lined with hypoallergenic, ventilated fabric at the support points. Modular sternum support pad, lined with fabric
- Rear plastic support pad, detachable, to be placed in a vertical or horizontal position
- Locking mechanism with buckle for easy application and safety
- Horizontal aluminum bar in the pelvis area that allows sitting and semi-sitting position
- Height adjustment on the sides of the splint

Indications:

Vertebral Compression Fractures. Osteoporosis. Dislocation. Spinal Arthritis, Compression, Vertebral Fracture T10-L2, Burst Vertebral Fracture T10-L2, Rotation Vertebral Fracture T10-L2, Osteoporosis Vertebral Fracture T10-L2 Spinal Fusion, Laminectomy (Thoracolumbar), Discectomy (Thoracolumbar).

REF	SIZE/CIRCUMFERENCE	HEIGHT
T2404	S 60 - 75cm	39 - 46cm
	M 75 - 90cm	42 - 49cm
	L 90 - 105cm	45 - 52cm
	XL 105 - 115cm	49 - 56cm







leg



CALF SUPPORT



- With velcro fixation straps for uniform compression on soft tissues
- Ambidextrous

Indications:

- Conservative rehabilitation of injuries
- Gastrocnemius muscle strains
- Preventive use to avoid injuries, gastrocnemius muscle strains

COMING SOON

REF SIZE

L2469 ONE SIZE





THIGH SUPPORT



- Features velcro fixation straps for uniform compression on soft tissues
- Ideal for athletic activity
- Ambidextrous

Indications:

- Conservative rehabilitation of injuries
- Thigh muscle strains,
- Preventive use to avoid injuries, thigh muscle strains

COMING SOON

REF SIZE

L2470 ONE SIZE





ELASTIC KNEE SUPPORT

Indications:

- Mild lateral ligament injuries
- Mild knee injuries and sprains
- Mild knee conditions such as swelling, inflammation, tendonitis, arthritis

- Made from high-strength elastic fabric with 3D knitting
- Provides gradual compression to soft tissues and muscles, increasing blood flow
- Ideal for use while working, exercising or walking
- Ambidextrous

COMING SOON

REF	SIZE/ KNEE CIRCUMFERENCE	HEIGHT
L2481	S 37 - 43cm	25.5cm
	M 43 - 49cm	25.5cm
	L 49 - 55cm	25.5cm
	XL 55 - 61cm	25.5cm
	XXL 61 - 67cm	25.5cm





NEOPRENE KNEE SUPPORT



- Made of neoprene fabric to prevent moisture
- Features an integrated foam ring to stabilize the patella and adhesive straps for gradual, uniform compression of soft tissues and muscles
- Ambidextrous

Indications:

- Patellar chondropathy
- Mild knee instability
- Return to sport after injury

COMING SOON

REF	SIZE	HEIGHT
L2475	ONE SIZE	27cm





PATELLA STRAP WITH SILICONE

- Made of neoprene fabric with an interior silicone pad that fits the shinbone
- With Velcro closure for easy application
- Ideal for sports activities

Indications:

- Recovery in knee ligament tendinitis (Jumper's knee)
- Osteochondritis (Osgood Schlatter)

REF

SIZE

L2429 ONE SIZE





ELASTIC KNEE SUPPORT WITH SILICONE RING AND STAYS

Indications:

- Patellar chondropathy
- Mild knee instability
- Rehabilitation of mild lateral ligament injuries

- Made of elastic, anti-sweat and breathable fabric
- Features flexible panels
- Features a silicone pad-ring in the patella area
- Provides uniform compression to tissues and soft tissues
- Ambidextrous

COMING SOON

REF	SIZE/ KNEE CIRCUMFERENCE	HEIGHT	
L2476	S	30 - 34cm	30cm
	M	34 - 38cm	30cm
	L	38 - 42cm	30cm
	XL	42 - 46cm	30cm
	XXL	46 - 50cm	30cm





NEOPRENE KNEE SUPPORT WITH STAYS SUPPORT AND STRAPS

- Made of neoprene fabric to prevent moisture
- It has six spiral stainless steel bars
- It has a special ring to stabilize and limit the mobility of the patella and two adhesive stabilization straps with velcro
- Suitable for exercise and daily activities
- Ambidextrous

Indications:

- Knee-Patella Injuries
- Patellar Chondropathy
- Lateral Ligament Injuries

COMING SOON

REF	SIZE/ KNEE CIRCUMFERENCE	HEIGHT
L2466	S/M 34 - 39cm	26cm
	L/XL 39 - 46cm	26cm





NEOPRENE KNEE SUPPORT WITH STAYS SUPPORT

Indications:

- Lateral ligament injuries
- Patellar chondropathy
- Mild osteoarthritis of the knee
- Knee joint instability

- Made of neoprene fabric to prevent moisture
- Features four spiral stainless steel rods
- Features a special ring to stabilize and limit the mobility of the patella and 2 adhesive stabilization straps with velcro
- Suitable for exercise and daily activities
- Ambidextrous

COMING SOON

REF	SIZE/ KNEE CIRCUMFERENCE	HEIGHT
L2465	S 32 - 34cm	31cm
	M 34 - 37cm	31cm
	L 37 - 40cm	31cm
	XL 40 - 43cm	31cm
	XXL 43 - 46cm	31cm





NEOPRENE KNEE SUPPORT WITH HINGE



- Made of neoprene fabric to prevent moisture
- Features two articulated stainless steel bars
- Features a special ring to stabilize and limit the mobility of the patella and two adhesive stabilization straps with velcro
- Ambidextrous

Indications:

- Meniscus injuries
1st and 2nd degree
- Knee arthritis
- Progressive rheumatoid arthritis
- Lateral knee ligament injuries
- Anterior – posterior cruciate ligament injuries

COMING SOON

REF	SIZE/ KNEE CIRCUMFERENCE	HEIGHT
L2464	S 32 – 34cm	31cm
	M 34 – 37cm	31cm
	L 37 – 40cm	31cm
	XL 40 – 43cm	31cm
	XXL 43 – 46cm	31cm



TRI-PANEL KNEE BRACE MOBIAK

- Made of hypoallergenic foam divided into three sections (1 back and 2 on the sides) which have rigid aluminum panels
- The three sections are joined by Velcro to shape the circumference
- It has straps that immobilize the knee

Indications:

- Patella dislocation
- Conservative rehabilitation in patella fracture
- Anterior (acl) or posterior (pcl) cruciate ligament injury
- Postoperative rehabilitation

REF	SIZE	HEIGHT
L2426		40cm
L2427	ONE SIZE	50cm
L2428		60cm





ROM KNEE BRACE

- Made of lightweight metal material
- It has foam pads, non-slip and fabric stabilization straps that offer comfort and easy application
- It has two multicenter mechanisms (right and left) that have an adjustment range of 0 to 90 degrees in extension and 0 to 120 degrees in flexion

Indications:

- Conservative rehabilitation in stable tibial plateau fracture
- Conservative rehabilitation in stable tibial-femoral condyle fracture
- Postoperative use in total patellar tendon rupture
- Quadriceps rupture
- Displaced patella fracture
- Anterior posterior (pcl) and lateral

REF	SIZE	HEIGHT
L2430	ONE SIZE	40cm
L2431		50cm





FUNCTIONAL KNEE BRACE WITH ROM

- Made of lightweight metal material
- It has foam pads, non-slip and fabric stabilization straps that offer comfort and ease of application
- It has two multicenter mechanisms (right and left) that have an adjustment range from 0 to 90 degrees in extension and 0 to 120 degrees in bending and quick lock system
- The length is easily adjusted with a button

Indications:

Conservative treatment for:

- High condyle and tibial condyle non-displaced stable fractures
- Tibial plateau non-displaced stable fractures
- Non-displaced stable patellar fractures

Post-operative treatment for:

- Patella tendon rupture, ACL, PCL, MCL, LCL injuries
- High tibial osteotomy
- Quadriceps rupture, Osteochondral condyle fractures

REF	SIZE	EXTENTION	FLEXION	LENGTH
L2418	ONE SIZE	0° - 90°	0° - 120°	50 - 70cm



TELESCOPIC ROM KNEE BRACE, UNIVERSAL

- Adjustable Range 58.5 to 73cm
- Features an Immobilisation mechanism and Articulated Range of Motion

REF	SIZE
0807669	ONE SIZE (adjustable)





ELASTIC ANKLE SUPPORT

Indications:

- Mild ankle sprains
- Minor ankle ligament injuries

- Made from high-strength elastic fabric with 3D knitting
- Ideal for minor sprains and ligament problems
- Provides support and stability to the ankle through gentle effective compression
- Ambidextrous

COMING SOON

REF	SIZE/ ANKLE CIRCUMFERENCE	HEIGHT
L2482	S 18 - 21cm	25.5cm
	M 21 - 24cm	25.5cm
	L 24 - 27cm	25.5cm
	XL 27 - 30cm	25.5cm





AEROPRENE ANKLE SUPPORT



- Made of Aeroprene fabric
- Easy to use
- Suitable for sports

Indications:

- Mild ankle sprains
- Ligament injuries

COMING SOON

REF	SIZE	HEIGHT
-----	------	--------

L2506	ONE SIZE	19cm
--------------	-----------------	------





ELASTIC ANKLE SUPPORT WITH STRAP

Indications:

- Minor ankle ligament injuries
- Mild ankle sprains
- Rheumatoid arthritis of the ankle

- Made of high-strength elastic fabric with 3D knitting
- Features a bandage strap
- Ideal for mild sprains and ligament problems
- Provides support and stability to the ankle through gentle effective compression
- Ambidextrous

COMING SOON

REF	SIZE/ ANKLE CIRCUMFERENCE	HEIGHT
L2483	S 18 – 21cm	26cm
	M 21 – 24cm	26cm
	L 24 – 27cm	26cm
	XL 27 – 30cm	26cm





LACE UP ANKLE BRACE



- Made of perforated fabric reinforced with plastic lining and lace closure
- Can completely replace air-cast braces
- Comfortable due to its design and material
- Ideal for sports activities and daily life
- Can be used on either the left or the right leg (ambidextrous)

Indications:

- Hammer injuries
- Type I and II sprains
- Chronic ankle instability
- Preventive use

REF	SIZE	CIRCUMFERENCE OF HEEL-TARSO
	S	27 - 30cm
L2412	M	30 - 33cm
	L	33 - 36cm





AIR/GEL ANKLE BRACE

- Made of soft plastic material divided into two sections
- Pads with air (upper section) and gel (bottom section)
- Easy application with two straps and Velcro, attached to the side parts of the brace
- Can be used on either the left or the right leg (ambidextrous)

Indications:

- Ankle sprain
- Ankle ligament injuries

REF	SIZE
L2425	ONE SIZE





AIR CAST BRACE

- Two braces providing mediolateral ankle stabilisation holding and maintaining the ankle joint in a functional position offering better adjustment and comfort to the user
- Easy to adjust thanks to external Velcro straps
- The lateral sections are lined with an adjustable air cushion

REF	SIZE
0807789	ONE SIZE





ANKLE REST BRACE NIGHT – SPLINT

- Made of durable plastic lined with, ventilated, antiperspirant, hypoallergenic fabric
- It has three support straps with cushions for comfort
- Two side straps with Velcro to adjust (increase) the dorsal flexion
- The brace can be used on either the left or the right leg

Indications:

- Foot drop from any pathological cause
- Plantar fasciitis
- Achilles tendonitis
- Ankle immobilization

REF	SHOE SIZE
L2414	S No <38
	M No 38 – 42
	L No 42 – 46
	XL No 46+





DROP FOOT BRACE

Indications:

Conservative treatment for:

- Flaccid paresis – hypopodia
- Ankle drop
- Foot support

- Made of thermoplastic material
- Easy to apply with a velcro

COMING SOON

REF

L2507 Right

L2508 Left

SIZE ΠΑΠΟΥΤΣΙΟΥ

S/M

L/XL

34 – 42

42 – 48





CHILD WALKER



- Made of plastic material
- It has an inner lining of antiperspirant, hypoallergenic fabric
- It has plastic side panels, for support and retention of the ankle, which have straps with Velcro for easy application
- Non-slip sole, with rocker bottom system for the correct walking of the patient
- The brace can be used on either the left or the right leg

Indications:

Conservative in stable non-displaced fibular fracture, non-displaced medial malleolus fracture, metatarsal fracture, stable non-displaced bifibular fracture, tibiofibular ligament rupture, postoperative in fibula fracture, medial fracture, medial malleolus fracture, bifibular fracture, ligament-tendon sutures.

REF	SHOE SIZE
	S No <25
L2416	M No 25 - 29
	L No 29 - 33



STD FRACTURE WALKER SHORT

REF	SIZING		
0806941	S	Sole 25cm	No 36 – 39
0806942	M	Sole 29cm	No 40 – 43
0806943	L	Sole 33cm	No 44 – 46
0807684	XL	Sole 36cm	No 47 +





LOW STANDARD WALKER

- Made of very lightweight plastic material
- Inner lining made of foam, antiperspirant and hypoallergenic
- All immobilizer straps are mounted on the splint making it easy to install
- Rocker bottom system on the foot and non-slip high density sole (Rubber)
- The brace can be used on either the left or the right leg

Indications:

- Stable fracture of the 5th metatarsal
- Transverse fracture of the 5th metatarsal
- Fractures of metatarsophalangeal joints
- Fracture of tarsal bones

REF	SHOE SIZE
L2411	XS No 34 – 36
	S No 36 – 39
	M No 40 – 43
	L No 44 – 46
	XL No 47 – 49





STD FRACTURE WALKER TALL

REF		SIZING	
0809091	S	Sole 25cm	No 36 – 39
0809092	M	Sole 29cm	No 40 – 43
0809093	L	Sole 33cm	No 44 – 46
0808081	XL	Sole 36cm	No 47 +





TALL STANDARD WALKER



- Made of very lightweight plastic material
- Inner lining made of foam, antiperspirant and hypoallergenic
- All immobilizer straps are mounted on the splint making it easy to install
- Rocker bottom system on the foot and non-slip high density sole (rubber)
- The brace can be used on either the left or the right leg

Indications:

- Conservative and postoperative in stable fracture of medial malleus
- Stable bicuspid fracture
- Type II and III sprain
- Stable tibial fibula fracture
- Fractured heel

REF	SHOE SIZE	
L2410	XS	No 34 - 36
	S	No 36 - 39
	M	No 40 - 43
	L	No 44 - 46
	XL	No 47 - 49



TALL PNEUMATIC WALKER "STABILITY"

- Made of lightweight plastic material
- It has built-in air chambers that offer maximum immobilization of the ankle in combination with the panels available in its front and rear part
- The lining (foam) is made of antiperspirant hypoallergenic material
- It has a rocker bottom system for better walking of the patient as well as a high density anti-slip sole (rubber)
- Easy air inflate-deflate system
- The brace can be used on either the left or the right leg

Indications:

- Conservatives in stable irregular fracture
- Irregular fracture of the hamlet
- Irregular stable fracture of the albulation
- Stable fracture of the talic head
- Tiofibular joint rupture
- Weber a type fracture kai b1
- Rupture of posterior tibial tendon
- Sublubric dislocation
- Postoperative for fracture of fibum
- Fracture of hamlet
- Bifire fracture
- Ligament sutures
- Median tibial-fibal membrane rupture
- Posterior tibial tendon rupture
- Pilon type fracture, maisonneuve, boswort

REF	SHOE SIZE	
L2436	S	No 36 – 39
	M	No 40 – 43
	L	No 44 – 47



ROM FRACTURE WALKER WITH AIR

- Features cells that inflate with air for further support when swelling subsides
- Equipped with mechanism for immobilization and adjustment of the range of motion, graduated in degrees (Range Of Motion)

REF	SHOE SIZE		
0809097	S	Sole 25cm	No 36 - 39
0809098	M	Sole 29cm	No 40 - 43
0809099	L	Sole 33cm	No 44 - 46



TALL PNEUMATIC WALKER WITH ROM

- It has aluminum side panels, antiperspirant and hypoallergenic lining with built-in air pump and valve, for maximum immobilization of the ankle
- Rocker bottom system with anti-slip high density rubber sole
- Goniometer to adjust the range of flexion and extension from 0 to 45 degrees
- The brace can be used on either the left or the right leg

Indications:

- Preservatives and post operative in achilles tendon suturing
- Amphisphia fractures
- Heel fracture
- Tibia fibula fracture
- Achilles tendon lengthening
- Suturing of medial and external collateral ligaments
- Weber b and c type fractures

REF	SHOE SIZE	
	S	No 36 – 39
L2409	M	No 40 – 43
	L	No 44 – 47





POST-OP SHOE

Indications:

Post-operatively in:

- Metatarsal osteotomies
- Finger osteotomies
- Morton's neuroma

- Made of ventilated non-elastic fabric.
- It has Velcro for easy closure and non-slip rigid sole.
- The brace can be used on either the left or the right leg

REF

SHOE SIZE

L2441

XS	No 34 - 36
S	No 37 - 39
M	No 40 - 42
L	No 41 - 43
XL	No 43 - 45
XXL	No 45 - 48





METATARSAL PRESSURE – RELIEF SHOE

Indications:

- Post-operative on the forefoot and toes
- Legs in plaster or with immobilizing gauze

- Construction of non-elastic, ventilated, hypoallergenic fabric
- Rigid, non-slip sole that leaves the metatarsal area free to discharge
- Closure with self-adhesive detachable Velcro
- The brace can be used on either the left or the right leg

REF	SHOE SIZE
L2417	XS No 35 – 37
	S No 37 – 39
	M No 40 – 42
	L No 42 – 44
	XL No 44 – 46





HEEL WEDGE PRESSURE- RELIEF SHOE



Indications:

- Heel soft tissue injuries
- Inflammatory condition in the heel area
- Fractured heel fracture
- Use after heel surgery

- Construction of non-elastic, ventilated, hypoallergenic fabric
- Rigid non-slip sole that leaves the heel area free to discharge
- Closure with self-adhesive detachable Velcro
- The brace can be used on either the left or the right leg

REF	SHOE SIZE	
L2415	XS	No 35 - 37
	S	No 37 - 39
	M	No 40 - 42
	L	No 42 - 44
	XL	No 44 - 46





WATERPROOF LEG COVER (ADULT)

REF	LENGTH
0808513	58.4cm

- Easy to apply
- Reusable waterproof PVC covering with plastic and elastic ring
- Convenient for use in shower or bath
- Intended for single patient use



WATERPROOF LEG COVER (CHILD)

REF	LENGTH
0808514	45.72cm



pressure therapy system

COMPRESSOR, ANALOGUE (4 PORTS)

Ref. 0806210



Specifications:

- Adjustable Pressure (Pressure range: 30 – 240 mmHg)
- Selection of Compression Therapy Duration
- Purchase of Suitable Pad (p. 113)
- One or two Pads can be Connected to the same Compressor (if it is suggested by the Doctor)
- To Connect to the Compressor two Pads, it is Necessary to Buy Ref: 0806221
- 2-year Warranty

COMPRESSOR, DIGITAL (4 PORTS)

Ref. 0806211



Specifications:

- Adjustable Pressure (Pressure range 20 – 250 mmHg)
- User-Selected Compression Therapy Duration
- 6 Programmes Available
- Purchase of Suitable Pads Necessary (p. 113)
- One or two Pads can be Connected to the same Compressor (if it is suggested by the Doctor)
- To Connect to the Compressor two Pads, it is Necessary to Buy Ref: 0806221
- 2-year Warranty

2-PAD CONNECTOR (4 CHAMBERS)

Ref. 0806221

Used to Connect two pads to the same Compressor if instructed by the Doctor.



pressure therapy system



REF WAIST COMPRESSION PAD

0806218 M: W: 24cm / L: 116cm

0806219 L: W: 36cm / L: 126cm

REF ARM COMPRESSION PAD

0806216 M: W: 19cm / L: 87cm

0806217 L: W: 25cm / L: 87cm

REF LEG COMPRESSION PAD

0806213 M: W: 24cm / L: 67cm / C: 62cm

0806214 L: W: 40cm / L: 67cm / C: 68cm

0806225 XL: W: 40cm / L: 86cm / C: 70cm





**foot
care**



SMALL TOE SEPARATOR WITH PROTECTIVE GEL




Indications:

- Bunionette /Tailor's Bunion

- Special design that allows for proper fit
- Hypoallergenic
- Minimizes friction
- Helps align the fingers
- Relieves the pain of the little finger bump
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2484	S	7.2cm	2cm
	L	8.5cm	2.5cm





HAMMERTOE SEPARATORS WITH PROTECTIVE GEL

Indications:


- Hammertoe



- Special design that allows for proper fit
- Provides comfort and relief from pressure on the toes and metatarsal heads
- Hypoallergenic
- Moisturizes the skin and massages the irritated area
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2485	ONE SIZE - UNIVERSAL	7.6cm	2.7cm






TOE RELAXING SEPARATORS WITH PROTECTIVE GEL

Indications:


- Toe relaxing



- Provides comfort and relief to the toes
- Easy to adjust
- Protection from friction
- Hypoallergenic
- Available in pairs

COMING SOON

REF	SIZE	A	B	C
L2487	S	6.8cm	3cm	1.9cm
	L	8.8cm	3.8cm	2.4cm





TOE RING SEPARATORS WITH GEL

Indications:

- Hallux valgus

Small




Large



- Provides a stable fit, preventing slipping during movement
- Helps align fingers in cases of deviation
- Reduces friction on fingers
- Hypoallergenic
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2488	S (4th- 5th toe)	3.3cm	1.4cm
	L (1st-2nd toe)	4.9cm	2.3cm





1st-2nd TOE RING SEPARATORS WITH PROTECTIVE GEL ONE SIZE




Indications:

- Hallux valgus

- Made of high quality and durable Gel
- Helps align the fingers
- Relieves the irritated area
- Hypoallergenic
- Ambidextrous
- One-size
- Available in pairs

COMING SOON

REF	SIZE	A	B	C
L2489	ONE SIZE – UNIVERSAL	6.8cm	4.7cm	2cm





1st-2nd TOE SEPARATORS WITH PROTECTIVE GEL

Indications:


- Hallux valgus



- Hypoallergenic
- Made of high quality Gel polymer, to separate the big and second toe
- Reduces friction and irritation
- Available in pairs

COMING SOON

REF	SIZE	A	B	C
L2490	S	3cm	1.8cm	2.2cm
	L	3.5cm	2.2cm	2.4cm





BIG TOE SEPARATOR WITH PROTECTIVE GEL




Indications:

- Hallux valgus

- Protects against calluses by minimizing friction
- Helps align the toes and relieves the pain of the big toe bump
- Special design for proper fit
- Hypoallergenic
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2502	ONE SIZE – UNIVERSAL	8cm	2.5cm





PROTECTIVE SOCK WITH GEL FOR HALLUX VALGUS

Indications:

- Hallux valgus



- Made of elastic and soft fabric
- Reduces friction
- Hypoallergenic
- Comfortable fit on all foot sizes
- Washable
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2491	S	9.5cm	8cm
	L	10.2cm	9cm





METATARSAL PROTECTIVE GEL

Indications:


- Metatarsalgia



- Special design that allows for proper fit
- Provides complete comfort when walking on the metatarsal
- Absorbs shock
- Hypoallergenic
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2486	ONE SIZE - UNIVERSAL	8.7cm	10cm





METATARSAL SOCK WITH GEL PAD


Indications:

- Metatarsalgia

- Hypoallergenic
- Protects metatarsal heads
- Absorbs shock
- Reduces friction
- Washable
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2493	S	9cm	7.8cm
	M	10cm	7.8cm





PROTECTIVE TOE TUBES WITH PROTECTIVE GEL


Indications:

- Protection of toes from calluses

- Can be cut to desired length
- Suitable for interdigital calluses, calluses and hammertoes
- Hypoallergenic
- Washable
- Easy to apply
- Available per piece


COMING SOON

REF	SIZE	A	B
L2492	S	10cm	1.5cm
	L	10cm	1.7cm





HALLUX VALGUS BRACE




Indications:

- Hallux valgus

- Made of high quality, versatile plastic material
- High-strength inner foam lining with anti-allergic, anti-perspirant, non-toxic action
- Washable
- Stabilizing adjustable adhesive fastening strap
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2503	ONE SIZE - UNIVERSAL	13.5cm	5cm



HALLUX VALGUS BRACE WITH PROTECTIVE GEL




Indications:

- Hallux valgus

- Hypoallergenic
- Corrects deformities of the big toe on the foot
- Ideal for those with bunions, calluses, glued or crossed toes
- Stretches the big toe and helps it until its normal shape is restored
- With soft Gel in contact with the skin
- Made with an ergonomic design, it is easily fixed and does not press on the toe
- Functional articulated design, according to the shape of the foot that helps restore the toe to its normal position
- Cleaned using lukewarm water and hand soap
- Available per piece

COMING SOON

REF	SIZE	A	B
L2494	ONE SIZE – UNIVERSAL	15cm	adjustable





ORTHOTIC - ANATOMICAL SILICONE INSOLES WITH GEL

Indications:


- Morton's neuroma
- Flat feet
- Heel spurs
- Metatarsalgia



- High elasticity & durability
- They offer comfort & provide uniform pressure distribution
- Balancing the forces exerted on the body while walking
- Facilitating daily activities
- Made of Silicone Gel for shock absorption
- Perforated surface for smooth air flow
- Available in pairs

COMING SOON

REF	SIZE	A	B
	S (35-38)	23cm	8cm
L2495	M (39-42)	28cm	9.5cm
	L (43-45)	29.5cm	10cm



FLATFOOT INSOLES WITH FOAM

Indications:


- Flatfoot



- They are suitable for use in various types of footwear, such as sports shoes or work shoes.
- They have memory foam on the upper part for better shock absorption
- Perforated surface on the bottom for smooth air flow
- They have arch support to relieve pain and correct the position of the foot in cases of flat feet
- They offer comfort & provide uniform pressure distribution
- Available in pairs

COMING SOON

REF	SIZE	A	B
	XS (36-38)	23.5cm	8.7cm
	S (39-40)	26.5cm	9.4cm
L2498	M (41-42)	28cm	9.8cm
	L (43-45)	30cm	10.3cm
	XL (46-47)	30.7cm	10.3cm





ANATOMICAL INSOLES WITH FOAM


Indications:

- Shock absorption
- Flatfoot

- They have memory foam on the upper for better shock absorption
- They offer support to the arch of the foot
- Made of antibacterial materials that allow the product to breathe
- Ideal for walking and extended standing
- Hypoallergenic
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2500	XS (36-38)	25cm	9cm
	S (39-40)	26.5cm	9.3cm
	M (41-43)	28.3cm	9.7cm
	L (44-45)	30.5cm	10.2cm



ATHLETIC INSOLES WITH GEL AND FOAM

Indications:


- Shock absorption



- Made of thermoplastic, elastic GEL layer, anatomically designed and highly shock-absorbing
- They have absorbent pads that offer continuous micro-massages
- Specially designed in the heel area
- Improves pain
- Can be cut to the desired length
- Use during sports, work
- Washable
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2499	S (35-40)	27.7cm	10cm
	L (41-46)	30.5cm	10.5cm





COMFORT INSOLES WITH FOAM


Indications:

- Shock absorption

- They are suitable for use in various types of footwear, such as sports shoes or work shoes.
- They have memory foam on the upper for better shock absorption
- Perforated surface on the bottom for smooth air flow
- They have arch support to relieve pain and correct the position of the foot
- They offer comfort & provide even pressure distribution
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2501	S (38-41)	26.5cm	9.5cm
	L (42-45)	29.5cm	10cm





SILICONE HEEL CUSHION PADS WITH REMOVABLE DOT GEL


Indications:

- Heel spur

- Suitable for all types of shoes ideal for hiking, walking, running, playing and other outdoor sports
- Hand washable & reusable
- Made of silicone gel for shock absorption
- Features a soft removable insert for bunions
- Protects the heel and provides even pressure distribution
- Available in pairs

COMING SOON

REF	SIZE	A	B
	S	8.7cm	5.8cm
L2496	M	9.3cm	6cm
	L	10.2cm	6.4cm





GEL HEEL CUSHION PADS


Indications:

- Heel spur

- Suitable for all types of shoes ideal for hiking, walking, running, playing and other outdoor sports
- Hand washable & reusable
- Made of Gel for shock absorption
- Protects the heel and provides even pressure distribution
- Available in pairs

COMING SOON

REF	SIZE	A	B
L2497	S	10.4cm	6.3cm
	L	11cm	7cm





PROTECTIVE SOCK WITH GEL FOR ACHILLES AND HEEL

Indications:

- Protection of tendons in the Achilles part

- Easy to apply
- Pressure relief
- Hypoallergenic
- Available per piece
- Washable

COMING SOON

REF	SIZE	A	B
L2510	S	20cm	11cm
	M	20cm	11.5cm







Kinesiology



SPORT TAPE

Indications:

- These products are used as sports bandages for athletes, suitable for additional support, but also for the prevention of injuries to joints, tendons and ligaments, mainly during sports activities.

- Adhesive, long-lasting
- Made of cotton
- Particularly durable and non-elastic
- Permeable to water vapor and air, so it does not easily come off due to sweating
- Easily cut
- Hypoallergenic
- Latex-free
- Can be used directly on the skin or over foam under wrap sport tape (spider)

REF	DIMENSIONS (WxL)
0815750	2.5cm x 13.7m
0815751	3.8cm x 13.7m
0815752	5cm x 13.7m



UNDER WRAP SPORT TAPE (SPIDER)

Indications:

- They are used as sports tapes for athletes, suitable for use under adhesive sports tapes to prevent skin irritation from repeated use of adhesive tapes. They also help in the painless removal of adhesive sports tapes. Recommended for use by people with sensitive skin.

- Made of 100% polyurethane foam (PU)
- It fits perfectly on the skin without containing glue and leaves no residue on the skin when removed
- It is used under the bandage with adhesive sports tape & helps in its painless detachment
- High quality construction & comfortable fit
- It is easily & accurately cut to the desired length
- Does not contain Latex
- Permeable to water vapor and air

REF	DIMENSIONS (WxL)	THICKNESS
0815753	7cm x 27.4m	0.6mm



KINESIOLOGY TAPE

Indications:

• These products are intended to promote functional movement of the body at an early stage during or after injuries, as clinical studies have shown that this leads to faster and better healing. They are used to speed up the healing of muscle injuries and support them. They are strategically applied to the body to provide support, reduce pain and swelling from injuries and improve performance. For applications in rehabilitation, sports and in cases of edema.

- Suitable for a wide range of applications in the fields of sports medicine and physiotherapy
- Adhesive
- Particularly durable and elastic
- Hypoallergenic
- Made of cotton
- Does not contain Latex

REF	DIMENSIONS (WxL)	COLOUR
0815754	5cm x 5m	BLACK
0815755	5cm x 5m	BEIGE
0815756	5cm x 5m	BLUE
0815757	5cm x 5m	PINK



www.mobiakcare.com

E: medicalexports@mobiak.com